



Hydrophobic Interaction Chromatography (HIC) of Proteins on a New Series of HIC Ligands that Retain Samples with Less Salt for Direct Analysis via Mass Spectrometry (MS) with Native Structures Intact.

Increased HIC Ligand Size Reduces Salt Required for Retention



Overview

Our new series of materials for HIC permit the retention of proteins using concentrations of ammonium acetate compatible with direct analysis via MS. Proteins can be separated and eluted with their native structures intact. This is useful for top-down proteomics and permits the analysis of some proteins that are not compatible with conditions of reversed-phase chromatography. The example to the left illustrates the reduction in salt achieved by increasing the ligand size in our new materials. Here we effectively reduced the salt required for retention from 1M to 0.4M Am-OAc.




Parameters:

Buffer A: 1.0M Ammonium Acetate
Buffer B: 25mM Ammonium Acetate, 50% ACN
Gradient: 0-100%B, 15min

HIC-MS of Antibodies


Sample: NIST mAb standard 8671 (IgG1 κ)
Column: PolyPENTYL A™ capillary, 100x0.2-mm, 2- μ m, 1000- Å
Mobile Phase: Decreasing gradient of NH₄-Oac
Flow rate: 2.6 μ l/min
MS: maXis II QTOF (Bruker)

PolyPENTYL A™
Base Peak Chromatogram



Theoretical: 148037.4 Da
Experimental: 148037.5 Da

MS



The composition of this antibody's glycoform variants is consistent with the literature: M. Hilliard *et al.*, *mAbs* 9 (2017) 1349.

The mass spectrum is that of a protein with native structure.

Common	Hexose	Hexosamine
G0	Fuc	GalNAc
G1	Rha	GlcNAc
G2	Qui	ManNAc
G0F	Gal	
G1F	Glc	
G2F	Man	

■ -203 Da (-GlcNAc) (data courtesy of Bifan Chen, Ziqing Lin, & Ying Ge, Univ. Wisconsin-Madison)
* +162 Da