

BioPureSPN™ PolyHYDROXYETHYL A™, PolySULFOETHYL A™, PolyWAX LP™, PolyCAT A® & Primary-Secondary Amine (PSA™ ERLIC) MACRO SPE Spin Columns (≈0.2-1mg maximum IEX capacity, ≈35-350µg max HILIC capacity) [Protocols](#) MACRO HILIC or IEX

Nest P/N	Price	For Bench Top Centrifuges	PolyLC Mfg P/N
HMM HIL	\$890	BioPureSPN MACRO , 100mg, 300Å, PolyHYDROXYETHYL A HILIC 35-350µg capacity, 50-200µL, 100/bx Protocols MACRO HILIC	SPNA-1CMG-HEA2
HMM HIL.20	\$250	BioPureSPN MACRO , 100mg, 300Å, PolyHYDROXYETHYL A HILIC 35-350µg capacity, 50-200µL, 20/bx	
HNS HIL-L	\$760	BioPureSPN 96-Well MACRO Plate , 100mg, 300Å, PolyHYDROXYETHYL A HILIC 35-350µg capacity, 50-200µL, ea. Protocols MACRO HILIC	96WA-1CMG-HEA2
HNS HIL60-L	\$760	BioPureSPN 96-Well MACRO Combinatorial Library Screening Plate , 100mg, 60Å, PolyHYDROXYETHYL A SEC 10-20µL, ea. Protocols SEC-60	96WC-1CMG-HIL6
HMM HIL-SCX	\$890	BioPureSPN MACRO , 100mg, 300Å, PolySULFOETHYL A HILIC-SCX 0.2-1mg/35-350µg capacity, 50-200µL, 100/bx	SPNA-1CMG-SCX2
HMM HIL-SCX.20	\$250	BioPureSPN MACRO , 100mg, 300Å, PolySULFOETHYL A HILIC-SCX 0.2-1mg/35-350µg capacity, 50-200µL, 20/bx	
HNS HIL-SCX-L	\$760	BioPureSPN 96-Well MACRO Plate , 100mg, 300Å, PolySULFOETHYL A HILIC-SCX 0.2-1mg/35-350µg capacity, 50-200µL, ea.	96WA-1CMG-SCX2
HMM HIL-PSA	\$490	BioPureSPN MACRO , 100mg, 120Å, "PSA" (Primary Secondary Amine) ERLIC-WAX 0.2-1mg/35-350µg capacity, 50-200µL, 100/bx	SPNC-1CMG-PSA2
HMM HIL-PSA.20	\$150	BioPureSPN MACRO , 100mg, 120Å, "PSA" (Primary Secondary Amine) ERLIC-WAX 0.2-1mg/35-350µg capacity, 50-200µL, 20/bx	
HNS HIL-PSA-L	\$460	BioPureSPN 96-Well MACRO Plate , 100mg, 120Å, "PSA" (Primary Secondary Amine) ERLIC-WAX 0.2-1mg/35-350µg capacity, 50-200µL, ea.	96WC-1CMG-PSA2
HMM HIL-WAX	\$890	BioPureSPN MACRO , 100mg, 300Å, PolyWAX LP HILIC-WAX 0.2-1mg/35-350µg capacity, 50-200µL, 100/bx	SPNA-1CMG-WAX2
HMM HIL-WAX.20	\$250	BioPureSPN MACRO , 100mg, 300Å, PolyWAX LP HILIC-WAX 0.2-1mg/35-350µg capacity, 50-200µL, 20/bx	
HNS HIL-WAX-L	\$890	BioPureSPN 96-Well MACRO Plate , 100mg, 300Å, PolyWAX LP HILIC-WAX 0.2-1mg/35-350µg capacity, 50-200µL, ea.	96WA-1CMG-WAX2
HMM HIL-WCX	\$890	BioPureSPN MACRO , 100mg, 300Å, PolyCAT A HILIC-WCX 0.2-1mg/35-350µg capacity, 50-200µL, 100/bx	SPNA-1CMG-WCX2
HMM HIL-WCX.20	\$250	BioPureSPN MACRO , 100mg, 300Å, PolyCAT A HILIC-WCX 0.2-1mg/35-350µg capacity, 50-200µL, 20/bx	
HNS HIL-WCX-L	\$760	BioPureSPN 96-Well MACRO Plate , 100mg, 300Å, PolyCAT A HILIC-WCX 0.2-1mg/35-350µg capacity, 50-200µL, ea.	96WA-1CMG-WCX2

BioPureSPN™ PolyHYDROXYETHYL A™, PolySULFOETHYL A™, PolyWAX LP™, & PolyCAT A® & Primary-Secondary Amine (PSA™ ERLIC) MIDI SPE Spin Columns (≈100-500µg maximum IEX capacity, ≈17-170µg max HILIC capacity) [Protocols](#) MIDI HILIC or IEX

Nest P/N	Price	For Bench Top Centrifuges	PolyLC Mfg P/N
HEM HIL	\$780	BioPureSPN Midi , 50mg, 300Å, PolyHYDROXYETHYL A HILIC 17-170µg capacity, 50-200µL, 100/bx Protocols MIDI HILIC	SPNA-50MG-HEA2
HEM HIL.20	\$220	BioPureSPN Midi , 50mg, 300Å, PolyHYDROXYETHYL A HILIC 17-170µg capacity, 50-200µL, 20/bx	
HNS HIL-M	\$550	BioPureSPE 96-Well Midi Plate , 50mg, 300Å, PolyHYDROXYETHYL A HILIC 17-170µg capacity, 50-200µL, ea. Protocols MIDI HILIC	96WA-50MG-HEA2
HEM HIL-SCX	\$780	BioPureSPN Midi , 50mg, 300Å, PolySULFOETHYL A HILIC-SCX 100-500µg/17-170µg capacity, 50-200µL, 100/bx	SPNA-50MG-SCX2
HEM HIL-SCX.20	\$220	BioPureSPN Midi , 50mg, 300Å, PolySULFOETHYL A HILIC-SCX 100-500µg/17-170µg capacity, 50-200µL, 20/bx	
HNS HIL-SCX-M	\$550	BioPureSPN 96-Well Midi Plate , 50mg, 300Å, PolySULFOETHYL A HILIC-SCX 100-500µg/17-170µg capacity, 50-200µL, ea.	96WA-50MG-SCX2
HEM HIL-PSA	\$460	BioPureSPN Midi , 50mg, 120Å, "PSA" (Primary Secondary Amine) ERLIC-WAX 100-500µg/17-170µg capacity, 50-200µL, 100/bx	SPNA-50MG-PSA2
HEM HIL-PSA.20	\$145	BioPureSPN Midi , 50mg, 120Å, "PSA" (Primary Secondary Amine) ERLIC-WAX 100-500µg/17-170µg capacity, 50-200µL, 20/bx	
HNS HIL-PSA-M	\$430	BioPureSPN Midi Plate , 50mg, 120Å, "PSA" (Primary Secondary Amine) ERLIC-WAX 100-500µg/17-170µg capacity, 50-200µL, ea.	96WA-50MG-PSA2
HEM HIL-WAX	\$780	BioPureSPN Midi , 50mg, 300Å, PolyWAX LP HILIC-WAX 100-500µg/17-170µg capacity, 50-200µL, 100/bx	SPNA-50MG-WAX2
HEM HIL-WAX.20	\$220	BioPureSPN Midi , 50mg, 300Å, PolyWAX LP HILIC-WAX 100-500µg/17-170µg capacity, 50-200µL, 20/bx	
HNS HIL-WAX-M	\$550	BioPureSPN 96-Well Midi Plate , 50mg, 300Å, PolyWAX LP HILIC-WAX 100-500µg/17-170µg capacity, 50-200µL, ea.	96WA-50MG-WAX2
HEM HIL-WCX	\$780	BioPureSPN Midi , 50mg, 300Å, PolyCAT A HILIC-WCX 100-500µg/17-170µg capacity, 50-200µL, 100/bx	SPNA-50MG-WCX2
HEM HIL-WCX.20	\$220	BioPureSPN Midi , 50mg, 300Å, PolyCAT A HILIC-WCX 100-500µg/17-170µg capacity, 50-200µL, 20/bx	
HNS HIL-WCX-M	\$550	BioPureSPN 96-Well Midi Plate , 50mg, 300Å, PolyCAT A HILIC-WCX 100-500µg/17-170µg capacity, 50-200µL, ea.	96WA-50MG-WCX2

BioPureSPN™ PolyHYDROXYETHYL A™, PolySULFOETHYL A™, PolyWAX LP™, & PolyCAT A® & Primary-Secondary Amine (PSA™ ERLIC) Mini SPE Spin Columns (≈40-200μg maximum IEX capacity, ≈7-70μg max HILIC capacity) [Protocols](#) Mini HILIC or IEX

Nest P/N	Price	For Bench Top Centrifuges	PolyLC Mfg P/N
HUM HIL	\$700	BioPureSPN Mini , 20mg, 300Å, PolyHYDROXYETHYL A HILIC 7-70μg capacity, 20-100μL, 100/bx Protocols Mini HILIC	SPNA-20MG-HEA2
HUM HIL.20	\$210	BioPureSPN Mini , 20mg, 300Å, PolyHYDROXYETHYL A HILIC 7-70μg capacity, 20-100μL, 20/bx	
HNS HIL	\$525	BioPureSPE 96-Well Mini Plate , 20mg, 300Å, PolyHYDROXYETHYL A HILIC 7-70μg capacity, 20-100μL, ea. Protocols Mini HILIC	96WA-20MG-HEA2
HUM HIL-SCX	\$700	BioPureSPN Mini , 20mg, 300Å, PolySULFOETHYL A HILIC-SCX 40-200μg/7-70μg capacity, 20-100μL, 100/bx	SPNA-20MG-SCX2
HUM HIL-SCX.20	\$210	BioPureSPN Mini , 20mg, 300Å, PolySULFOETHYL A HILIC-SCX 40-200μg/7-70μg capacity, 20-100μL, 20/bx	
HNS HIL-SCX	\$525	BioPureSPN 96-Well Mini Plate , 20mg, 300Å, PolySULFOETHYL A HILIC-SCX 40-200μg/7-70μg capacity, 20-100μL, ea	96WA-20MG-SCX2
HUM HIL-PSA	\$440	BioPureSPN Mini , 20mg, 120Å, "PSA" (Primary Secondary Amine) ERLIC-WAX 40-200μg/7-70μg capacity, 20-100μL, 100/bx	SPNA-20MG-PSA2
HUM HIL-PSA.20	\$140	BioPureSPN Mini , 20mg, 120Å, "PSA" (Primary Secondary Amine) ERLIC-WAX 40-200μg/7-70μg capacity, 20-100μL, 20/bx	
HNS HIL-PSA	\$410	BioPureSPN Mini Plate , 20mg, 120Å, "PSA" (Primary Secondary Amine) ERLIC-WAX 40-200μg/7-70μg capacity, 20-100μL, ea.	96WA-20MG-PSA2
HUM HIL-WAX	\$700	BioPureSPN Mini , 20mg, 300Å, PolyWAX LP ERLIC-WAX 40-200μg/7-70μg capacity, 20-100μL, 100/bx	SPNA-20MG-WAX2
HUM HIL-WAX.20	\$210	BioPureSPN Mini , 20mg, 300Å, PolyWAX LP ERLIC-WAX 40-200μg/7-70μg capacity, 20-100μL, 20/bx	
HNS HIL-WAX	\$700	BioPureSPN 96-Well Mini Plate , 20mg, 300Å, PolyWAX LP ERLIC-WAX 40-200μg/7-70μg capacity, 20-100μL, ea.	96WA-20MG-WAX2
HUM HIL-WCX	\$700	BioPureSPN Mini , 20mg, 300Å, PolyCAT A HILIC-WCX 40-200μg/7-70μg capacity, 20-100μL, 100/bx	SPNA-20MG-WCX2
HUM HIL-WCX.20	\$210	BioPureSPN Mini , 20mg, 300Å, PolyCAT A HILIC-WCX 40-200μg/7-70μg capacity, 20-100μL, 20/bx	
HNS HIL-WCX	\$700	BioPureSPN 96-Well Mini Plate , 20mg, 300Å, PolyCAT A HILIC-WCX 40-200μg/7-70μg capacity, 20-100μL, ea	96WA-20MG-WCX2