

2026 Product Pricelist

Specializing in manufacturing alternatives to reversed-phase HPLC for proteins, peptides, and other biomolecules since 1985, with a focus on the most challenging separations.

PolyLC Inc.
9151 Rumsey Rd.
Suites 175 & 180
Columbia, MD 21045
(410) 992-5400
www.polylc.com



2026 PolyLC Inc. Retail Pricelist



PolyCAT A®

(Weak-Cation Exchange)



Offering high-performance liquid chromatography media, columns, capillaries, cartridges, and instrumentation. We specialize in advanced alternatives to reversed-phase for proteins, peptides, and other biomolecules, supported by deep application expertise and customized solutions for the most demanding separations.

HPLC Columns

PolyCAT A® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore size options.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
051CT0203	PolyCAT A® Column, Glass Lined SS, 50 x 1.0mm, 2µm,300Å	WCX	Analytical	GLT	10-32 coned	50	1.0	2	300	615
051CT0210	PolyCAT A® Column, Glass Lined SS, 50 x 1.0mm, 2µm,1000Å	WCX	Analytical	GLT	10-32 coned	50	1.0	2	1000	615
051CT0215	PolyCAT A® Column, Glass Lined SS, 50 x 1.0mm, 2µm,1500Å	WCX	Analytical	GLT	10-32 coned	50	1.0	2	1500	615
101CT0203	PolyCAT A® Column, Glass Lined SS, 100 x 1.0mm, 2µm,300Å	WCX	Analytical	GLT	10-32 coned	100	1.0	2	300	750
101CT0210	PolyCAT A® Column, Glass Lined SS, 100 x 1.0mm, 2µm,1000Å	WCX	Analytical	GLT	10-32 coned	100	1.0	2	1000	750
101CT0215	PolyCAT A® Column, Glass Lined SS, 100 x 1.0mm, 2µm,1500Å	WCX	Analytical	GLT	10-32 coned	100	1.0	2	1500	750
151CT0203	PolyCAT A® Column, Glass Lined SS, 150 x 1.0mm, 2µm,300Å	WCX	Analytical	GLT	10-32 coned	150	1.0	2	300	800
151CT0210	PolyCAT A® Column, Glass Lined SS, 150 x 1.0mm, 2µm,1000Å	WCX	Analytical	GLT	10-32 coned	150	1.0	2	1000	800
151CT0215	PolyCAT A® Column, Glass Lined SS, 150 x 1.0mm, 2µm,1500Å	WCX	Analytical	GLT	10-32 coned	150	1.0	2	1500	800
3.52CT0203	PolyCAT A® Column, SS, 35 x 2.1mm, 2µm,300Å	WCX	Analytical	SS	10-32 coned	35	2.1	2	300	570
3.52CT0210	PolyCAT A® Column, SS, 35 x 2.1mm, 2µm,1000Å	WCX	Analytical	SS	10-32 coned	35	2.1	2	1000	570
3.52CT0215	PolyCAT A® Column, SS, 35 x 2.1mm, 2µm,1500Å	WCX	Analytical	SS	10-32 coned	35	2.1	2	1500	570
052CT0203	PolyCAT A® Column, SS, 50 x 2.1mm, 2µm,300Å	WCX	Analytical	SS	10-32 coned	50	2.1	2	300	580
052CT0210	PolyCAT A® Column, SS, 50 x 2.1mm, 2µm,1000Å	WCX	Analytical	SS	10-32 coned	50	2.1	2	1000	580
052CT0215	PolyCAT A® Column, SS, 50 x 2.1mm, 2µm,1500Å	WCX	Analytical	SS	10-32 coned	50	2.1	2	1500	580
102CT0203	PolyCAT A® Column, SS, 100 x 2.1mm, 2µm,300Å	WCX	Analytical	SS	10-32 coned	100	2.1	2	300	750
102CT0210	PolyCAT A® Column, SS, 100 x 2.1mm, 2µm,1000Å	WCX	Analytical	SS	10-32 coned	100	2.1	2	1000	750
102CT0215	PolyCAT A® Column, SS, 100 x 2.1mm, 2µm,1500Å	WCX	Analytical	SS	10-32 coned	100	2.1	2	1500	750
3.54CT0203	PolyCAT A® Column, SS, 35 x 4.6mm, 2µm,300Å	WCX	Analytical	SS	10-32 coned	35	4.6	2	300	570
3.54CT0210	PolyCAT A® Column, SS, 35 x 4.6mm, 2µm,1000Å	WCX	Analytical	SS	10-32 coned	35	4.6	2	1000	570
3.54CT0215	PolyCAT A® Column, SS, 35 x 4.6mm, 2µm,1500Å	WCX	Analytical	SS	10-32 coned	35	4.6	2	1500	570
054CT0203	PolyCAT A® Column, SS, 50 x 4.6mm, 2µm,300Å	WCX	Analytical	SS	10-32 coned	50	4.6	2	300	580
054CT0210	PolyCAT A® Column, SS, 50 x 4.6mm, 2µm,1000Å	WCX	Analytical	SS	10-32 coned	50	4.6	2	1000	580
054CT0215	PolyCAT A® Column, SS, 50 x 4.6mm, 2µm,1500Å	WCX	Analytical	SS	10-32 coned	50	4.6	2	1500	580
104CT0203	PolyCAT A® Column, SS, 100 x 4.6mm, 2µm,300Å	WCX	Analytical	SS	10-32 coned	100	4.6	2	300	750
104CT0210	PolyCAT A® Column, SS, 100 x 4.6mm, 2µm,1000Å	WCX	Analytical	SS	10-32 coned	100	4.6	2	1000	750
104CT0215	PolyCAT A® Column, SS, 100 x 4.6mm, 2µm,1500Å	WCX	Analytical	SS	10-32 coned	100	4.6	2	1500	750
154CT0203	PolyCAT A® Column, SS, 150 x 4.6mm, 2µm,300Å	WCX	Analytical	SS	10-32 coned	150	4.6	2	300	790
154CT0210	PolyCAT A® Column, SS, 150 x 4.6mm, 2µm,1000Å	WCX	Analytical	SS	10-32 coned	150	4.6	2	1000	790
154CT0215	PolyCAT A® Column, SS, 150 x 4.6mm, 2µm,1500Å	WCX	Analytical	SS	10-32 coned	150	4.6	2	1500	790
204CT0203	PolyCAT A® Column, SS, 200 x 4.6mm, 2µm,300Å	WCX	Analytical	SS	10-32 coned	200	4.6	2	300	810
204CT0210	PolyCAT A® Column, SS, 200 x 4.6mm, 2µm,1000Å	WCX	Analytical	SS	10-32 coned	200	4.6	2	1000	810
204CT0215	PolyCAT A® Column, SS, 200 x 4.6mm, 2µm,1500Å	WCX	Analytical	SS	10-32 coned	200	4.6	2	1500	810

PolyCAT A® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore size options.) Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
109CT0203	PolyCAT A® Column, SS, 100 x 9.4mm, 2µm,300Å	WCX	Semi-Prep	SS	10-32 coned	100	9.4	2	300	1525
109CT0210	PolyCAT A® Column, SS, 100 x 9.4mm, 2µm,1000Å	WCX	Semi-Prep	SS	10-32 coned	100	9.4	2	1000	1525
109CT0215	PolyCAT A® Column, SS, 100 x 9.4mm, 2µm,1500Å	WCX	Semi-Prep	SS	10-32 coned	100	9.4	2	1500	1525
209CT0203	PolyCAT A® Column, SS, 200 x 9.4mm, 2µm,300Å	WCX	Semi-Prep	SS	10-32 coned	200	9.4	2	300	2450
209CT0210	PolyCAT A® Column, SS, 200 x 9.4mm, 2µm,1000Å	WCX	Semi-Prep	SS	10-32 coned	200	9.4	2	1000	2450
209CT0215	PolyCAT A® Column, SS, 200 x 9.4mm, 2µm,1500Å	WCX	Semi-Prep	SS	10-32 coned	200	9.4	2	1500	2450
051CT0303	PolyCAT A® Column, Glass Lined SS, 50 x 1.0mm, 3µm,300Å	WCX	Analytical	GLT	10-32 coned	50	1.0	3	300	570
051CT0310	PolyCAT A® Column, Glass Lined SS, 50 x 1.0mm, 3µm,1000Å	WCX	Analytical	GLT	10-32 coned	50	1.0	3	1000	570
051CT0315	PolyCAT A® Column, Glass Lined SS, 50 x 1.0mm, 3µm,1500Å	WCX	Analytical	GLT	10-32 coned	50	1.0	3	1500	570
101CT0303	PolyCAT A® Column, Glass Lined SS, 100 x 1.0mm, 3µm,300Å	WCX	Analytical	GLT	10-32 coned	100	1.0	3	300	660
101CT0310	PolyCAT A® Column, Glass Lined SS, 100 x 1.0mm, 3µm,1000Å	WCX	Analytical	GLT	10-32 coned	100	1.0	3	1000	660
101CT0315	PolyCAT A® Column, Glass Lined SS, 100 x 1.0mm, 3µm,1500Å	WCX	Analytical	GLT	10-32 coned	100	1.0	3	1500	660
151CT0303	PolyCAT A® Column, Glass Lined SS, 150 x 1.0mm, 3µm,300Å	WCX	Analytical	GLT	10-32 coned	150	1.0	3	300	740
151CT0310	PolyCAT A® Column, Glass Lined SS, 150 x 1.0mm, 3µm,1000Å	WCX	Analytical	GLT	10-32 coned	150	1.0	3	1000	740
151CT0315	PolyCAT A® Column, Glass Lined SS, 150 x 1.0mm, 3µm,1500Å	WCX	Analytical	GLT	10-32 coned	150	1.0	3	1500	740
3.52CT0303	PolyCAT A® Column, SS, 35 x 2.1mm, 3µm,300Å	WCX	Analytical	SS	10-32 coned	35	2.1	3	300	550
3.52CT0310	PolyCAT A® Column, SS, 35 x 2.1mm, 3µm,1000Å	WCX	Analytical	SS	10-32 coned	35	2.1	3	1000	550
3.52CT0315	PolyCAT A® Column, SS, 35 x 2.1mm, 3µm,1500Å	WCX	Analytical	SS	10-32 coned	35	2.1	3	1500	550
052CT0303	PolyCAT A® Column, SS, 50 x 2.1mm, 3µm,300Å	WCX	Analytical	SS	10-32 coned	50	2.1	3	300	560
052CT0310	PolyCAT A® Column, SS, 50 x 2.1mm, 3µm,1000Å	WCX	Analytical	SS	10-32 coned	50	2.1	3	1000	560
052CT0315	PolyCAT A® Column, SS, 50 x 2.1mm, 3µm,1500Å	WCX	Analytical	SS	10-32 coned	50	2.1	3	1500	560
102CT0303	PolyCAT A® Column, SS, 100 x 2.1mm, 3µm,300Å	WCX	Analytical	SS	10-32 coned	100	2.1	3	300	710
102CT0310	PolyCAT A® Column, SS, 100 x 2.1mm, 3µm,1000Å	WCX	Analytical	SS	10-32 coned	100	2.1	3	1000	710
102CT0315	PolyCAT A® Column, SS, 100 x 2.1mm, 3µm,1500Å	WCX	Analytical	SS	10-32 coned	100	2.1	3	1500	710
024.OCT0315	PolyCAT A® Column, SS, 20 x 4.6mm, 3µm,300Å (Ultra fast Hb Variant Screening)	WCX	Analytical	SS	10-32 coned	35	4.0	3	1500	350
3.54CT0303	PolyCAT A® Column, SS, 35 x 4.6mm, 3µm,300Å	WCX	Analytical	SS	10-32 coned	35	4.6	3	300	545
3.54CT0310	PolyCAT A® Column, SS, 35 x 4.6mm, 3µm,1000Å	WCX	Analytical	SS	10-32 coned	35	4.6	3	1000	545
3.54CT0315	PolyCAT A® Column, SS, 35 x 4.6mm, 3µm,1500Å	WCX	Analytical	SS	10-32 coned	35	4.6	3	1500	545
054CT0303	PolyCAT A® Column, SS, 50 x 4.6mm, 3µm,300Å	WCX	Analytical	SS	10-32 coned	50	4.6	3	300	560
054CT0310	PolyCAT A® Column, SS, 50 x 4.6mm, 3µm,1000Å	WCX	Analytical	SS	10-32 coned	50	4.6	3	1000	560
054CT0315	PolyCAT A® Column, SS, 50 x 4.6mm, 3µm,1500Å	WCX	Analytical	SS	10-32 coned	50	4.6	3	1500	560
104CT0303	PolyCAT A® Column, SS, 100 x 4.6mm, 3µm,300Å	WCX	Analytical	SS	10-32 coned	100	4.6	3	300	710
104CT0310	PolyCAT A® Column, SS, 100 x 4.6mm, 3µm,1000Å	WCX	Analytical	SS	10-32 coned	100	4.6	3	1000	710
104CT0315	PolyCAT A® Column, SS, 100 x 4.6mm, 3µm,1500Å	WCX	Analytical	SS	10-32 coned	100	4.6	3	1500	710
154CT0303	PolyCAT A® Column, SS, 150 x 4.6mm, 3µm,300Å	WCX	Analytical	SS	10-32 coned	150	4.6	3	300	740
154CT0310	PolyCAT A® Column, SS, 150 x 4.6mm, 3µm,1000Å	WCX	Analytical	SS	10-32 coned	150	4.6	3	1000	740
154CT0315	PolyCAT A® Column, SS, 150 x 4.6mm, 3µm,1500Å	WCX	Analytical	SS	10-32 coned	150	4.6	3	1500	740
204CT0303	PolyCAT A® Column, SS, 200 x 4.6mm, 3µm,300Å	WCX	Analytical	SS	10-32 coned	200	4.6	3	300	760
204CT0310	PolyCAT A® Column, SS, 200 x 4.6mm, 3µm,1000Å	WCX	Analytical	SS	10-32 coned	200	4.6	3	1000	760
204CT0315	PolyCAT A® Column, SS, 200 x 4.6mm, 3µm,1500Å	WCX	Analytical	SS	10-32 coned	200	4.6	3	1500	760
254CT0303	PolyCAT A® Column, SS, 250 x 4.6mm, 3µm,300Å	WCX	Analytical	SS	10-32 coned	250	4.6	3	300	780
254CT0310	PolyCAT A® Column, SS, 250 x 4.6mm, 3µm,1000Å	WCX	Analytical	SS	10-32 coned	250	4.6	3	1000	780
254CT0315	PolyCAT A® Column, SS, 250 x 4.6mm, 3µm,1500Å	WCX	Analytical	SS	10-32 coned	250	4.6	3	1500	780
109CT0303	PolyCAT A® Column, SS, 100 x 9.4mm, 3µm,300Å	WCX	Analytical	SS	10-32 coned	100	9.4	3	300	1400

PolyCAT A® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore size options.) Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
109CT0310	PolyCAT A® Column, SS, 100 x 9.4mm, 3µm,1000Å	WCX	Analytical	SS	10-32 coned	100	9.4	3	1000	1400
109CT0315	PolyCAT A® Column, SS, 100 x 9.4mm, 3µm,1500Å	WCX	Analytical	SS	10-32 coned	100	9.4	3	1500	1400
209CT0303	PolyCAT A® Column, SS, 200 x 9.4mm, 3µm,300Å	WCX	Analytical	SS	10-32 coned	200	9.4	3	300	1620
209CT0310	PolyCAT A® Column, SS, 200 x 9.4mm, 3µm,1000Å	WCX	Analytical	SS	10-32 coned	200	9.4	3	1000	1620
209CT0315	PolyCAT A® Column, SS, 200 x 9.4mm, 3µm,1500Å	WCX	Analytical	SS	10-32 coned	200	9.4	3	1500	1620
259CT0303	PolyCAT A® Column, SS, 250 x 9.4mm, 3µm,300Å	WCX	Analytical	SS	10-32 coned	250	9.4	3	300	1750
259CT0310	PolyCAT A® Column, SS, 250 x 9.4mm, 3µm,1000Å	WCX	Analytical	SS	10-32 coned	250	9.4	3	1000	1750
259CT0315	PolyCAT A® Column, SS, 250 x 9.4mm, 3µm,1500Å	WCX	Analytical	SS	10-32 coned	250	9.4	3	1500	1750
309CT0303	PolyCAT A® Column, SS, 300 x 9.4mm, 3µm,300Å	WCX	Analytical	SS	10-32 coned	300	9.4	3	300	1850
309CT0310	PolyCAT A® Column, SS, 300 x 9.4mm, 3µm,1000Å	WCX	Analytical	SS	10-32 coned	300	9.4	3	1000	1850
309CT0315	PolyCAT A® Column, SS, 300 x 9.4mm, 3µm,1500Å	WCX	Analytical	SS	10-32 coned	300	9.4	3	1500	1850
051CT0503	PolyCAT A® Column, Glass Lined SS, 50 x 1.0mm, 5µm,300Å	WCX	Analytical	GLT	10-32 coned	50	1.0	5	300	570
051CT0510	PolyCAT A® Column, Glass Lined SS, 50 x 1.0mm, 5µm,1000Å	WCX	Analytical	GLT	10-32 coned	50	1.0	5	1000	570
101CT0503	PolyCAT A® Column, Glass Lined SS, 100 x 1.0mm, 5µm,300Å	WCX	Analytical	GLT	10-32 coned	100	1.0	5	300	660
101CT0510	PolyCAT A® Column, Glass Lined SS, 100 x 1.0mm, 5µm,1000Å	WCX	Analytical	GLT	10-32 coned	100	1.0	5	1000	660
151CT0503	PolyCAT A® Column, Glass Lined SS, 150 x 1.0mm, 5µm,300Å	WCX	Analytical	GLT	10-32 coned	150	1.0	5	300	740
151CT0510	PolyCAT A® Column, Glass Lined SS, 150 x 1.0mm, 5µm,1000Å	WCX	Analytical	GLT	10-32 coned	150	1.0	5	1000	740
3.52CT0503	PolyCAT A® Column, SS, 35 x 2.1mm, 5µm,300Å	WCX	Analytical	SS	10-32 coned	35	2.1	5	300	500
3.52CT0510	PolyCAT A® Column, SS, 35 x 2.1mm, 5µm,1000Å	WCX	Analytical	SS	10-32 coned	35	2.1	5	1000	500
052CT0503	PolyCAT A® Column, SS, 50 x 2.1mm, 5µm,300Å	WCX	Analytical	SS	10-32 coned	50	2.1	5	300	520
052CT0510	PolyCAT A® Column, SS, 50 x 2.1mm, 5µm,1000Å	WCX	Analytical	SS	10-32 coned	50	2.1	5	1000	520
102CT0503	PolyCAT A® Column, SS, 100 x 2.1mm, 5µm,300Å	WCX	Analytical	SS	10-32 coned	100	2.1	5	300	600
102CT0510	PolyCAT A® Column, SS, 100 x 2.1mm, 5µm,1000Å	WCX	Analytical	SS	10-32 coned	100	2.1	5	1000	600
152CT0503	PolyCAT A® Column, SS, 150 x 2.1mm, 5µm,300Å	WCX	Analytical	SS	10-32 coned	150	2.1	5	300	650
152CT0510	PolyCAT A® Column, SS, 150 x 2.1mm, 5µm,1000Å	WCX	Analytical	SS	10-32 coned	150	2.1	5	1000	650
202CT0503	PolyCAT A® Column, SS, 200 x 2.1mm, 5µm,300Å	WCX	Analytical	SS	10-32 coned	200	2.1	5	300	710
202CT0510	PolyCAT A® Column, SS, 200 x 2.1mm, 5µm,1000Å	WCX	Analytical	SS	10-32 coned	200	2.1	5	1000	710
3.54CT0503	PolyCAT A® Column, SS, 35 x 4.6mm, 5µm,300Å	WCX	Analytical	SS	10-32 coned	35	4.6	5	300	500
3.54CT0510	PolyCAT A® Column, SS, 35 x 4.6mm, 5µm,1000Å	WCX	Analytical	SS	10-32 coned	35	4.6	5	1000	500
054CT0503	PolyCAT A® Column, SS, 50 x 4.6mm, 5µm,300Å	WCX	Analytical	SS	10-32 coned	50	4.6	5	300	520
054CT0510	PolyCAT A® Column, SS, 50 x 4.6mm, 5µm,1000Å	WCX	Analytical	SS	10-32 coned	50	4.6	5	1000	520
104CT0503	PolyCAT A® Column, SS, 100 x 4.6mm, 5µm,300Å	WCX	Analytical	SS	10-32 coned	100	4.6	5	300	600
104CT0510	PolyCAT A® Column, SS, 100 x 4.6mm, 5µm,1000Å	WCX	Analytical	SS	10-32 coned	100	4.6	5	1000	600
154CT0503	PolyCAT A® Column, SS, 150 x 4.6mm, 5µm,300Å	WCX	Analytical	SS	10-32 coned	150	4.6	5	300	650
154CT0510	PolyCAT A® Column, SS, 150 x 4.6mm, 5µm,1000Å	WCX	Analytical	SS	10-32 coned	150	4.6	5	1000	650
204CT0503	PolyCAT A® Column, SS, 200 x 4.6mm, 5µm,300Å	WCX	Analytical	SS	10-32 coned	200	4.6	5	300	710
204CT0510	PolyCAT A® Column, SS, 200 x 4.6mm, 5µm,1000Å	WCX	Analytical	SS	10-32 coned	200	4.6	5	1000	710
254CT0503	PolyCAT A® Column, SS, 250 x 4.6mm, 5µm,300Å	WCX	Analytical	SS	10-32 coned	250	4.6	5	300	740
254CT0510	PolyCAT A® Column, SS, 250 x 4.6mm, 5µm,1000Å	WCX	Analytical	SS	10-32 coned	250	4.6	5	1000	740
109CT0503	PolyCAT A® Column, SS, 100 x 9.4mm, 5µm,300Å	WCX	Semi-Prep	SS	10-32 coned	100	9.4	5	300	1110
109CT0510	PolyCAT A® Column, SS, 100 x 9.4mm, 5µm,1000Å	WCX	Semi-Prep	SS	10-32 coned	100	9.4	5	1000	1110
209CT0503	PolyCAT A® Column, SS, 200 x 9.4mm, 5µm,300Å	WCX	Semi-Prep	SS	10-32 coned	200	9.4	5	300	1400
209CT0510	PolyCAT A® Column, SS, 200 x 9.4mm, 5µm,1000Å	WCX	Semi-Prep	SS	10-32 coned	200	9.4	5	1000	1400
259CT0503	PolyCAT A® Column, SS, 250 x 9.4mm, 5µm,300Å	WCX	Semi-Prep	SS	10-32 coned	250	9.4	5	300	1500
259CT0510	PolyCAT A® Column, SS, 250 x 9.4mm, 5µm,1000Å	WCX	Semi-Prep	SS	10-32 coned	250	9.4	5	1000	1500

PolyCAT A® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore size options. Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
309CT0503	PolyCAT A® Column, SS, 300 x 9.4mm, 5µm,300Å	WCX	Semi-Prep	SS	10-32 coned	300	9.4	5	300	1600
309CT0510	PolyCAT A® Column, SS, 300 x 9.4mm, 5µm,1000Å	WCX	Semi-Prep	SS	10-32 coned	300	9.4	5	1000	1600
2521CT0503	PolyCAT A® Column, SS, 250 x 21mm, 5µm,300Å	WCX	Semi-Prep	SS	10-32 coned	250	21.0	5	300	4500
2521CT0510	PolyCAT A® Column, SS, 250 x 21mm, 5µm,1000Å	WCX	Semi-Prep	SS	10-32 coned	250	21.0	5	1000	4500
204CT1203	PolyCAT A® Column, SS, 200 x 4.6mm, 12µm,300Å	WCX	Analytical	SS	10-32 coned	200	4.6	12	300	420
204CT1215	PolyCAT A® Column, SS, 200 x 4.6mm, 12µm,1500Å	WCX	Analytical	SS	10-32 coned	200	4.6	12	1500	420
254CT1203	PolyCAT A® Column, SS, 250 x 4.6mm, 12µm,300Å	WCX	Analytical	SS	10-32 coned	250	4.6	12	300	500
254CT1215	PolyCAT A® Column, SS, 250 x 4.6mm, 12µm,1500Å	WCX	Analytical	SS	10-32 coned	250	4.6	12	1500	500
209CT1203	PolyCAT A® Column, SS, 200 x 9.4mm, 12µm,300Å	WCX	Semi-Prep	SS	10-32 coned	200	9.4	12	300	1000
209CT1215	PolyCAT A® Column, SS, 200 x 9.4mm, 12µm,1500Å	WCX	Semi-Prep	SS	10-32 coned	200	9.4	12	1500	1000
259CT1203	PolyCAT A® Column, SS, 250 x 9.4mm, 12µm,300Å	WCX	Semi-Prep	SS	10-32 coned	250	9.4	12	300	1100
259CT1215	PolyCAT A® Column, SS, 250 x 9.4mm, 12µm,1500Å	WCX	Semi-Prep	SS	10-32 coned	250	9.4	12	1500	1100
309CT1203	PolyCAT A® Column, SS, 300 x 9.4mm, 12µm,300Å	WCX	Semi-Prep	SS	10-32 coned	300	9.4	12	300	1300
309CT1215	PolyCAT A® Column, SS, 300 x 9.4mm, 12µm,1500Å	WCX	Semi-Prep	SS	10-32 coned	300	9.4	12	1500	1300
2521CT1203	PolyCAT A® Column, SS, 250 x 21mm, 12µm,300Å	WCX	Semi-Prep	SS	10-32 coned	250	21.0	12	300	3200
2521CT1215	PolyCAT A® Column, SS, 250 x 21mm, 12µm,1500Å	WCX	Semi-Prep	SS	10-32 coned	250	21.0	12	1500	3200
2550CT1203	PolyCAT A® Column, SS, 250 x 50.8mm, 12µm,300Å	WCX	Prep	SS	10-32 coned	250	50.8	12	300	12250
2550CT1215	PolyCAT A® Column, SS, 250 x 50.8mm, 12µm,1500Å	WCX	Prep	SS	10-32 coned	250	50.8	12	1500	12250



HPLC Capillaries

PolyCAT A® Capillaries (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore size options.) Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
050.075CT0203	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 2µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.075	2	300	520
050.075CT0210	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 2µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.075	2	1000	520
050.075CT0215	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 2µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.075	2	1500	520
100.075CT0203	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 2µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.075	2	300	580
100.075CT0210	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 2µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.075	2	1000	580
100.075CT0215	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 2µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.075	2	1500	580
150.075CT0203	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 2µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.075	2	300	640
150.075CT0210	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 2µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.075	2	1000	640
150.075CT0215	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 2µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.075	2	1500	640
200.075CT0203	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 2µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.075	2	300	700
200.075CT0210	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 2µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.075	2	1000	700
200.075CT0215	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 2µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.075	2	1500	700
050.15CT0203	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 2µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.150	2	300	520
050.15CT0210	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 2µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.150	2	1000	520
050.15CT0215	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 2µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.150	2	1500	520

PolyCAT A® Capillaries (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore size options.) Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
100.15CT0203	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 2µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.150	2	300	580
100.15CT0210	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 2µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.150	2	1000	580
100.15CT0215	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 2µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.150	2	1500	580
150.15CT0203	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 2µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.150	2	300	640
150.15CT0210	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 2µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.150	2	1000	640
150.15CT0215	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 2µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.150	2	1500	640
200.15CT0203	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 2µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.150	2	300	700
200.15CT0210	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 2µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.150	2	1000	700
200.15CT0215	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 2µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.150	2	1500	700
050.30CT0203	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 2µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.300	2	300	520
050.30CT0210	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 2µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.300	2	1000	520
050.30CT0215	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 2µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.300	2	1500	520
100.30CT0203	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 2µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.300	2	300	580
100.30CT0210	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 2µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.300	2	1000	580
100.30CT0215	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 2µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.300	2	1500	580
150.30CT0203	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 2µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.300	2	300	640
150.30CT0210	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 2µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.300	2	1000	640
150.30CT0215	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 2µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.300	2	1500	640
200.30CT0203	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 2µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.300	2	300	700
200.30CT0210	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 2µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.300	2	1000	700
200.30CT0215	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 2µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.300	2	1500	700
050.075CT0303	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.075	3	300	520
050.075CT0310	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.075	3	1000	520
050.075CT0315	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.075	3	1500	520
100.075CT0303	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.075	3	300	580
100.075CT0310	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.075	3	1000	580
100.075CT0315	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.075	3	1500	580
150.075CT0303	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.075	3	300	640
150.075CT0310	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.075	3	1000	640
150.075CT0315	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.075	3	1500	640
200.075CT0303	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.075	3	300	700
200.075CT0310	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.075	3	1000	700
200.075CT0315	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.075	3	1500	700
050.15CT0303	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.150	3	300	520
050.15CT0310	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.150	3	1000	520
050.15CT0315	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.150	3	1500	520
100.15CT0303	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.150	3	300	580
100.15CT0310	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.150	3	1000	580
100.15CT0315	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.150	3	1500	580
150.15CT0303	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.150	3	300	640
150.15CT0310	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.150	3	1000	640
150.15CT0315	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.150	3	1500	640
200.15CT0303	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.150	3	300	700
200.15CT0310	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.150	3	1000	700
200.15CT0315	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.150	3	1500	700

PolyCAT A® Capillaries (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore size options.) Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
050.30CT0303	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.300	3	300	520
050.30CT0310	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.300	3	1000	520
050.30CT0315	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.300	3	1500	520
100.30CT0303	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.300	3	300	580
100.30CT0310	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.300	3	1000	580
100.30CT0315	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.300	3	1500	580
150.30CT0303	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.300	3	300	640
150.30CT0310	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.300	3	1000	640
150.30CT0315	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.300	3	1500	640
200.30CT0303	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.300	3	300	700
200.30CT0310	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.300	3	1000	700
200.30CT0315	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,1500Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.300	3	1500	700
050.075CT0503	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 5µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.075	5	300	520
050.075CT0510	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 5µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.075	5	1000	520
100.075CT0503	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 5µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.075	5	300	580
100.075CT0510	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 5µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.075	5	1000	580
150.075CT0503	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 5µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.075	5	300	640
150.075CT0510	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 5µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.075	5	1000	640
200.075CT0503	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 5µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.075	5	300	700
200.075CT0510	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 5µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.075	5	1000	700
050.15CT0503	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 5µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.150	5	300	520
050.15CT0510	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 5µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.150	5	1000	520
100.15CT0503	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 5µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.150	5	300	580
100.15CT0510	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 5µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.150	5	1000	580
150.15CT0503	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 5µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.150	5	300	640
150.15CT0510	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 5µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.150	5	1000	640
200.15CT0503	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 5µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.150	5	300	700
200.15CT0510	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 5µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.150	5	1000	700
050.30CT0503	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 5µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.300	5	300	520
050.30CT0510	PolyCAT A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 5µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	50	0.300	5	1000	520
100.30CT0503	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 5µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.300	5	300	580
100.30CT0510	PolyCAT A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 5µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	100	0.300	5	1000	580
150.30CT0503	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 5µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.300	5	300	640
150.30CT0510	PolyCAT A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 5µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	150	0.300	5	1000	640
200.30CT0503	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 5µm,300Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.300	5	300	700
200.30CT0510	PolyCAT A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 5µm,1000Å	WCX	Analytical	PEEK SIL	10-32 coned	200	0.300	5	1000	700

Guard Cartridges



PolyCAT A° Guard Cartridges (Available in 2, 3, 5, 12, 20, and 30µm particles with pore size options.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
GC1-CT0203	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 1mm, 2µm, 300Å	WCX	Analytical	SS/PEEK	10-32 coned	10	1.0	2	300	120
GC1-CT0210	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 1mm, 2µm, 1000Å	WCX	Analytical	SS/PEEK	10-32 coned	10	1.0	2	1000	120
GC1-CT0215	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 1mm, 2µm, 1500Å	WCX	Analytical	SS/PEEK	10-32 coned	10	1.0	2	1500	120
GC2-CT0203	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 2µm, 300Å	WCX	Analytical	SS/PEEK	10-32 coned	10	2.0	2	300	120
GC2-CT0210	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 2µm, 1000Å	WCX	Analytical	SS/PEEK	10-32 coned	10	2.0	2	1000	120
GC2-CT0215	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 2µm, 1500Å	WCX	Analytical	SS/PEEK	10-32 coned	10	2.0	2	1500	120
GC3-CT0203	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 2µm, 300Å	WCX	Analytical	SS/PEEK	10-32 coned	10	3.0	2	300	120
GC3-CT0210	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 2µm, 1000Å	WCX	Analytical	SS/PEEK	10-32 coned	10	3.0	2	1000	120
GC3-CT0215	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 2µm, 1500Å	WCX	Analytical	SS/PEEK	10-32 coned	10	3.0	2	1500	120
GC4-CT0203	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 2µm, 300Å	WCX	Analytical	SS/PEEK	10-32 coned	10	4.0	2	300	120
GC4-CT0210	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 2µm, 1000Å	WCX	Analytical	SS/PEEK	10-32 coned	10	4.0	2	1000	120
GC4-CT0215	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 2µm, 1500Å	WCX	Analytical	SS/PEEK	10-32 coned	10	4.0	2	1500	120
GC1-CT0303	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 300Å	WCX	Analytical	SS/PEEK	10-32 coned	10	1.0	3	300	100
GC1-CT0310	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 1000Å	WCX	Analytical	SS/PEEK	10-32 coned	10	1.0	3	1000	100
GC1-CT0315	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 1500Å	WCX	Analytical	SS/PEEK	10-32 coned	10	1.0	3	1500	100
GC2-CT0303	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 3µm, 300Å	WCX	Analytical	SS/PEEK	10-32 coned	10	2.0	3	300	100
GC2-CT0310	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 3µm, 1000Å	WCX	Analytical	SS/PEEK	10-32 coned	10	2.0	3	1000	100
GC2-CT0315	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 3µm, 1500Å	WCX	Analytical	SS/PEEK	10-32 coned	10	2.0	3	1500	100
GC3-CT0303	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 3µm, 300Å	WCX	Analytical	SS/PEEK	10-32 coned	10	3.0	3	300	100
GC3-CT0310	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 3µm, 1000Å	WCX	Analytical	SS/PEEK	10-32 coned	10	3.0	3	1000	100
GC3-CT0315	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 3µm, 1500Å	WCX	Analytical	SS/PEEK	10-32 coned	10	3.0	3	1500	100
GC4-CT0303	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 3µm, 300Å	WCX	Analytical	SS/PEEK	10-32 coned	10	4.0	3	300	100
GC4-CT0310	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 3µm, 1000Å	WCX	Analytical	SS/PEEK	10-32 coned	10	4.0	3	1000	100
GC4-CT0315	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 3µm, 1500Å	WCX	Analytical	SS/PEEK	10-32 coned	10	4.0	3	1500	100
GC1-CT0503	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 1mm, 5µm, 300Å	WCX	Analytical	SS/PEEK	10-32 coned	10	1.0	5	300	90
GC1-CT0510	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 1mm, 5µm, 1000Å	WCX	Analytical	SS/PEEK	10-32 coned	10	1.0	5	1000	90
GC2-CT0503	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 5µm, 300Å	WCX	Analytical	SS/PEEK	10-32 coned	10	2.0	5	300	90
GC2-CT0510	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 5µm, 1000Å	WCX	Analytical	SS/PEEK	10-32 coned	10	2.0	5	1000	90
GC3-CT0503	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 5µm, 300Å	WCX	Analytical	SS/PEEK	10-32 coned	10	3.0	5	300	90
GC3-CT0510	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 5µm, 1000Å	WCX	Analytical	SS/PEEK	10-32 coned	10	3.0	5	1000	90
GC4-CT0503	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 5µm, 300Å	WCX	Analytical	SS/PEEK	10-32 coned	10	4.0	5	300	90
GC4-CT0510	PolyCAT A° Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 5µm, 1000Å	WCX	Analytical	SS/PEEK	10-32 coned	10	4.0	5	1000	90
DGCH10	Direct Guard Cartridge Holder, 10mm (For all analytical guards, REUSABLE)	N/A	Analytical	SS/PEEK	10-32 coned	10	N/A	N/A	N/A	75
IGCH10	Inline Guard Cartridge Holder, 10mm (For all analytical guards, REUSABLE)	N/A	Analytical	SS/PEEK	10-32 coned	10	N/A	N/A	N/A	75

*** NOTE: DGCH10 for Direct Connection to Column (female inlet/male outlet) IGCH10 for Connection to Line (two female ports)

GC10-CT0303	PolyCAT A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 300Å	WCX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	300	225
GC10-CT0310	PolyCAT A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 1000Å	WCX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	1000	225
GC10-CT0315	PolyCAT A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 1500Å	WCX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	1500	225
GC10-CT0503	PolyCAT A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 5µm, 300Å	WCX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	5	300	185
GC10-CT0510	PolyCAT A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 5µm, 1000Å	WCX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	5	1000	185
CH10	Semi-Prep Guard Cartridge Holder, 10mm (For 9.4mm columns, REUSABLE)	WCX	Semi-Prep	SS/PEEK	10-32 coned	10	N/A	N/A	N/A	175

*** NOTE: CH10 Holder Required for Semi-Prep Guard Cartridges (10 x 10mm) Connection to Line (two female ports)

Solid Phase Extraction Cartridges

PolyCAT A® SPE Cartridges

Custom combinations available upon request info@polylc.c

PN	Product Description	Phase	Type	Hardware	Connection	V (ml)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
SPECT1203	PolyCAT A® SPE Cartridge 1mL, PE, 12µm, 300Å (Pack of 10 cartridges)	WCX	SPE	PE	Lure	1	N/A	12	300	130
SPECT1215	PolyCAT A® SPE Cartridge 1mL, PE, 12µm, 1500Å (Pack of 10 cartridges)	WCX	SPE	PE	Lure	1	N/A	12	1500	130
SPECT2003	PolyCAT A® SPE Cartridge 1mL, PE, 20µm, 300Å (Pack of 10 cartridges)	WCX	SPE	PE	Lure	1	N/A	20	300	130
SPECT2015	PolyCAT A® SPE Cartridge 1mL, PE, 20µm, 1500Å (Pack of 10 cartridges)	WCX	SPE	PE	Lure	1	N/A	20	1500	130

Spin Columns

PolyCAT A® BioPureSPN Columns

Custom combinations available upon request info@polylc.c

PN	Product Description	Phase	Type	Bed (mg)	Pore (Å)	WCX Capacity	HILIC Capacity	2026 Price (USD)
HMM HIL-WCX	PolyCAT A® BioPureSPN MACRO for WCX or HILIC 50-200µl (Box of 100 each)	WCX	Spin	100	300	0.2 -1mg	35 - 350µg	940
HMM HIL-WCX.20	PolyCAT A® BioPureSPN MACRO for WCX or HILIC 50-200µl (Pack of 20 each)	WCX	Spin	100	300	0.2 -1mg	35 - 350µg	265
HMM HIL-WCX-L	PolyCAT A® BioPureSPN MACRO PLATE 96 Well for WCX or HILIC 50-200µl	WCX	Spin	100	300	0.2 -1mg	35 - 350µg	800
HEM HIL-WCX	PolyCAT A® BioPureSPN MIDI for WCX/HILIC 50-200µl (Box of 100 each)	WCX	Spin	50	300	100-500µg	17 - 170µg	820
HEM HIL-WCX.20	PolyCAT A® BioPureSPN MIDI for WCX/HILIC 50-200µl (Pack of 20 each)	WCX	Spin	50	300	100-500µg	17 - 170µg	240
HNS HIL-WCX-M	PolyCAT A® BioPureSPN MIDI PLATE 96 Well for WCX/HILIC 50-200µl	WCX	Spin	50	300	100-500µg	17 - 170µg	575
HUM HIL-WCX	PolyCAT A® BioPureSPN MINI for WCX/HILIC 20-100µl (Box of 100 each)	WCX	Spin	20	300	40-200µg	7 - 70µg	735
HUM HIL-WCX.20	PolyCAT A® BioPureSPN MINI for WCX/HILIC 20-100µl (Pack of 20 each)	WCX	Spin	20	300	40-200µg	7 - 70µg	225
HNS HIL-WCX-M	PolyCAT A® BioPureSPN MINI PLATE 96 Well for WCX/HILIC 20-100µl	WCX	Spin	20	300	40-200µg	7 - 70µg	735

Bulk Material



PolyCAT A® Bulk Material (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore size options.) Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	QTY (g)	Particle (µm)	Pore (Å)	2026 Price (USD)
BMCT0203	PolyCAT A® Bulk Material, 2µm, 300Å	WCX	BM	N/A	0.5	2	300	100
BMCT0210	PolyCAT A® Bulk Material, 2µm, 1000Å	WCX	BM	N/A	0.5	2	1000	100
BMCT0215	PolyCAT A® Bulk Material, 2µm, 1500Å	WCX	BM	N/A	0.5	2	1500	100
BMCT0303	PolyCAT A® Bulk Material, 3µm, 300Å	WCX	BM	N/A	1	3	300	130
BMCT0310	PolyCAT A® Bulk Material, 3µm, 1000Å	WCX	BM	N/A	1	3	1000	130
BMCT0315	PolyCAT A® Bulk Material, 3µm, 1500Å	WCX	BM	N/A	1	3	1500	130
BMCT0503	PolyCAT A® Bulk Material, 5µm, 300Å	WCX	BM	N/A	1	5	300	85
BMCT0510	PolyCAT A® Bulk Material, 5µm, 1000Å	WCX	BM	N/A	1	5	1000	85

PolyCAT A® Bulk Material (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore size options.) Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	QTY (g)	Particle (µm)	Pore (Å)	2026 Price (USD)
BMCT1203	PolyCAT A® Bulk Material, 12µm, 300Å	WCX	BM	N/A	1	12	300	32
BMCT1215	PolyCAT A® Bulk Material, 12µm, 1000Å	WCX	BM	N/A	1	12	1000	32
BMCT2003	PolyCAT A® Bulk Material, 20µm, 300Å	WCX	BM	N/A	1	20	300	32
BMCT2015	PolyCAT A® Bulk Material, 20µm, 1500Å	WCX	BM	N/A	1	20	1500	32
BMCT3003	PolyCAT A® Bulk Material, 30µm, 300Å	WCX	BM	N/A	1	30	300	32
BMCT3010	PolyCAT A® Bulk Material, 30µm, 1000Å	WCX	BM	N/A	1	30	1000	32

Please contact info@polylc.com with any questions or for information on any product combination not listed in the table above

PolyCAT A®

(Weak-Cation Exchange)

2026 PolyLC Inc. Retail Pricelist

PolyLC Inc., 9151 Rumsey Rd., Suite 180, Columbia, MD 21045 410-992-5400 www.polylc.com



2026 PolyLC Inc. Retail Pricelist

PolyHYDROXYETHYL A®

(HILIC and SEC Separations)



Offering high-performance liquid chromatography media, columns, capillaries, cartridges, and instrumentation. We specialize in advanced alternatives to reversed-phase for proteins, peptides, and other biomolecules, supported by deep application expertise and customized solutions for the most demanding separations.

HPLC Columns

PolyHYDROXYETHYL A® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore size options.)

Custom combinations available: info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
051HY0201	PolyHYDROXYETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 2µm,100Å	HILIC / SEC	Analytical	GLT	10-32 coned	50	1.0	2	100	615
051HY0203	PolyHYDROXYETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 2µm,300Å	HILIC / SEC	Analytical	GLT	10-32 coned	50	1.0	2	300	615
051HY0210	PolyHYDROXYETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 2µm,1000Å	HILIC / SEC	Analytical	GLT	10-32 coned	50	1.0	2	1000	615
101HY0201	PolyHYDROXYETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 2µm,100Å	HILIC / SEC	Analytical	GLT	10-32 coned	100	1.0	2	100	750
101HY0203	PolyHYDROXYETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 2µm,300Å	HILIC / SEC	Analytical	GLT	10-32 coned	100	1.0	2	300	750
101HY0210	PolyHYDROXYETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 2µm,1000Å	HILIC / SEC	Analytical	GLT	10-32 coned	100	1.0	2	1000	750
151HY0201	PolyHYDROXYETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 2µm,100Å	HILIC / SEC	Analytical	GLT	10-32 coned	150	1.0	2	100	800
151HY0203	PolyHYDROXYETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 2µm,300Å	HILIC / SEC	Analytical	GLT	10-32 coned	150	1.0	2	300	800
151HY0210	PolyHYDROXYETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 2µm,1000Å	HILIC / SEC	Analytical	GLT	10-32 coned	150	1.0	2	1000	800
3.52HY0201	PolyHYDROXYETHYL A® Column, SS, 35 x 2.1mm, 2µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	35	2.1	2	100	570
3.52HY0203	PolyHYDROXYETHYL A® Column, SS, 35 x 2.1mm, 2µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	35	2.1	2	300	570
3.52HY0210	PolyHYDROXYETHYL A® Column, SS, 35 x 2.1mm, 2µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	35	2.1	2	1000	570
052HY0201	PolyHYDROXYETHYL A® Column, SS, 50 x 2.1mm, 2µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	50	2.1	2	100	580
052HY0203	PolyHYDROXYETHYL A® Column, SS, 50 x 2.1mm, 2µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	50	2.1	2	300	580
052HY0210	PolyHYDROXYETHYL A® Column, SS, 50 x 2.1mm, 2µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	50	2.1	2	1000	580
102HY0201	PolyHYDROXYETHYL A® Column, SS, 100 x 2.1mm, 2µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	100	2.1	2	100	750
102HY0203	PolyHYDROXYETHYL A® Column, SS, 100 x 2.1mm, 2µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	100	2.1	2	300	750
102HY0210	PolyHYDROXYETHYL A® Column, SS, 100 x 2.1mm, 2µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	100	2.1	2	1000	750
3.54HY0201	PolyHYDROXYETHYL A® Column, SS, 35 x 4.6mm, 2µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	35	4.6	2	100	570
3.54HY0203	PolyHYDROXYETHYL A® Column, SS, 35 x 4.6mm, 2µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	35	4.6	2	300	570
3.54HY0210	PolyHYDROXYETHYL A® Column, SS, 35 x 4.6mm, 2µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	35	4.6	2	1000	570
054HY0201	PolyHYDROXYETHYL A® Column, SS, 50 x 4.6mm, 2µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	50	4.6	2	100	580
054HY0203	PolyHYDROXYETHYL A® Column, SS, 50 x 4.6mm, 2µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	50	4.6	2	300	580
054HY0210	PolyHYDROXYETHYL A® Column, SS, 50 x 4.6mm, 2µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	50	4.6	2	1000	580
104HY0201	PolyHYDROXYETHYL A® Column, SS, 100 x 4.6mm, 2µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	100	4.6	2	100	750
104HY0203	PolyHYDROXYETHYL A® Column, SS, 100 x 4.6mm, 2µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	100	4.6	2	300	750
104HY0210	PolyHYDROXYETHYL A® Column, SS, 100 x 4.6mm, 2µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	100	4.6	2	1000	750
154HY0201	PolyHYDROXYETHYL A® Column, SS, 150 x 4.6mm, 2µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	150	4.6	2	100	790
154HY0203	PolyHYDROXYETHYL A® Column, SS, 150 x 4.6mm, 2µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	150	4.6	2	300	790
154HY0210	PolyHYDROXYETHYL A® Column, SS, 150 x 4.6mm, 2µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	150	4.6	2	1000	790
204HY0201	PolyHYDROXYETHYL A® Column, SS, 200 x 4.6mm, 2µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	200	4.6	2	100	810
204HY0203	PolyHYDROXYETHYL A® Column, SS, 200 x 4.6mm, 2µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	200	4.6	2	300	810
204HY0210	PolyHYDROXYETHYL A® Column, SS, 200 x 4.6mm, 2µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	200	4.6	2	1000	810

PolyHYDROXYETHYL A® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore size options.) Custom combinations available: info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
109HY0201	PolyHYDROXYETHYL A® Column, SS, 100 x 9.4mm, 2µm,100Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	100	9.4	2	100	1525
109HY0203	PolyHYDROXYETHYL A® Column, SS, 100 x 9.4mm, 2µm,300Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	100	9.4	2	300	1525
109HY0210	PolyHYDROXYETHYL A® Column, SS, 100 x 9.4mm, 2µm,1000Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	100	9.4	2	1000	1525
209HY0201	PolyHYDROXYETHYL A® Column, SS, 200 x 9.4mm, 2µm,100Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	200	9.4	2	100	2450
209HY0203	PolyHYDROXYETHYL A® Column, SS, 200 x 9.4mm, 2µm,300Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	200	9.4	2	300	2450
209HY0210	PolyHYDROXYETHYL A® Column, SS, 200 x 9.4mm, 2µm,1000Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	200	9.4	2	1000	2450
051HY03006	PolyHYDROXYETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 3µm,60Å	HILIC / SEC	Analytical	GLT	10-32 coned	50	1.0	3	60	570
051HY0301	PolyHYDROXYETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 3µm,100Å	HILIC / SEC	Analytical	GLT	10-32 coned	50	1.0	3	100	570
051HY0303	PolyHYDROXYETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 3µm,300Å	HILIC / SEC	Analytical	GLT	10-32 coned	50	1.0	3	300	570
051HY0305	PolyHYDROXYETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 3µm,500Å	HILIC / SEC	Analytical	GLT	10-32 coned	50	1.0	3	500	570
051HY0310	PolyHYDROXYETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 3µm,1000Å	HILIC / SEC	Analytical	GLT	10-32 coned	50	1.0	3	1000	570
051HY0315	PolyHYDROXYETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 3µm,1500Å	HILIC / SEC	Analytical	GLT	10-32 coned	50	1.0	3	1500	570
101HY03006	PolyHYDROXYETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 3µm,60Å	HILIC / SEC	Analytical	GLT	10-32 coned	100	1.0	3	60	660
101HY0301	PolyHYDROXYETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 3µm,100Å	HILIC / SEC	Analytical	GLT	10-32 coned	100	1.0	3	100	660
101HY0303	PolyHYDROXYETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 3µm,300Å	HILIC / SEC	Analytical	GLT	10-32 coned	100	1.0	3	300	660
101HY0305	PolyHYDROXYETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 3µm,500Å	HILIC / SEC	Analytical	GLT	10-32 coned	100	1.0	3	500	660
101HY0310	PolyHYDROXYETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 3µm,1000Å	HILIC / SEC	Analytical	GLT	10-32 coned	100	1.0	3	1000	660
101HY0315	PolyHYDROXYETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 3µm,1500Å	HILIC / SEC	Analytical	GLT	10-32 coned	100	1.0	3	1500	660
151HY03006	PolyHYDROXYETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 3µm,60Å	HILIC / SEC	Analytical	GLT	10-32 coned	150	1.0	3	60	740
151HY0301	PolyHYDROXYETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 3µm,100Å	HILIC / SEC	Analytical	GLT	10-32 coned	150	1.0	3	100	740
151HY0303	PolyHYDROXYETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 3µm,300Å	HILIC / SEC	Analytical	GLT	10-32 coned	150	1.0	3	300	740
151HY0305	PolyHYDROXYETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 3µm,500Å	HILIC / SEC	Analytical	GLT	10-32 coned	150	1.0	3	500	740
151HY0310	PolyHYDROXYETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 3µm,1000Å	HILIC / SEC	Analytical	GLT	10-32 coned	150	1.0	3	1000	740
151HY0315	PolyHYDROXYETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 3µm,1500Å	HILIC / SEC	Analytical	GLT	10-32 coned	150	1.0	3	1500	740
3.52HY03006	PolyHYDROXYETHYL A® Column, SS, 35 x 2.1mm, 3µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	35	2.1	3	60	550
3.52HY0301	PolyHYDROXYETHYL A® Column, SS, 35 x 2.1mm, 3µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	35	2.1	3	100	550
3.52HY0303	PolyHYDROXYETHYL A® Column, SS, 35 x 2.1mm, 3µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	35	2.1	3	300	550
3.52HY0305	PolyHYDROXYETHYL A® Column, SS, 35 x 2.1mm, 3µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	50	2.1	3	500	550
3.52HY0310	PolyHYDROXYETHYL A® Column, SS, 35 x 2.1mm, 3µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	50	2.1	3	1000	550
3.52HY0315	PolyHYDROXYETHYL A® Column, SS, 35 x 2.1mm, 3µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	50	2.1	3	1500	550
052HY03006	PolyHYDROXYETHYL A® Column, SS, 50 x 2.1mm, 3µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	100	2.1	3	60	560
052HY0301	PolyHYDROXYETHYL A® Column, SS, 50 x 2.1mm, 3µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	100	2.1	3	100	560
052HY0303	PolyHYDROXYETHYL A® Column, SS, 50 x 2.1mm, 3µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	100	2.1	3	300	560
052HY0305	PolyHYDROXYETHYL A® Column, SS, 50 x 2.1mm, 3µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	35	2.1	3	500	560
052HY0310	PolyHYDROXYETHYL A® Column, SS, 50 x 2.1mm, 3µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	35	2.1	3	1000	560
052HY0315	PolyHYDROXYETHYL A® Column, SS, 50 x 2.1mm, 3µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	35	2.1	3	1500	560
102HY03006	PolyHYDROXYETHYL A® Column, SS, 100 x 2.1mm, 3µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	50	2.1	3	60	710
102HY0301	PolyHYDROXYETHYL A® Column, SS, 100 x 2.1mm, 3µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	50	2.1	3	100	710
102HY0303	PolyHYDROXYETHYL A® Column, SS, 100 x 2.1mm, 3µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	50	2.1	3	300	710
102HY0305	PolyHYDROXYETHYL A® Column, SS, 100 x 2.1mm, 3µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	100	2.1	3	500	710
102HY0310	PolyHYDROXYETHYL A® Column, SS, 100 x 2.1mm, 3µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	100	2.1	3	1000	710
102HY0315	PolyHYDROXYETHYL A® Column, SS, 100 x 2.1mm, 3µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	100	2.1	3	1500	710
3.54HY03006	PolyHYDROXYETHYL A® Column, SS, 35 x 4.6mm, 3µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	35	4.6	3	60	550
3.54HY0301	PolyHYDROXYETHYL A® Column, SS, 35 x 4.6mm, 3µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	35	4.6	3	100	550
3.54HY0303	PolyHYDROXYETHYL A® Column, SS, 35 x 4.6mm, 3µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	35	4.6	3	300	550
3.54HY0305	PolyHYDROXYETHYL A® Column, SS, 35 x 4.6mm, 3µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	35	4.6	3	500	550
3.54HY0310	PolyHYDROXYETHYL A® Column, SS, 35 x 4.6mm, 3µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	35	4.6	3	1000	550

PolyHYDROXYETHYL A® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore size options.) Custom combinations available: info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
3.54HY0315	PolyHYDROXYETHYL A® Column, SS, 35 x 4.6mm, 3µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	35	4.6	3	1500	550
054HY03006	PolyHYDROXYETHYL A® Column, SS, 50 x 4.6mm, 3µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	50	4.6	3	60	560
054HY0301	PolyHYDROXYETHYL A® Column, SS, 50 x 4.6mm, 3µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	50	4.6	3	100	560
054HY0303	PolyHYDROXYETHYL A® Column, SS, 50 x 4.6mm, 3µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	50	4.6	3	300	560
054HY0305	PolyHYDROXYETHYL A® Column, SS, 50 x 4.6mm, 3µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	50	4.6	3	500	560
054HY0310	PolyHYDROXYETHYL A® Column, SS, 50 x 4.6mm, 3µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	50	4.6	3	1000	560
054HY0315	PolyHYDROXYETHYL A® Column, SS, 50 x 4.6mm, 3µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	50	4.6	3	1500	560
104HY03006	PolyHYDROXYETHYL A® Column, SS, 100 x 4.6mm, 3µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	100	4.6	3	60	710
104HY0301	PolyHYDROXYETHYL A® Column, SS, 100 x 4.6mm, 3µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	100	4.6	3	100	710
104HY0303	PolyHYDROXYETHYL A® Column, SS, 100 x 4.6mm, 3µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	100	4.6	3	300	710
104HY0305	PolyHYDROXYETHYL A® Column, SS, 100 x 4.6mm, 3µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	100	4.6	3	500	710
104HY0310	PolyHYDROXYETHYL A® Column, SS, 100 x 4.6mm, 3µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	100	4.6	3	1000	710
104HY0315	PolyHYDROXYETHYL A® Column, SS, 100 x 4.6mm, 3µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	100	4.6	3	1500	710
154HY03006	PolyHYDROXYETHYL A® Column, SS, 150 x 4.6mm, 3µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	150	4.6	3	60	740
154HY0301	PolyHYDROXYETHYL A® Column, SS, 150 x 4.6mm, 3µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	150	4.6	3	100	740
154HY0303	PolyHYDROXYETHYL A® Column, SS, 150 x 4.6mm, 3µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	150	4.6	3	300	740
154HY0305	PolyHYDROXYETHYL A® Column, SS, 150 x 4.6mm, 3µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	150	4.6	3	500	740
154HY0310	PolyHYDROXYETHYL A® Column, SS, 150 x 4.6mm, 3µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	150	4.6	3	1000	740
154HY0315	PolyHYDROXYETHYL A® Column, SS, 150 x 4.6mm, 3µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	150	4.6	3	1500	740
204HY03006	PolyHYDROXYETHYL A® Column, SS, 200 x 4.6mm, 3µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	200	4.6	3	60	760
204HY0301	PolyHYDROXYETHYL A® Column, SS, 200 x 4.6mm, 3µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	200	4.6	3	100	760
204HY0303	PolyHYDROXYETHYL A® Column, SS, 200 x 4.6mm, 3µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	200	4.6	3	300	760
204HY0305	PolyHYDROXYETHYL A® Column, SS, 200 x 4.6mm, 3µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	200	4.6	3	500	760
204HY0310	PolyHYDROXYETHYL A® Column, SS, 200 x 4.6mm, 3µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	200	4.6	3	1000	760
204HY0315	PolyHYDROXYETHYL A® Column, SS, 200 x 4.6mm, 3µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	200	4.6	3	1500	760
254HY03006	PolyHYDROXYETHYL A® Column, SS, 250 x 4.6mm, 3µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	250	4.6	3	60	780
254HY0301	PolyHYDROXYETHYL A® Column, SS, 250 x 4.6mm, 3µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	250	4.6	3	100	780
254HY0303	PolyHYDROXYETHYL A® Column, SS, 250 x 4.6mm, 3µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	250	4.6	3	300	780
254HY0305	PolyHYDROXYETHYL A® Column, SS, 250 x 4.6mm, 3µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	250	4.6	3	500	780
254HY0310	PolyHYDROXYETHYL A® Column, SS, 250 x 4.6mm, 3µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	250	4.6	3	1000	780
254HY0315	PolyHYDROXYETHYL A® Column, SS, 250 x 4.6mm, 3µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	250	4.6	3	1500	780
109HY03006	PolyHYDROXYETHYL A® Column, SS, 100 x 9.4mm, 3µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	100	9.4	3	60	1400
109HY0301	PolyHYDROXYETHYL A® Column, SS, 100 x 9.4mm, 3µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	100	9.4	3	100	1400
109HY0303	PolyHYDROXYETHYL A® Column, SS, 100 x 9.4mm, 3µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	100	9.4	3	300	1400
109HY0305	PolyHYDROXYETHYL A® Column, SS, 100 x 9.4mm, 3µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	100	9.4	3	500	1400
109HY0310	PolyHYDROXYETHYL A® Column, SS, 100 x 9.4mm, 3µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	100	9.4	3	1000	1400
109HY0315	PolyHYDROXYETHYL A® Column, SS, 100 x 9.4mm, 3µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	100	9.4	3	1500	1400
209HY03006	PolyHYDROXYETHYL A® Column, SS, 200 x 9.4mm, 3µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	200	9.4	3	60	1620
209HY0301	PolyHYDROXYETHYL A® Column, SS, 200 x 9.4mm, 3µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	200	9.4	3	100	1620
209HY0303	PolyHYDROXYETHYL A® Column, SS, 200 x 9.4mm, 3µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	200	9.4	3	300	1620
209HY0305	PolyHYDROXYETHYL A® Column, SS, 200 x 9.4mm, 3µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	200	9.4	3	500	1620
209HY0310	PolyHYDROXYETHYL A® Column, SS, 200 x 9.4mm, 3µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	200	9.4	3	1000	1620
209HY0315	PolyHYDROXYETHYL A® Column, SS, 200 x 9.4mm, 3µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	200	9.4	3	1500	1620
259HY03006	PolyHYDROXYETHYL A® Column, SS, 250 x 9.4mm, 3µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	250	9.4	3	60	1750
259HY0301	PolyHYDROXYETHYL A® Column, SS, 250 x 9.4mm, 3µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	250	9.4	3	100	1750
259HY0303	PolyHYDROXYETHYL A® Column, SS, 250 x 9.4mm, 3µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	250	9.4	3	300	1750
259HY0305	PolyHYDROXYETHYL A® Column, SS, 250 x 9.4mm, 3µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	250	9.4	3	500	1750
259HY0310	PolyHYDROXYETHYL A® Column, SS, 250 x 9.4mm, 3µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	250	9.4	3	1000	1750

PolyHYDROXYETHYL A® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore size options.) Custom combinations available: info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
259HY0315	PolyHYDROXYETHYL A® Column, SS, 250 x 9.4mm, 3µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	250	9.4	3	1500	1750
309HY03006	PolyHYDROXYETHYL A® Column, SS, 300 x 9.4mm, 3µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	300	9.4	3	60	1850
309HY0301	PolyHYDROXYETHYL A® Column, SS, 300 x 9.4mm, 3µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	300	9.4	3	100	1850
309HY0303	PolyHYDROXYETHYL A® Column, SS, 300 x 9.4mm, 3µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	300	9.4	3	300	1850
309HY0305	PolyHYDROXYETHYL A® Column, SS, 300 x 9.4mm, 3µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	300	9.4	3	500	1850
309HY0310	PolyHYDROXYETHYL A® Column, SS, 300 x 9.4mm, 3µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	300	9.4	3	1000	1850
309HY0315	PolyHYDROXYETHYL A® Column, SS, 300 x 9.4mm, 3µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	300	9.4	3	1500	1850
051HY05006	PolyHYDROXYETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 5µm,60Å	HILIC / SEC	Analytical	GLT	10-32 coned	50	1.0	5	60	570
051HY0501	PolyHYDROXYETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 5µm,100Å	HILIC / SEC	Analytical	GLT	10-32 coned	50	1.0	5	100	570
051HY0502	PolyHYDROXYETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 5µm,200Å	HILIC / SEC	Analytical	GLT	10-32 coned	50	1.0	5	200	570
051HY0503	PolyHYDROXYETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 5µm,300Å	HILIC / SEC	Analytical	GLT	10-32 coned	50	1.0	5	300	570
051HY0505	PolyHYDROXYETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 5µm,500Å	HILIC / SEC	Analytical	GLT	10-32 coned	50	1.0	5	500	570
051HY0510	PolyHYDROXYETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 5µm,1000Å	HILIC / SEC	Analytical	GLT	10-32 coned	50	1.0	5	1000	570
051HY0515	PolyHYDROXYETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 5µm,1500Å	HILIC / SEC	Analytical	GLT	10-32 coned	50	1.0	5	1500	570
101HY05006	PolyHYDROXYETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 5µm,60Å	HILIC / SEC	Analytical	GLT	10-32 coned	100	1.0	5	60	660
101HY0501	PolyHYDROXYETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 5µm,100Å	HILIC / SEC	Analytical	GLT	10-32 coned	100	1.0	5	100	660
101HY0502	PolyHYDROXYETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 5µm,200Å	HILIC / SEC	Analytical	GLT	10-32 coned	100	1.0	5	200	660
101HY0503	PolyHYDROXYETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 5µm,300Å	HILIC / SEC	Analytical	GLT	10-32 coned	100	1.0	5	300	660
101HY0505	PolyHYDROXYETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 5µm,500Å	HILIC / SEC	Analytical	GLT	10-32 coned	100	1.0	5	500	660
101HY0510	PolyHYDROXYETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 5µm,1000Å	HILIC / SEC	Analytical	GLT	10-32 coned	100	1.0	5	1000	660
101HY0515	PolyHYDROXYETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 5µm,1500Å	HILIC / SEC	Analytical	GLT	10-32 coned	100	1.0	5	1500	660
151HY05006	PolyHYDROXYETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 5µm,60Å	HILIC / SEC	Analytical	GLT	10-32 coned	150	1.0	5	60	740
151HY0501	PolyHYDROXYETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 5µm,100Å	HILIC / SEC	Analytical	GLT	10-32 coned	150	1.0	5	100	740
151HY0502	PolyHYDROXYETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 5µm,200Å	HILIC / SEC	Analytical	GLT	10-32 coned	150	1.0	5	200	740
151HY0503	PolyHYDROXYETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 5µm,300Å	HILIC / SEC	Analytical	GLT	10-32 coned	150	1.0	5	300	740
151HY0505	PolyHYDROXYETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 5µm,500Å	HILIC / SEC	Analytical	GLT	10-32 coned	150	1.0	5	500	740
151HY0510	PolyHYDROXYETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 5µm,1000Å	HILIC / SEC	Analytical	GLT	10-32 coned	150	1.0	5	1000	740
151HY0515	PolyHYDROXYETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 5µm,1500Å	HILIC / SEC	Analytical	GLT	10-32 coned	150	1.0	5	1500	740
3.52HY05006	PolyHYDROXYETHYL A® Column, SS, 35 x 2.1mm, 5µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	35	2.1	5	60	500
3.52HY0501	PolyHYDROXYETHYL A® Column, SS, 35 x 2.1mm, 5µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	35	2.1	5	100	500
3.52HY0502	PolyHYDROXYETHYL A® Column, SS, 35 x 2.1mm, 5µm,200Å	HILIC / SEC	Analytical	SS	10-32 coned	35	2.1	5	200	500
3.52HY0503	PolyHYDROXYETHYL A® Column, SS, 35 x 2.1mm, 5µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	35	2.1	5	300	500
3.52HY0505	PolyHYDROXYETHYL A® Column, SS, 35 x 2.1mm, 5µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	35	2.1	5	500	500
3.52HY0510	PolyHYDROXYETHYL A® Column, SS, 35 x 2.1mm, 5µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	35	2.1	5	1000	500
3.52HY0515	PolyHYDROXYETHYL A® Column, SS, 35 x 2.1mm, 5µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	35	2.1	5	1500	500
052HY05006	PolyHYDROXYETHYL A® Column, SS, 50 x 2.1mm, 5µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	50	2.1	5	60	520
052HY0501	PolyHYDROXYETHYL A® Column, SS, 50 x 2.1mm, 5µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	50	2.1	5	100	520
052HY0502	PolyHYDROXYETHYL A® Column, SS, 50 x 2.1mm, 5µm,200Å	HILIC / SEC	Analytical	SS	10-32 coned	50	2.1	5	200	520
052HY0503	PolyHYDROXYETHYL A® Column, SS, 50 x 2.1mm, 5µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	50	2.1	5	300	520
052HY0505	PolyHYDROXYETHYL A® Column, SS, 50 x 2.1mm, 5µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	50	2.1	5	500	520
052HY0510	PolyHYDROXYETHYL A® Column, SS, 50 x 2.1mm, 5µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	50	2.1	5	1000	520
052HY0515	PolyHYDROXYETHYL A® Column, SS, 50 x 2.1mm, 5µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	50	2.1	5	1500	520
102HY05006	PolyHYDROXYETHYL A® Column, SS, 100 x 2.1mm, 5µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	100	2.1	5	60	600
102HY0501	PolyHYDROXYETHYL A® Column, SS, 100 x 2.1mm, 5µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	100	2.1	5	100	600
102HY0502	PolyHYDROXYETHYL A® Column, SS, 100 x 2.1mm, 5µm,200Å	HILIC / SEC	Analytical	SS	10-32 coned	100	2.1	5	200	600
102HY0503	PolyHYDROXYETHYL A® Column, SS, 100 x 2.1mm, 5µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	100	2.1	5	300	600
102HY0505	PolyHYDROXYETHYL A® Column, SS, 100 x 2.1mm, 5µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	100	2.1	5	500	600

PolyHYDROXYETHYL A® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore size options.) Custom combinations available: info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
102HY0510	PolyHYDROXYETHYL A® Column, SS, 100 x 2.1mm, 5µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	100	2.1	5	1000	600
102HY0515	PolyHYDROXYETHYL A® Column, SS, 100 x 2.1mm, 5µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	100	2.1	5	1500	600
152HY05006	PolyHYDROXYETHYL A® Column, SS, 150 x 2.1mm, 5µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	150	2.1	5	60	650
152HY0501	PolyHYDROXYETHYL A® Column, SS, 150 x 2.1mm, 5µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	150	2.1	5	100	650
152HY0502	PolyHYDROXYETHYL A® Column, SS, 150 x 2.1mm, 5µm,200Å	HILIC / SEC	Analytical	SS	10-32 coned	150	2.1	5	200	650
152HY0503	PolyHYDROXYETHYL A® Column, SS, 150 x 2.1mm, 5µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	150	2.1	5	300	650
152HY0505	PolyHYDROXYETHYL A® Column, SS, 150 x 2.1mm, 5µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	150	2.1	5	500	650
152HY0510	PolyHYDROXYETHYL A® Column, SS, 150 x 2.1mm, 5µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	150	2.1	5	1000	650
152HY0515	PolyHYDROXYETHYL A® Column, SS, 150 x 2.1mm, 5µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	150	2.1	5	1500	650
202HY05006	PolyHYDROXYETHYL A® Column, SS, 200 x 2.1mm, 5µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	200	2.1	5	60	710
202HY0501	PolyHYDROXYETHYL A® Column, SS, 200 x 2.1mm, 5µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	200	2.1	5	100	710
202HY0502	PolyHYDROXYETHYL A® Column, SS, 200 x 2.1mm, 5µm,200Å	HILIC / SEC	Analytical	SS	10-32 coned	200	2.1	5	200	710
202HY0503	PolyHYDROXYETHYL A® Column, SS, 200 x 2.1mm, 5µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	200	2.1	5	300	710
202HY0505	PolyHYDROXYETHYL A® Column, SS, 200 x 2.1mm, 5µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	200	2.1	5	500	710
202HY0510	PolyHYDROXYETHYL A® Column, SS, 200 x 2.1mm, 5µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	200	2.1	5	1000	710
202HY0515	PolyHYDROXYETHYL A® Column, SS, 200 x 2.1mm, 5µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	200	2.1	5	1500	710
252HY05006	PolyHYDROXYETHYL A® Column, SS, 250 x 2.1mm, 5µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	250	2.1	5	60	740
252HY0501	PolyHYDROXYETHYL A® Column, SS, 250 x 2.1mm, 5µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	250	2.1	5	100	740
252HY0502	PolyHYDROXYETHYL A® Column, SS, 250 x 2.1mm, 5µm,200Å	HILIC / SEC	Analytical	SS	10-32 coned	250	2.1	5	200	740
252HY0503	PolyHYDROXYETHYL A® Column, SS, 250 x 2.1mm, 5µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	250	2.1	5	300	740
252HY0505	PolyHYDROXYETHYL A® Column, SS, 250 x 2.1mm, 5µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	250	2.1	5	500	740
252HY0510	PolyHYDROXYETHYL A® Column, SS, 250 x 2.1mm, 5µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	250	2.1	5	1000	740
252HY0515	PolyHYDROXYETHYL A® Column, SS, 250 x 2.1mm, 5µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	250	2.1	5	1500	740
3.54HY05006	PolyHYDROXYETHYL A® Column, SS, 35 x 4.6mm, 5µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	35	4.6	5	60	500
3.54HY0501	PolyHYDROXYETHYL A® Column, SS, 35 x 4.6mm, 5µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	35	4.6	5	100	500
3.54HY0502	PolyHYDROXYETHYL A® Column, SS, 35 x 4.6mm, 5µm,200Å	HILIC / SEC	Analytical	SS	10-32 coned	35	4.6	5	200	500
3.54HY0503	PolyHYDROXYETHYL A® Column, SS, 35 x 4.6mm, 5µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	35	4.6	5	300	500
3.54HY0505	PolyHYDROXYETHYL A® Column, SS, 35 x 4.6mm, 5µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	35	4.6	5	500	500
3.54HY0510	PolyHYDROXYETHYL A® Column, SS, 35 x 4.6mm, 5µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	35	4.6	5	1000	500
3.54HY0515	PolyHYDROXYETHYL A® Column, SS, 35 x 4.6mm, 5µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	35	4.6	5	1500	500
054HY05006	PolyHYDROXYETHYL A® Column, SS, 50 x 4.6mm, 5µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	50	4.6	5	60	520
054HY0501	PolyHYDROXYETHYL A® Column, SS, 50 x 4.6mm, 5µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	50	4.6	5	100	520
054HY0502	PolyHYDROXYETHYL A® Column, SS, 50 x 4.6mm, 5µm,200Å	HILIC / SEC	Analytical	SS	10-32 coned	50	4.6	5	200	520
054HY0503	PolyHYDROXYETHYL A® Column, SS, 50 x 4.6mm, 5µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	50	4.6	5	300	520
054HY0505	PolyHYDROXYETHYL A® Column, SS, 50 x 4.6mm, 5µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	50	4.6	5	500	520
054HY0510	PolyHYDROXYETHYL A® Column, SS, 50 x 4.6mm, 5µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	50	4.6	5	1000	520
054HY0515	PolyHYDROXYETHYL A® Column, SS, 50 x 4.6mm, 5µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	50	4.6	5	1500	520
104HY05006	PolyHYDROXYETHYL A® Column, SS, 100 x 4.6mm, 5µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	100	4.6	5	60	600
104HY0501	PolyHYDROXYETHYL A® Column, SS, 100 x 4.6mm, 5µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	100	4.6	5	100	600
104HY0502	PolyHYDROXYETHYL A® Column, SS, 100 x 4.6mm, 5µm,200Å	HILIC / SEC	Analytical	SS	10-32 coned	100	4.6	5	200	600
104HY0503	PolyHYDROXYETHYL A® Column, SS, 100 x 4.6mm, 5µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	100	4.6	5	300	600
104HY0505	PolyHYDROXYETHYL A® Column, SS, 100 x 4.6mm, 5µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	100	4.6	5	500	600
104HY0510	PolyHYDROXYETHYL A® Column, SS, 100 x 4.6mm, 5µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	100	4.6	5	1000	600
104HY0515	PolyHYDROXYETHYL A® Column, SS, 100 x 4.6mm, 5µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	100	4.6	5	1500	600
154HY05006	PolyHYDROXYETHYL A® Column, SS, 150 x 4.6mm, 5µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	150	4.6	5	60	650
154HY0501	PolyHYDROXYETHYL A® Column, SS, 150 x 4.6mm, 5µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	150	4.6	5	100	650
154HY0502	PolyHYDROXYETHYL A® Column, SS, 150 x 4.6mm, 5µm,200Å	HILIC / SEC	Analytical	SS	10-32 coned	150	4.6	5	200	650
154HY0503	PolyHYDROXYETHYL A® Column, SS, 150 x 4.6mm, 5µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	150	4.6	5	300	650

PolyHYDROXYETHYL A® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore size options.) Custom combinations available: info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
154HY0505	PolyHYDROXYETHYL A® Column, SS, 150 x 4.6mm, 5µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	150	4.6	5	500	650
154HY0510	PolyHYDROXYETHYL A® Column, SS, 150 x 4.6mm, 5µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	150	4.6	5	1000	650
154HY0515	PolyHYDROXYETHYL A® Column, SS, 150 x 4.6mm, 5µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	150	4.6	5	1500	650
204HY05006	PolyHYDROXYETHYL A® Column, SS, 200 x 4.6mm, 5µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	200	4.6	5	60	710
204HY0501	PolyHYDROXYETHYL A® Column, SS, 200 x 4.6mm, 5µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	200	4.6	5	100	710
204HY0502	PolyHYDROXYETHYL A® Column, SS, 200 x 4.6mm, 5µm,200Å	HILIC / SEC	Analytical	SS	10-32 coned	200	4.6	5	200	710
204HY0503	PolyHYDROXYETHYL A® Column, SS, 200 x 4.6mm, 5µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	200	4.6	5	300	710
204HY0505	PolyHYDROXYETHYL A® Column, SS, 200 x 4.6mm, 5µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	200	4.6	5	500	710
204HY0510	PolyHYDROXYETHYL A® Column, SS, 200 x 4.6mm, 5µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	200	4.6	5	1000	710
204HY0515	PolyHYDROXYETHYL A® Column, SS, 200 x 4.6mm, 5µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	200	4.6	5	1500	710
254HY05006	PolyHYDROXYETHYL A® Column, SS, 250 x 4.6mm, 5µm,60Å	HILIC / SEC	Analytical	SS	10-32 coned	250	4.6	5	60	740
254HY0501	PolyHYDROXYETHYL A® Column, SS, 250 x 4.6mm, 5µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	250	4.6	5	100	740
254HY0502	PolyHYDROXYETHYL A® Column, SS, 250 x 4.6mm, 5µm,200Å	HILIC / SEC	Analytical	SS	10-32 coned	250	4.6	5	200	740
254HY0503	PolyHYDROXYETHYL A® Column, SS, 250 x 4.6mm, 5µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	250	4.6	5	300	740
254HY0505	PolyHYDROXYETHYL A® Column, SS, 250 x 4.6mm, 5µm,500Å	HILIC / SEC	Analytical	SS	10-32 coned	250	4.6	5	500	740
254HY0510	PolyHYDROXYETHYL A® Column, SS, 250 x 4.6mm, 5µm,1000Å	HILIC / SEC	Analytical	SS	10-32 coned	250	4.6	5	1000	740
254HY0515	PolyHYDROXYETHYL A® Column, SS, 200 x 4.6mm, 5µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	250	4.6	5	1500	740
109HY05006	PolyHYDROXYETHYL A® Column, SS, 100 x 9.4mm, 5µm,60Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	100	9.4	5	60	1110
109HY0501	PolyHYDROXYETHYL A® Column, SS, 100 x 9.4mm, 5µm,100Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	100	9.4	5	100	1110
109HY0502	PolyHYDROXYETHYL A® Column, SS, 100 x 9.4mm, 5µm,200Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	100	9.4	5	200	1110
109HY0503	PolyHYDROXYETHYL A® Column, SS, 100 x 9.4mm, 5µm,300Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	100	9.4	5	300	1110
109HY0505	PolyHYDROXYETHYL A® Column, SS, 100 x 9.4mm, 5µm,500Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	100	9.4	5	500	1110
109HY0510	PolyHYDROXYETHYL A® Column, SS, 100 x 9.4mm, 5µm,1000Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	100	9.4	5	1000	1110
109HY0515	PolyHYDROXYETHYL A® Column, SS, 100 x 9.4mm, 5µm,1500Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	100	9.4	5	1500	1110
209HY05006	PolyHYDROXYETHYL A® Column, SS, 200 x 9.4mm, 5µm,60Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	200	9.4	5	60	1400
209HY0501	PolyHYDROXYETHYL A® Column, SS, 200 x 9.4mm, 5µm,100Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	200	9.4	5	100	1400
209HY0502	PolyHYDROXYETHYL A® Column, SS, 200 x 9.4mm, 5µm,200Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	200	9.4	5	200	1400
209HY0503	PolyHYDROXYETHYL A® Column, SS, 200 x 9.4mm, 5µm,300Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	200	9.4	5	300	1400
209HY0505	PolyHYDROXYETHYL A® Column, SS, 200 x 9.4mm, 5µm,500Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	200	9.4	5	500	1400
209HY0510	PolyHYDROXYETHYL A® Column, SS, 200 x 9.4mm, 5µm,1000Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	200	9.4	5	1000	1400
209HY0515	PolyHYDROXYETHYL A® Column, SS, 200 x 9.4mm, 5µm,1500Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	200	9.4	5	1500	1400
259HY05006	PolyHYDROXYETHYL A® Column, SS, 250 x 9.4mm, 5µm,60Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	250	9.4	5	60	1500
259HY0501	PolyHYDROXYETHYL A® Column, SS, 250 x 9.4mm, 5µm,100Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	250	9.4	5	100	1500
259HY0502	PolyHYDROXYETHYL A® Column, SS, 250 x 9.4mm, 5µm,200Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	250	9.4	5	200	1500
259HY0503	PolyHYDROXYETHYL A® Column, SS, 250 x 9.4mm, 5µm,300Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	250	9.4	5	300	1500
259HY0505	PolyHYDROXYETHYL A® Column, SS, 250 x 9.4mm, 5µm,500Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	250	9.4	5	500	1500
259HY0510	PolyHYDROXYETHYL A® Column, SS, 250 x 9.4mm, 5µm,1000Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	250	9.4	5	1000	1500
259HY0515	PolyHYDROXYETHYL A® Column, SS, 250 x 9.4mm, 5µm,1500Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	250	9.4	5	1500	1500
309HY05006	PolyHYDROXYETHYL A® Column, SS, 300 x 9.4mm, 5µm,60Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	300	9.4	5	60	1600
309HY0501	PolyHYDROXYETHYL A® Column, SS, 300 x 9.4mm, 5µm,100Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	300	9.4	5	100	1600
309HY0502	PolyHYDROXYETHYL A® Column, SS, 300 x 9.4mm, 5µm,200Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	300	9.4	5	200	1600
309HY0503	PolyHYDROXYETHYL A® Column, SS, 300 x 9.4mm, 5µm,300Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	300	9.4	5	300	1600
309HY0505	PolyHYDROXYETHYL A® Column, SS, 300 x 9.4mm, 5µm,500Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	300	9.4	5	500	1600
309HY0510	PolyHYDROXYETHYL A® Column, SS, 300 x 9.4mm, 5µm,1000Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	300	9.4	5	1000	1600
309HY0515	PolyHYDROXYETHYL A® Column, SS, 300 x 9.4mm, 5µm,1500Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	300	9.4	5	1500	1600
2521HY05006	PolyHYDROXYETHYL A® Column, SS, 250x 21.0mm, 5µm,60Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	250	21.0	5	60	4500
2521HY0501	PolyHYDROXYETHYL A® Column, SS, 250x 21.0mm, 5µm,100Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	250	21.0	5	100	4500

PolyHYDROXYETHYL A® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore size options.) Custom combinations available: info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
2521HY0502	PolyHYDROXYETHYL A® Column, SS, 250x 21.0mm, 5µm,200Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	250	21.0	5	200	4500
2521HY0503	PolyHYDROXYETHYL A® Column, SS, 250x 21.0mm, 5µm,300Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	250	21.0	5	300	4500
2521HY0505	PolyHYDROXYETHYL A® Column, SS, 250x 21.0mm, 5µm,500Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	250	21.0	5	500	4500
2521HY0510	PolyHYDROXYETHYL A® Column, SS, 250x 21.0mm, 5µm,1000Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	250	21.0	5	1000	4500
2521HY0515	PolyHYDROXYETHYL A® Column, SS, 250x 21.0mm, 5µm,1500Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	250	21.0	5	1500	4500
204HY1201	PolyHYDROXYETHYL A® Column, SS, 200x 4.6mm, 12µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	200	4.6	12	100	420
204HY1203	PolyHYDROXYETHYL A® Column, SS, 200x 4.6mm, 12µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	200	4.6	12	300	420
204HY1215	PolyHYDROXYETHYL A® Column, SS, 200x 4.6mm, 12µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	200	4.6	12	1500	420
254HY1201	PolyHYDROXYETHYL A® Column, SS, 250x 4.6mm, 12µm,100Å	HILIC / SEC	Analytical	SS	10-32 coned	250	4.6	12	100	500
254HY1203	PolyHYDROXYETHYL A® Column, SS, 250x 4.6mm, 12µm,300Å	HILIC / SEC	Analytical	SS	10-32 coned	250	4.6	12	300	500
254HY1215	PolyHYDROXYETHYL A® Column, SS, 250x 4.6mm, 12µm,1500Å	HILIC / SEC	Analytical	SS	10-32 coned	250	4.6	12	1500	500
209HY1201	PolyHYDROXYETHYL A® Column, SS, 200x 9.4mm, 12µm,100Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	200	9.4	12	100	1000
209HY1203	PolyHYDROXYETHYL A® Column, SS, 200x 9.4mm, 12µm,300Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	200	9.4	12	300	1000
209HY1215	PolyHYDROXYETHYL A® Column, SS, 200x 9.4mm, 12µm,1500Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	200	9.4	12	1500	1000
259HY1201	PolyHYDROXYETHYL A® Column, SS, 250x 9.4mm, 12µm,100Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	250	9.4	12	100	1100
259HY1203	PolyHYDROXYETHYL A® Column, SS, 250x 9.4mm, 12µm,300Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	250	9.4	12	300	1100
259HY1215	PolyHYDROXYETHYL A® Column, SS, 250x 9.4mm, 12µm,1500Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	250	9.4	12	1500	1100
309HY1201	PolyHYDROXYETHYL A® Column, SS, 300x 9.4mm, 12µm,100Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	300	9.4	12	100	1300
309HY1203	PolyHYDROXYETHYL A® Column, SS, 300x 9.4mm, 12µm,300Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	300	9.4	12	300	1300
309HY1215	PolyHYDROXYETHYL A® Column, SS, 300x 9.4mm, 12µm,1500Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	300	9.4	12	1500	1300
2521HY1201	PolyHYDROXYETHYL A® Column, SS, 250x 21.0mm, 12µm,100Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	250	21.0	12	100	3200
2521HY1203	PolyHYDROXYETHYL A® Column, SS, 250x 21.0mm, 12µm,300Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	250	21.0	12	300	3200
2521HY1215	PolyHYDROXYETHYL A® Column, SS, 250x 21.0mm, 12µm,1500Å	HILIC / SEC	Semi-Prep	SS	10-32 coned	250	21.0	12	1500	3200
2550HY1201	PolyHYDROXYETHYL A® Column, SS, 250x 50.8mm, 12µm,100Å	HILIC / SEC	Prep	SS	10-32 coned	250	50.8	12	100	12250
2550HY1203	PolyHYDROXYETHYL A® Column, SS, 250x 50.8mm, 12µm,300Å	HILIC / SEC	Prep	SS	10-32 coned	250	50.8	12	300	12250
2550HY1215	PolyHYDROXYETHYL A® Column, SS, 250x 80.5mm, 12µm,1500Å	HILIC / SEC	Prep	SS	10-32 coned	250	50.8	12	1500	12250



HPLC Capillaries

PolyHYDROXYETHYL A® Capillaries (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes. Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
050.075HY0201	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 2µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.075	2	100	520
050.075HY0203	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 2µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.075	2	300	520
050.075HY0210	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 2µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.075	2	1000	520
100.075HY0201	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 2µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.075	2	100	580
100.075HY0203	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 2µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.075	2	300	580
100.075HY0210	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 2µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.075	2	1000	580
150.075HY0201	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 2µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.075	2	100	640
150.075HY0203	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 2µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.075	2	300	640
150.075HY0210	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 2µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.075	2	1000	640

PolyHYDROXYETHYL A® Capillaries (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes. Custom combinations available upon request info@polylc.com)

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
200.075HY0201	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 2µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.075	2	100	700
200.075HY0203	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 2µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.075	2	300	700
200.075HY0210	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 2µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.075	2	1000	700
050.15HY0201	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 2µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.150	2	100	520
050.15HY0203	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 2µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.150	2	300	520
050.15HY0210	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 2µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.150	2	1000	520
100.15HY0201	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 2µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.150	2	100	580
100.15HY0203	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 2µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.150	2	300	580
100.15HY0210	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 2µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.150	2	1000	580
150.15HY0201	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 2µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.150	2	100	640
150.15HY0203	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 2µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.150	2	300	640
150.15HY0210	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 2µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.150	2	1000	640
200.15HY0201	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 2µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.150	2	100	700
200.15HY0203	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 2µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.150	2	300	700
200.15HY0210	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 2µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.150	2	1000	700
050.30HY0201	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 2µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.300	2	100	520
050.30HY0203	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 2µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.300	2	300	520
050.30HY0210	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 2µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.300	2	1000	520
100.30HY0201	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 2µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.300	2	100	580
100.30HY0203	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 2µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.300	2	300	580
100.30HY0210	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 2µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.300	2	1000	580
150.30HY0201	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 2µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.300	2	100	640
150.30HY0203	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 2µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.300	2	300	640
150.30HY0210	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 2µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.300	2	1000	640
200.30HY0201	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 2µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.300	2	100	700
200.30HY0203	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 2µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.300	2	300	700
200.30HY0210	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 2µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.300	2	1000	700
050.075HY03006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.075	3	60	520
050.075HY0301	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.075	3	100	520
050.075HY0302	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.075	3	200	520
050.075HY0303	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.075	3	300	520
050.075HY0305	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.075	3	500	520
050.075HY0310	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.075	3	1000	520
050.075HY0315	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.075	3	1500	520
100.075HY03006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.075	3	60	580
100.075HY0301	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.075	3	100	580
100.075HY0302	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.075	3	200	580
100.075HY0303	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.075	3	300	580
100.075HY0305	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.075	3	500	580
100.075HY0310	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.075	3	1000	580
100.075HY0315	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.075	3	1500	580
150.075HY03006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.075	3	60	640
150.075HY0301	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.075	3	100	640
150.075HY0302	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.075	3	200	640
150.075HY0303	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.075	3	300	640
150.075HY0305	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.075	3	500	640
150.075HY0310	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.075	3	1000	640

PolyHYDROXYETHYL A® Capillaries (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes. Custom combinations available upon request info@polylc.com)

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
150.075HY0315	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.075	3	1500	640
200.075HY03006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.075	3	60	700
200.075HY0301	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.075	3	100	700
200.075HY0302	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.075	3	200	700
200.075HY0303	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.075	3	300	700
200.075HY0305	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.075	3	500	700
200.075HY0310	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.075	3	1000	700
200.075HY0315	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.075	3	1500	700
050.15HY03006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.150	3	60	520
050.15HY0301	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.150	3	100	520
050.15HY0302	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.150	3	200	520
050.15HY0303	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.150	3	300	520
050.15HY0305	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.150	3	500	520
050.15HY0310	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.150	3	1000	520
050.15HY0315	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.150	3	1500	520
100.15HY03006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.150	3	60	580
100.15HY0301	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.150	3	100	580
100.15HY0302	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.150	3	200	580
100.15HY0303	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.150	3	300	580
100.15HY0305	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.150	3	500	580
100.15HY0310	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.150	3	1000	580
100.15HY0315	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.150	3	1500	580
150.15HY03006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.150	3	60	640
150.15HY0301	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.150	3	100	640
150.15HY0302	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.150	3	200	640
150.15HY0303	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.150	3	300	640
150.15HY0305	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.150	3	500	640
150.15HY0310	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.150	3	1000	640
150.15HY0315	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.150	3	1500	640
200.15HY03006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.150	3	60	700
200.15HY0301	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.150	3	100	700
200.15HY0302	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.150	3	200	700
200.15HY0303	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.150	3	300	700
200.15HY0305	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.150	3	500	700
200.15HY0310	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.150	3	1000	700
200.15HY0315	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.150	3	1500	700
050.30HY03006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.300	3	60	520
050.30HY0301	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.300	3	100	520
050.30HY0302	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.300	3	200	520
050.30HY0303	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.300	3	300	520
050.30HY0305	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.300	3	500	520
050.30HY0310	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.300	3	1000	520
050.30HY0315	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.300	3	1500	520
100.30HY03006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.300	3	60	580
100.30HY0301	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.300	3	100	580
100.30HY0302	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.300	3	200	580
100.30HY0303	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.300	3	300	580

PolyHYDROXYETHYL A® Capillaries (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes. Custom combinations available upon request info@polylc.com)

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
100.30HY0305	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.300	3	500	580
100.30HY0310	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.300	3	1000	580
100.30HY0315	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.300	3	1500	580
150.30HY03006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.300	3	60	640
150.30HY0301	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.300	3	100	640
150.30HY0302	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.300	3	200	640
150.30HY0303	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.300	3	300	640
150.30HY0305	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.300	3	500	640
150.30HY0310	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.300	3	1000	640
150.30HY0315	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.300	3	1500	640
200.30HY03006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.300	3	60	700
200.30HY0301	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.300	3	100	700
200.30HY0302	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.300	3	200	700
200.30HY0303	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.300	3	300	700
200.30HY0305	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.300	3	500	700
200.30HY0310	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.300	3	1000	700
200.30HY0315	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.300	3	1500	700
050.075HY05006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 5µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.075	5	60	520
050.075HY0501	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 5µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.075	5	100	520
050.075HY0502	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 5µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.075	5	200	520
050.075HY0503	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 5µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.075	5	300	520
050.075HY0505	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 5µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.075	5	500	520
050.075HY0510	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 5µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.075	5	1000	520
050.075HY0515	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 5µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.075	5	1500	520
100.075HY05006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 5µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.075	5	60	580
100.075HY0501	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 5µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.075	5	100	580
100.075HY0502	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 5µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.075	5	200	580
100.075HY0503	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 5µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.075	5	300	580
100.075HY0505	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 5µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.075	5	500	580
100.075HY0510	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 5µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.075	5	1000	580
100.075HY0515	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 5µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.075	5	1500	580
150.075HY05006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 5µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.075	5	60	640
150.075HY0501	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 5µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.075	5	100	640
150.075HY0502	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 5µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.075	5	200	640
150.075HY0503	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 5µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.075	5	300	640
150.075HY0505	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 5µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.075	5	500	640
150.075HY0510	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 5µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.075	5	1000	640
150.075HY0515	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 5µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.075	5	1500	640
200.075HY05006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 5µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.075	5	60	700
200.075HY0501	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 5µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.075	5	100	700
200.075HY0502	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 5µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.075	5	200	700
200.075HY0503	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 5µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.075	5	300	700
200.075HY0505	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 5µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.075	5	500	700
200.075HY0510	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 5µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.075	5	1000	700
200.075HY0515	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 5µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.075	5	1500	700
050.15HY05006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 5µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.150	5	60	520
050.15HY0501	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 5µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.150	5	100	520

PolyHYDROXYETHYL A® Capillaries (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes. Custom combinations available upon request info@polylc.com)

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
050.15HY0502	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 5µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.150	5	200	520
050.15HY0503	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 5µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.150	5	300	520
050.15HY0505	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 5µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.150	5	500	520
050.15HY0510	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 5µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.150	5	1000	520
050.15HY0515	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 5µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.150	5	1500	520
100.15HY05006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 5µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.150	5	60	580
100.15HY0501	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 5µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.150	5	100	580
100.15HY0502	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 5µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.150	5	200	580
100.15HY0503	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 5µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.150	5	300	580
100.15HY0505	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 5µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.150	5	500	580
100.15HY0510	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 5µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.150	5	1000	580
100.15HY0515	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 5µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.150	5	1500	580
150.15HY05006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 5µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.150	5	60	640
150.15HY0501	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 5µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.150	5	100	640
150.15HY0502	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 5µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.150	5	200	640
150.15HY0503	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 5µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.150	5	300	640
150.15HY0505	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 5µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.150	5	500	640
150.15HY0510	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 5µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.150	5	1000	640
150.15HY0515	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 5µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.150	5	1500	640
200.15HY05006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 5µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.150	5	60	700
200.15HY0501	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 5µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.150	5	100	700
200.15HY0502	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 5µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.150	5	200	700
200.15HY0503	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 5µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.150	5	300	700
200.15HY0505	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 5µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.150	5	500	700
200.15HY0510	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 5µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.150	5	1000	700
200.15HY0515	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 5µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.150	5	1500	700
050.30HY05006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 5µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.300	5	60	520
050.30HY0501	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 5µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.300	5	100	520
050.30HY0502	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 5µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.300	5	200	520
050.30HY0503	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 5µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.300	5	300	520
050.30HY0505	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 5µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.300	5	500	520
050.30HY0510	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 5µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.300	5	1000	520
050.30HY0515	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 5µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	50	0.300	5	1500	520
100.30HY05006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 5µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.300	5	60	580
100.30HY0501	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 5µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.300	5	100	580
100.30HY0502	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 5µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.300	5	200	580
100.30HY0503	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 5µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.300	5	300	580
100.30HY0505	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 5µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.300	5	500	580
100.30HY0510	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 5µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.300	5	1000	580
100.30HY0515	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 5µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	100	0.300	5	1500	580
150.30HY05006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 5µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.300	5	60	640
150.30HY0501	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 5µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.300	5	100	640
150.30HY0502	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 5µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.300	5	200	640
150.30HY0503	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 5µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.300	5	300	640
150.30HY0505	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 5µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.300	5	500	640
150.30HY0510	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 5µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.300	5	1000	640
150.30HY0515	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 5µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	150	0.300	5	1500	640

PolyHYDROXYETHYL A® Capillaries (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes. Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
200.30HY05006	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 5µm,60Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.300	5	60	700
200.30HY0501	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 5µm,100Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.300	5	100	700
200.30HY0502	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 5µm,200Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.300	5	200	700
200.30HY0503	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 5µm,300Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.300	5	300	700
200.30HY0505	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 5µm,500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.300	5	500	700
200.30HY0510	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 5µm,1000Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.300	5	1000	700
200.30HY0515	PolyHYDROXYETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 5µm,1500Å	HILIC / SEC	Analytical	PEEK SIL	10-32 coned	200	0.300	5	1500	700

Guard Cartridges



PolyHYDROXYETHYL A® Guard Cartridges (Available in 2, 3, 5, 12, 20, and 30µm particles with pore sizes.) Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
GC1-HY0201	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 2µm, 100Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	1.0	2	100	120
GC1-HY0203	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 2µm, 300Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	1.0	2	300	120
GC1-HY0210	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 2µm, 1000Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	1.0	2	1000	120
GC2-HY0201	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2mm, 2µm, 100Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	2.0	2	100	120
GC2-HY0203	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2mm, 2µm, 300Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	2.0	2	300	120
GC2-HY0210	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2mm, 2µm, 1000Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	2.0	2	1000	120
GC3-HY0201	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 2µm, 100Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	3.0	2	100	120
GC3-HY0203	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 2µm, 300Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	3.0	2	300	120
GC3-HY0210	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 2µm, 1000Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	3.0	2	1000	120
GC4-HY0201	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 4mm, 2µm, 100Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	4.0	2	100	120
GC4-HY0203	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 4mm, 2µm, 300Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	4.0	2	300	120
GC4-HY0210	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 4mm, 2µm, 1000Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	4.0	2	1000	120
GC1-HY03006	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 60Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	1.0	3	60	100
GC1-HY0301	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 100Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	1.0	3	100	100
GC1-HY0302	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 200Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	1.0	3	200	100
GC1-HY0303	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 300Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	1.0	3	300	100
GC1-HY0305	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 500Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	1.0	3	500	100
GC1-HY0310	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 1000Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	1.0	3	1000	100
GC1-HY0315	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 1500Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	1.0	3	1500	100
GC2-HY03006	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2mm, 3µm, 60Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	2.0	3	60	100
GC2-HY0301	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2mm, 3µm, 100Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	2.0	3	100	100
GC2-HY0302	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2mm, 3µm, 200Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	2.0	3	200	100
GC2-HY0303	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2mm, 3µm, 300Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	2.0	3	300	100
GC2-HY0305	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2mm, 3µm, 500Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	2.0	3	500	100
GC2-HY0310	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2mm, 3µm, 1000Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	2.0	3	1000	100
GC2-HY0315	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2mm, 3µm, 1500Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	2.0	3	1500	100
GC3-HY03006	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 3µm, 60Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	3.0	3	60	100

PolyHYDROXYETHYL A® Guard Cartridges (Available in 2, 3, 5, 12, 20, and 30µm particles with pore sizes.) Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
GC3-HY0301	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 3µm, 100Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	3.0	3	100	100
GC3-HY0302	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 3µm, 200Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	3.0	3	200	100
GC3-HY0303	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 3µm, 300Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	3.0	3	300	100
GC3-HY0305	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 3µm, 500Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	3.0	3	500	100
GC3-HY0310	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 3µm, 1000Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	3.0	3	1000	100
GC3-HY0315	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 3µm, 1500Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	3.0	3	1500	100
GC4-HY03006	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 4mm, 3µm, 60Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	4.0	3	60	100
GC4-HY0301	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 4mm, 3µm, 100Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	4.0	3	100	100
GC4-HY0302	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 4mm, 3µm, 200Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	4.0	3	200	100
GC4-HY0303	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 4mm, 3µm, 300Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	4.0	3	300	100
GC4-HY0305	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 4mm, 3µm, 500Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	4.0	3	500	100
GC4-HY0310	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 4mm, 3µm, 1000Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	4.0	3	1000	100
GC4-HY0315	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 4mm, 3µm, 1500Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	4.0	3	1500	100
GC1-HY05006	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 5µm, 60Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	1.0	5	60	90
GC1-HY0501	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 5µm, 100Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	1.0	5	100	90
GC1-HY0502	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 5µm, 200Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	1.0	5	200	90
GC1-HY0503	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 5µm, 300Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	1.0	5	300	90
GC1-HY0505	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 5µm, 500Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	1.0	5	500	90
GC1-HY0510	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 5µm, 1000Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	1.0	5	1000	90
GC1-HY0515	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 5µm, 1500Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	1.0	5	1500	90
GC2-HY05006	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2mm, 5µm, 60Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	2.0	5	60	90
GC2-HY0501	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2mm, 5µm, 100Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	2.0	5	100	90
GC2-HY0502	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2mm, 5µm, 200Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	2.0	5	200	90
GC2-HY0503	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2mm, 5µm, 300Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	2.0	5	300	90
GC2-HY0505	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2mm, 5µm, 500Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	2.0	5	500	90
GC2-HY0510	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2mm, 5µm, 1000Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	2.0	5	1000	90
GC2-HY0515	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2mm, 5µm, 1500Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	2.0	5	1500	90
GC3-HY05006	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 5µm, 60Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	3.0	5	60	90
GC3-HY0501	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 5µm, 100Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	3.0	5	100	90
GC3-HY0502	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 5µm, 200Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	3.0	5	200	90
GC3-HY0503	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 5µm, 300Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	3.0	5	300	90
GC3-HY0505	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 5µm, 500Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	3.0	5	500	90
GC3-HY0510	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 5µm, 1000Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	3.0	5	1000	90
GC3-HY0515	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 5µm, 1500Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	3.0	5	1500	90
GC4-HY05006	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 5µm, 60Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	4.0	5	60	90
GC4-HY0501	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 5µm, 100Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	4.0	5	100	90
GC4-HY0502	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 5µm, 200Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	4.0	5	200	90
GC4-HY0503	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 5µm, 300Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	4.0	5	300	90
GC4-HY0505	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 5µm, 500Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	4.0	5	500	90
GC4-HY0510	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 5µm, 1000Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	4.0	5	1000	90
GC4-HY0515	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3mm, 5µm, 1500Å	HILIC / SEC	Analytical	SS/PEEK	10-32 coned	10	4.0	5	1500	90
DGCH10	Direct Guard Cartridge Holder, 10mm (For all analytical guards, REUSABLE)	N/A	Analytical	SS/PEEK	10-32 coned	10	N/A	N/A	N/A	75
IGCH10	Inline Guard Cartridge Holder, 10mm (For all analytical guards, REUSABLE)	N/A	Analytical	SS/PEEK	10-32 coned	10	N/A	N/A	N/A	75
*** NOTE: DGCH10 for Direct Connection to Column (female inlet/male outlet) IGCH10 for Connection to Line (two female ports)										
GC10-HY03006	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 60Å	HILIC / SEC	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	60	225

PolyHYDROXYETHYL A® Guard Cartridges (Available in 2, 3, 5, 12, 20, and 30µm particles with pore sizes.) Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
GC10-HY0301	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 100Å	HILIC / SEC	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	100	225
GC10-HY0302	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 200Å	HILIC / SEC	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	200	225
GC10-HY0303	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 300Å	HILIC / SEC	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	300	225
GC10-HY0305	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 500Å	HILIC / SEC	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	500	225
GC10-HY0310	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 1000Å	HILIC / SEC	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	1000	225
GC10-HY0315	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 1500Å	HILIC / SEC	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	1500	225
GC10-HY05006	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 5µm, 60Å	HILIC / SEC	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	5	60	185
GC10-HY0501	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 5µm, 100Å	HILIC / SEC	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	5	100	185
GC10-HY0502	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 5µm, 200Å	HILIC / SEC	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	5	200	185
GC10-HY0503	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 5µm, 300Å	HILIC / SEC	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	5	300	185
GC10-HY0505	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 5µm, 500Å	HILIC / SEC	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	5	500	185
GC10-HY0510	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 5µm, 1000Å	HILIC / SEC	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	5	1000	185
GC10-HY0515	PolyHYDROXYETHYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 5µm, 1500Å	HILIC / SEC	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	5	1500	185
CH10	Semi-Prep Guard Cartridge Holder, 10mm (For 9.4mm columns, REUSABLE)	HILIC / SEC	Semi-Prep	SS/PEEK	10-32 coned	10	N/A	N/A	N/A	175

*** NOTE: CH10 Holder Required for Semi-Prep Guard Cartridges (10 x 10mm) Connection to Line (two female ports)

Solid Phase Extraction Cartridges




PolyHYDROXYETHYL A® SPE Cartridges

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	V (ml)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
SPEHY1201	PolyHYDROXYETHYL A® SPE Cartridge 1mL, PE, 12µm, 100Å (Pack of 10 cartridges)	HILIC / SEC	SPE	PE	Lure	1	N/A	12	100	130
SPEHY1203	PolyHYDROXYETHYL A® SPE Cartridge 1mL, PE, 12µm, 300Å (Pack of 10 cartridges)	HILIC / SEC	SPE	PE	Lure	1	N/A	12	300	130
SPEHY1215	PolyHYDROXYETHYL A® SPE Cartridge 1mL, PE, 12µm, 1500Å (Pack of 10 cartridges)	HILIC / SEC	SPE	PE	Lure	1	N/A	12	1500	130
SPEHY2001	PolyHYDROXYETHYL A® SPE Cartridge 1mL, PE, 20µm, 100Å (Pack of 10 cartridges)	HILIC / SEC	SPE	PE	Lure	1	N/A	20	100	130
SPEHY2003	PolyHYDROXYETHYL A® SPE Cartridge 1mL, PE, 20µm, 300Å (Pack of 10 cartridges)	HILIC / SEC	SPE	PE	Lure	1	N/A	20	300	130
SPEHY2015	PolyHYDROXYETHYL A® SPE Cartridge 1mL, PE, 20µm, 1500Å (Pack of 10 cartridges)	HILIC / SEC	SPE	PE	Lure	1	N/A	20	1500	130

Spin Columns



PolyHYDROXYETHYL A® BioPureSPN Columns

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Bed (mg)	Pore (Å)		HILIC Capacity	2026 Price (USD)
HMM HIL	PolyHYDROXYETHYL A® BioPureSPN MACRO for HILIC 50-200µl (Box of 100 each)	HILIC / SEC	Spin	100	300		35 - 350µg	940
HMM HIL.20	PolyHYDROXYETHYL A® BioPureSPN MACRO for HILIC 50-200µl (Pack of 20 each)	HILIC / SEC	Spin	100	300		35 - 350µg	265
HNS HIL-L	PolyHYDROXYETHYL A® BioPureSPN MACRO PLATE 96 Well for HILIC 50-200µl	HILIC / SEC	Spin	100	300		35 - 350µg	800
HNS HIL60-L	PolyHYDROXYETHYL A® BioPureSPN MIDI PLATE 96 Well for HILIC 50-200µl	HILIC / SEC	Spin	100	60		35 - 350µg	800
HEM HIL	PolyHYDROXYETHYL A® BioPureSPN MIDI for HILIC 50-200µl (Box of 100 each)	HILIC / SEC	Spin	50	300		17 - 170µg	820
HEM HIL.20	PolyHYDROXYETHYL A® BioPureSPN MIDI for HILIC 50-200µl (Pack of 20 each)	HILIC / SEC	Spin	50	300		17 - 170µg	240
HNS HIL-M	PolyHYDROXYETHYL A® BioPureSPN MIDI PLATE 96 Well for HILIC 50-200µl	HILIC / SEC	Spin	50	300		17 - 170µg	820
HUM HIL	PolyHYDROXYETHYL A® BioPureSPN MINI for HILIC 20-100µl (Box of 100 each)	HILIC / SEC	Spin	20	300		7 - 70µg	735
HUM HIL.20	PolyHYDROXYETHYL A® BioPureSPN MINI for HILIC 20-100µl (Pack of 20 each)	HILIC / SEC	Spin	20	300		7 - 70µg	225
HNS HIL	PolyHYDROXYETHYL A® BioPureSPN MINI PLATE 96 Well for HILIC 20-100µl	HILIC / SEC	Spin	20	300		7 - 70µg	560

Bulk Material



PolyHYDROXYETHYL A® Bulk Material (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.) Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	QTY (g)	Particle (µm)	Pore (Å)	2026 Price (USD)
BMHY0201	PolyHYDROXYETHYL A® Bulk Material, 2µm, 100Å	HILIC / SEC	BM	N/A	0.5	2	100	100
BMHY0203	PolyHYDROXYETHYL A® Bulk Material, 2µm, 300Å	HILIC / SEC	BM	N/A	0.5	2	300	100
BMHY0210	PolyHYDROXYETHYL A® Bulk Material, 2µm, 1000Å	HILIC / SEC	BM	N/A	0.5	2	1000	100
BMHY03006	PolyHYDROXYETHYL A® Bulk Material, 3µm, 60Å	HILIC / SEC	BM	N/A	1	3	60	130
BMHY0301	PolyHYDROXYETHYL A® Bulk Material, 3µm, 100Å	HILIC / SEC	BM	N/A	1	3	100	130
BMHY0302	PolyHYDROXYETHYL A® Bulk Material, 3µm, 200Å	HILIC / SEC	BM	N/A	1	3	200	130
BMHY0303	PolyHYDROXYETHYL A® Bulk Material, 3µm, 300Å	HILIC / SEC	BM	N/A	1	3	300	130
BMHY0304	PolyHYDROXYETHYL A® Bulk Material, 3µm, 500Å	HILIC / SEC	BM	N/A	1	3	500	130
BMHY0310	PolyHYDROXYETHYL A® Bulk Material, 3µm, 1000Å	HILIC / SEC	BM	N/A	1	3	1000	130
BMHY0315	PolyHYDROXYETHYL A® Bulk Material, 3µm, 1500Å	HILIC / SEC	BM	N/A	1	3	1500	130
BMHY05006	PolyHYDROXYETHYL A® Bulk Material, 5µm, 60Å	HILIC / SEC	BM	N/A	1	5	60	85
BMHY0501	PolyHYDROXYETHYL A® Bulk Material, 5µm, 100Å	HILIC / SEC	BM	N/A	1	5	100	85
BMHY0502	PolyHYDROXYETHYL A® Bulk Material, 5µm, 200Å	HILIC / SEC	BM	N/A	1	5	200	85
BMHY0503	PolyHYDROXYETHYL A® Bulk Material, 5µm, 300Å	HILIC / SEC	BM	N/A	1	5	300	85
BMHY0504	PolyHYDROXYETHYL A® Bulk Material, 5µm, 500Å	HILIC / SEC	BM	N/A	1	5	500	85
BMHY0510	PolyHYDROXYETHYL A® Bulk Material, 5µm, 1000Å	HILIC / SEC	BM	N/A	1	5	1000	85
BMCT0503	PolyHYDROXYETHYL A® Bulk Material, 5µm, 1500Å	HILIC / SEC	BM	N/A	1	5	1500	85
BMHY1201	PolyHYDROXYETHYL A® Bulk Material, 12µm, 100Å	HILIC / SEC	BM	N/A	1	12	100	32
BMHY1203	PolyHYDROXYETHYL A® Bulk Material, 12µm, 300Å	HILIC / SEC	BM	N/A	1	12	300	32
BMHY1215	PolyHYDROXYETHYL A® Bulk Material, 12µm, 1500Å	HILIC / SEC	BM	N/A	1	12	1500	32
BMHY2001	PolyHYDROXYETHYL A® Bulk Material, 20µm, 100Å	HILIC / SEC	BM	N/A	1	20	100	32
BMHY2003	PolyHYDROXYETHYL A® Bulk Material, 20µm, 300Å	HILIC / SEC	BM	N/A	1	20	300	32
BMHY2015	PolyHYDROXYETHYL A® Bulk Material, 20µm, 1500Å	HILIC / SEC	BM	N/A	1	20	1500	32
BMHY3003	PolyHYDROXYETHYL A® Bulk Material, 30µm, 300Å	HILIC / SEC	BM	N/A	1	30	300	32
BMHY3010	PolyHYDROXYETHYL A® Bulk Material, 30µm, 1000Å	HILIC / SEC	BM	N/A	1	30	1000	32

Please contact info@polylc.com with any questions or for information on any product combination not listed in the table above

PolyHYDROXYETHYL A®

(HILIC / SEC)

2026 PolyLC Inc. Retail Pricelist

PolyLC Inc., 9151 Rumsey Rd., Suite 180, Columbia, MD 21045 410-992-5400 www.polylc.com



2026 PolyLC Inc. Retail Pricelist



Offering high-performance liquid chromatography media, columns, capillaries, cartridges, and instrumentation. We specialize in advanced alternatives to reversed-phase for proteins, peptides, and other biomolecules, supported by deep application expertise and customized solutions for the most demanding separations.

PolySULFOETHYL A®

(Strong-Cation Exchange)



HPLC Columns

PolySULFOETHYL A® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.) Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
051SE0203	PolySULFOETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 2µm,300Å	SCX	Analytical	GLT	10-32 coned	50	1.0	2	300	615
051SE0210	PolySULFOETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 2µm,1000Å	SCX	Analytical	GLT	10-32 coned	50	1.0	2	1000	615
051SE0215	PolySULFOETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 2µm,1500Å	SCX	Analytical	GLT	10-32 coned	50	1.0	2	1500	615
101SE0203	PolySULFOETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 2µm,300Å	SCX	Analytical	GLT	10-32 coned	100	1.0	2	300	750
101SE0210	PolySULFOETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 2µm,1000Å	SCX	Analytical	GLT	10-32 coned	100	1.0	2	1000	750
101SE0215	PolySULFOETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 2µm,1500Å	SCX	Analytical	GLT	10-32 coned	100	1.0	2	1500	750
151SE0203	PolySULFOETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 2µm,300Å	SCX	Analytical	GLT	10-32 coned	150	1.0	2	300	800
151SE0210	PolySULFOETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 2µm,1000Å	SCX	Analytical	GLT	10-32 coned	150	1.0	2	1000	800
151SE0215	PolySULFOETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 2µm,1500Å	SCX	Analytical	GLT	10-32 coned	150	1.0	2	1500	800
3.52SE0203	PolySULFOETHYL A® Column, SS, 35 x 2.1mm, 2µm,300Å	SCX	Analytical	SS	10-32 coned	35	2.1	2	300	615
3.52SE0210	PolySULFOETHYL A® Column, SS, 35 x 2.1mm, 2µm,1000Å	SCX	Analytical	SS	10-32 coned	35	2.1	2	1000	615
3.52SE0215	PolySULFOETHYL A® Column, SS, 35 x 2.1mm, 2µm,1500Å	SCX	Analytical	SS	10-32 coned	35	2.1	2	1500	615
052SE0203	PolySULFOETHYL A® Column, SS, 50 x 2.1mm, 2µm,300Å	SCX	Analytical	SS	10-32 coned	50	2.1	2	300	625
052SE0210	PolySULFOETHYL A® Column, SS, 50 x 2.1mm, 2µm,1000Å	SCX	Analytical	SS	10-32 coned	50	2.1	2	1000	625
052SE0215	PolySULFOETHYL A® Column, SS, 50 x 2.1mm, 2µm,1500Å	SCX	Analytical	SS	10-32 coned	50	2.1	2	1500	625
102SE0203	PolySULFOETHYL A® Column, SS, 100 x 2.1mm, 2µm,300Å	SCX	Analytical	SS	10-32 coned	100	2.1	2	300	800
102SE0210	PolySULFOETHYL A® Column, SS, 100 x 2.1mm, 2µm,1000Å	SCX	Analytical	SS	10-32 coned	100	2.1	2	1000	800
102SE0215	PolySULFOETHYL A® Column, SS, 100 x 2.1mm, 2µm,1500Å	SCX	Analytical	SS	10-32 coned	100	2.1	2	1500	800
3.54SE0203	PolySULFOETHYL A® Column, SS, 35 x 4.6mm, 2µm,300Å	SCX	Analytical	SS	10-32 coned	35	4.6	2	300	615
3.54SE0210	PolySULFOETHYL A® Column, SS, 35 x 4.6mm, 2µm,1000Å	SCX	Analytical	SS	10-32 coned	35	4.6	2	1000	615
3.54SE0215	PolySULFOETHYL A® Column, SS, 35 x 4.6mm, 2µm,1500Å	SCX	Analytical	SS	10-32 coned	35	4.6	2	1500	615
054SE0203	PolySULFOETHYL A® Column, SS, 50 x 4.6mm, 2µm,300Å	SCX	Analytical	SS	10-32 coned	50	4.6	2	300	625
054SE0210	PolySULFOETHYL A® Column, SS, 50 x 4.6mm, 2µm,1000Å	SCX	Analytical	SS	10-32 coned	50	4.6	2	1000	625
054SE0215	PolySULFOETHYL A® Column, SS, 50 x 4.6mm, 2µm,1500Å	SCX	Analytical	SS	10-32 coned	50	4.6	2	1500	625
104SE0203	PolySULFOETHYL A® Column, SS, 100 x 4.6mm, 2µm,300Å	SCX	Analytical	SS	10-32 coned	100	4.6	2	300	800
104SE0210	PolySULFOETHYL A® Column, SS, 100 x 4.6mm, 2µm,1000Å	SCX	Analytical	SS	10-32 coned	100	4.6	2	1000	800
104SE0215	PolySULFOETHYL A® Column, SS, 100 x 4.6mm, 2µm,1500Å	SCX	Analytical	SS	10-32 coned	100	4.6	2	1500	800
154SE0203	PolySULFOETHYL A® Column, SS, 150 x 4.6mm, 2µm,300Å	SCX	Analytical	SS	10-32 coned	150	4.6	2	300	830
154SE0210	PolySULFOETHYL A® Column, SS, 150 x 4.6mm, 2µm,1000Å	SCX	Analytical	SS	10-32 coned	150	4.6	2	1000	830
154SE0215	PolySULFOETHYL A® Column, SS, 150 x 4.6mm, 2µm,1500Å	SCX	Analytical	SS	10-32 coned	150	4.6	2	1500	830
204SE0203	PolySULFOETHYL A® Column, SS, 200 x 4.6mm, 2µm,300Å	SCX	Analytical	SS	10-32 coned	200	4.6	2	300	860
204SE0210	PolySULFOETHYL A® Column, SS, 200 x 4.6mm, 2µm,1000Å	SCX	Analytical	SS	10-32 coned	200	4.6	2	1000	860
204SE0215	PolySULFOETHYL A® Column, SS, 200 x 4.6mm, 2µm,1500Å	SCX	Analytical	SS	10-32 coned	200	4.6	2	1500	860

PolySULFOETHYL A® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.) Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
109SE0203	PolySULFOETHYL A® Column, SS, 100 x 9.4mm, 2µm,300Å	SCX	Semi-Prep	SS	10-32 coned	100	9.4	2	300	1620
109SE0210	PolySULFOETHYL A® Column, SS, 100 x 9.4mm, 2µm,1000Å	SCX	Semi-Prep	SS	10-32 coned	100	9.4	2	1000	1620
109SE0215	PolySULFOETHYL A® Column, SS, 100 x 9.4mm, 2µm,1500Å	SCX	Semi-Prep	SS	10-32 coned	100	9.4	2	1500	1620
209SE0203	PolySULFOETHYL A® Column, SS, 200 x 9.4mm, 2µm,300Å	SCX	Semi-Prep	SS	10-32 coned	200	9.4	2	300	2550
209SE0210	PolySULFOETHYL A® Column, SS, 200 x 9.4mm, 2µm,1000Å	SCX	Semi-Prep	SS	10-32 coned	200	9.4	2	1000	2550
209SE0215	PolySULFOETHYL A® Column, SS, 200 x 9.4mm, 2µm,1500Å	SCX	Semi-Prep	SS	10-32 coned	200	9.4	2	1500	2550
051SE0302	PolySULFOETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 3µm,200Å	SCX	Analytical	GLT	10-32 coned	50	1.0	3	200	570
051SE0303	PolySULFOETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 3µm,300Å	SCX	Analytical	GLT	10-32 coned	50	1.0	3	300	570
051SE0310	PolySULFOETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 3µm,1000Å	SCX	Analytical	GLT	10-32 coned	50	1.0	3	1000	570
051SE0315	PolySULFOETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 3µm,1500Å	SCX	Analytical	GLT	10-32 coned	50	1.0	3	1500	570
101SE0302	PolySULFOETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 3µm,200Å	SCX	Analytical	GLT	10-32 coned	100	1.0	3	200	660
101SE0303	PolySULFOETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 3µm,300Å	SCX	Analytical	GLT	10-32 coned	100	1.0	3	300	660
101SE0310	PolySULFOETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 3µm,1000Å	SCX	Analytical	GLT	10-32 coned	100	1.0	3	1000	660
101SE0315	PolySULFOETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 3µm,1500Å	SCX	Analytical	GLT	10-32 coned	100	1.0	3	1500	660
151SE0302	PolySULFOETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 3µm,200Å	SCX	Analytical	GLT	10-32 coned	150	1.0	3	200	740
151SE0303	PolySULFOETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 3µm,300Å	SCX	Analytical	GLT	10-32 coned	150	1.0	3	300	740
151SE0310	PolySULFOETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 3µm,1000Å	SCX	Analytical	GLT	10-32 coned	150	1.0	3	1000	740
151SE0315	PolySULFOETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 3µm,1500Å	SCX	Analytical	GLT	10-32 coned	150	1.0	3	1500	740
3.52SE0302	PolySULFOETHYL A® Column, SS, 35 x 2.1mm, 3µm,200Å	SCX	Analytical	SS	10-32 coned	35	2.1	3	200	560
3.52SE0303	PolySULFOETHYL A® Column, SS, 35 x 2.1mm, 3µm,300Å	SCX	Analytical	SS	10-32 coned	35	2.1	3	300	560
3.52SE0310	PolySULFOETHYL A® Column, SS, 35 x 2.1mm, 3µm,1000Å	SCX	Analytical	SS	10-32 coned	35	2.1	3	1000	560
3.52SE0315	PolySULFOETHYL A® Column, SS, 35 x 2.1mm, 3µm,1500Å	SCX	Analytical	SS	10-32 coned	35	2.1	3	1500	560
052SE0302	PolySULFOETHYL A® Column, SS, 50 x 2.1mm, 3µm,200Å	SCX	Analytical	SS	10-32 coned	35	2.1	3	200	585
052SE0303	PolySULFOETHYL A® Column, SS, 50 x 2.1mm, 3µm,300Å	SCX	Analytical	SS	10-32 coned	35	2.1	3	300	585
052SE0310	PolySULFOETHYL A® Column, SS, 50 x 2.1mm, 3µm,1000Å	SCX	Analytical	SS	10-32 coned	35	2.1	3	1000	585
052SE0315	PolySULFOETHYL A® Column, SS, 50 x 2.1mm, 3µm,1500Å	SCX	Analytical	SS	10-32 coned	35	2.1	3	1500	585
102SE0302	PolySULFOETHYL A® Column, SS, 100 x 2.1mm, 3µm,200Å	SCX	Analytical	SS	10-32 coned	35	2.1	3	200	750
102SE0303	PolySULFOETHYL A® Column, SS, 100 x 2.1mm, 3µm,300Å	SCX	Analytical	SS	10-32 coned	100	2.1	3	300	750
102SE0310	PolySULFOETHYL A® Column, SS, 100 x 2.1mm, 3µm,1000Å	SCX	Analytical	SS	10-32 coned	100	2.1	3	1000	750
102SE0315	PolySULFOETHYL A® Column, SS, 100 x 2.1mm, 3µm,1500Å	SCX	Analytical	SS	10-32 coned	100	2.1	3	1500	750
3.54SE0302	PolySULFOETHYL A® Column, SS, 35 x 4.6mm, 3µm,200Å	SCX	Analytical	SS	10-32 coned	35	4.6	3	200	560
3.54SE0303	PolySULFOETHYL A® Column, SS, 35 x 4.6mm, 3µm,300Å	SCX	Analytical	SS	10-32 coned	35	4.6	3	300	560
3.54SE0310	PolySULFOETHYL A® Column, SS, 35 x 4.6mm, 3µm,1000Å	SCX	Analytical	SS	10-32 coned	35	4.6	3	1000	560
3.54SE0315	PolySULFOETHYL A® Column, SS, 35 x 4.6mm, 3µm,1500Å	SCX	Analytical	SS	10-32 coned	35	4.6	3	1500	560
054SE0302	PolySULFOETHYL A® Column, SS, 50 x 4.6mm, 3µm,200Å	SCX	Analytical	SS	10-32 coned	50	4.6	3	200	585
054SE0303	PolySULFOETHYL A® Column, SS, 50 x 4.6mm, 3µm,300Å	SCX	Analytical	SS	10-32 coned	50	4.6	3	300	585
054SE0310	PolySULFOETHYL A® Column, SS, 50 x 4.6mm, 3µm,1000Å	SCX	Analytical	SS	10-32 coned	50	4.6	3	1000	585
054SE0315	PolySULFOETHYL A® Column, SS, 50 x 4.6mm, 3µm,1500Å	SCX	Analytical	SS	10-32 coned	50	4.6	3	1500	585
104SE0302	PolySULFOETHYL A® Column, SS, 100 x 4.6mm, 3µm,200Å	SCX	Analytical	SS	10-32 coned	100	4.6	3	1500	750
104SE0303	PolySULFOETHYL A® Column, SS, 100 x 4.6mm, 3µm,300Å	SCX	Analytical	SS	10-32 coned	100	4.6	3	300	750
104SE0310	PolySULFOETHYL A® Column, SS, 100 x 4.6mm, 3µm,1000Å	SCX	Analytical	SS	10-32 coned	100	4.6	3	1000	750
104SE0315	PolySULFOETHYL A® Column, SS, 100 x 4.6mm, 3µm,1500Å	SCX	Analytical	SS	10-32 coned	100	4.6	3	1500	750
154SE0302	PolySULFOETHYL A® Column, SS, 150 x 4.6mm, 3µm,200Å	SCX	Analytical	SS	10-32 coned	150	4.6	3	200	810
154SE0303	PolySULFOETHYL A® Column, SS, 150 x 4.6mm, 3µm,300Å	SCX	Analytical	SS	10-32 coned	150	4.6	3	300	810
154SE0310	PolySULFOETHYL A® Column, SS, 150 x 4.6mm, 3µm,1000Å	SCX	Analytical	SS	10-32 coned	150	4.6	3	1000	810

PolySULFOETHYL A® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.) Custom combinations available upon request info@polylc.com

PN	ProduSE Description	Phase	Type	Hardware	ConneSEion	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
154SE0315	PolySULFOETHYL A® Column, SS, 150 x 4.6mm, 3µm,1500Å	SCX	Analytical	SS	10-32 coned	150	4.6	3	1500	810
204SE0302	PolySULFOETHYL A® Column, SS, 200 x 4.6mm, 3µm,200Å	SCX	Analytical	SS	10-32 coned	200	4.6	3	200	870
204SE0303	PolySULFOETHYL A® Column, SS, 200 x 4.6mm, 3µm,300Å	SCX	Analytical	SS	10-32 coned	200	4.6	3	300	870
204SE0310	PolySULFOETHYL A® Column, SS, 200 x 4.6mm, 3µm,1000Å	SCX	Analytical	SS	10-32 coned	200	4.6	3	1000	870
204SE0315	PolySULFOETHYL A® Column, SS, 200 x 4.6mm, 3µm,1500Å	SCX	Analytical	SS	10-32 coned	200	4.6	3	1500	870
254SE0302	PolySULFOETHYL A® Column, SS, 250 x 4.6mm, 3µm,200Å	SCX	Analytical	SS	10-32 coned	250	4.6	3	200	930
254SE0303	PolySULFOETHYL A® Column, SS, 250 x 4.6mm, 3µm,300Å	SCX	Analytical	SS	10-32 coned	250	4.6	3	300	930
254SE0310	PolySULFOETHYL A® Column, SS, 250 x 4.6mm, 3µm,1000Å	SCX	Analytical	SS	10-32 coned	250	4.6	3	1000	930
254SE0315	PolySULFOETHYL A® Column, SS, 250 x 4.6mm, 3µm,1500Å	SCX	Analytical	SS	10-32 coned	250	4.6	3	1500	930
109SE0302	PolySULFOETHYL A® Column, SS, 100 x 9.4mm, 3µm,200Å	SCX	Analytical	SS	10-32 coned	100	9.4	3	200	1490
109SE0303	PolySULFOETHYL A® Column, SS, 100 x 9.4mm, 3µm,300Å	SCX	Analytical	SS	10-32 coned	100	9.4	3	300	1490
109SE0310	PolySULFOETHYL A® Column, SS, 100 x 9.4mm, 3µm,1000Å	SCX	Analytical	SS	10-32 coned	100	9.4	3	1000	1490
109SE0315	PolySULFOETHYL A® Column, SS, 100 x 9.4mm, 3µm,1500Å	SCX	Analytical	SS	10-32 coned	100	9.4	3	1500	1490
209SE0302	PolySULFOETHYL A® Column, SS, 200 x 9.4mm, 3µm,200Å	SCX	Analytical	SS	10-32 coned	200	9.4	3	200	1700
209SE0303	PolySULFOETHYL A® Column, SS, 200 x 9.4mm, 3µm,300Å	SCX	Analytical	SS	10-32 coned	200	9.4	3	300	1700
209SE0310	PolySULFOETHYL A® Column, SS, 200 x 9.4mm, 3µm,1000Å	SCX	Analytical	SS	10-32 coned	200	9.4	3	1000	1700
209SE0315	PolySULFOETHYL A® Column, SS, 200 x 9.4mm, 3µm,1500Å	SCX	Analytical	SS	10-32 coned	200	9.4	3	1500	1700
051SE0502	PolySULFOETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 5µm,200Å	SCX	Analytical	GLT	10-32 coned	50	1.0	5	200	570
051SE0503	PolySULFOETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 5µm,300Å	SCX	Analytical	GLT	10-32 coned	50	1.0	5	300	570
051SE0510	PolySULFOETHYL A® Column, Glass Lined SS, 50 x 1.0mm, 5µm,1000Å	SCX	Analytical	GLT	10-32 coned	50	1.0	5	1000	570
101SE0502	PolySULFOETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 5µm,200Å	SCX	Analytical	GLT	10-32 coned	100	1.0	5	200	660
101SE0503	PolySULFOETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 5µm,300Å	SCX	Analytical	GLT	10-32 coned	100	1.0	5	300	660
101SE0510	PolySULFOETHYL A® Column, Glass Lined SS, 100 x 1.0mm, 5µm,1000Å	SCX	Analytical	GLT	10-32 coned	100	1.0	5	1000	660
151SE0502	PolySULFOETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 5µm,200Å	SCX	Analytical	GLT	10-32 coned	150	1.0	5	200	740
151SE0503	PolySULFOETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 5µm,300Å	SCX	Analytical	GLT	10-32 coned	150	1.0	5	300	740
151SE0510	PolySULFOETHYL A® Column, Glass Lined SS, 150 x 1.0mm, 5µm,1000Å	SCX	Analytical	GLT	10-32 coned	150	1.0	5	1000	740
3.52SE0502	PolySULFOETHYL A® Column, SS, 35 x 2.1mm, 5µm,200Å	SCX	Analytical	SS	10-32 coned	35	2.1	5	200	510
3.52SE0503	PolySULFOETHYL A® Column, SS, 35 x 2.1mm, 5µm,300Å	SCX	Analytical	SS	10-32 coned	35	2.1	5	300	510
3.52SE0510	PolySULFOETHYL A® Column, SS, 35 x 2.1mm, 5µm,1000Å	SCX	Analytical	SS	10-32 coned	35	2.1	5	1000	510
052SE0502	PolySULFOETHYL A® Column, SS, 50 x 2.1mm, 5µm,200Å	SCX	Analytical	SS	10-32 coned	50	2.1	5	200	540
052SE0503	PolySULFOETHYL A® Column, SS, 50 x 2.1mm, 5µm,300Å	SCX	Analytical	SS	10-32 coned	50	2.1	5	300	540
052SE0510	PolySULFOETHYL A® Column, SS, 50 x 2.1mm, 5µm,1000Å	SCX	Analytical	SS	10-32 coned	50	2.1	5	1000	540
102SE0502	PolySULFOETHYL A® Column, SS, 100 x 2.1mm, 5µm,200Å	SCX	Analytical	SS	10-32 coned	100	2.1	5	200	650
102SE0503	PolySULFOETHYL A® Column, SS, 100 x 2.1mm, 5µm,300Å	SCX	Analytical	SS	10-32 coned	100	2.1	5	300	650
102SE0510	PolySULFOETHYL A® Column, SS, 100 x 2.1mm, 5µm,1000Å	SCX	Analytical	SS	10-32 coned	100	2.1	5	1000	650
152SE0502	PolySULFOETHYL A® Column, SS, 150 x 2.1mm, 5µm,200Å	SCX	Analytical	SS	10-32 coned	150	2.1	5	200	700
152SE0503	PolySULFOETHYL A® Column, SS, 150 x 2.1mm, 5µm,300Å	SCX	Analytical	SS	10-32 coned	150	2.1	5	300	700
152SE0510	PolySULFOETHYL A® Column, SS, 150 x 2.1mm, 5µm,1000Å	SCX	Analytical	SS	10-32 coned	150	2.1	5	1000	700
202SE0502	PolySULFOETHYL A® Column, SS, 150 x 2.1mm, 5µm,200Å	SCX	Analytical	SS	10-32 coned	200	2.1	5	200	750
202SE0503	PolySULFOETHYL A® Column, SS, 200 x 2.1mm, 5µm,300Å	SCX	Analytical	SS	10-32 coned	200	2.1	5	300	750
202SE0510	PolySULFOETHYL A® Column, SS, 200 x 2.1mm, 5µm,1000Å	SCX	Analytical	SS	10-32 coned	200	2.1	5	1000	750
3.54SE0502	PolySULFOETHYL A® Column, SS, 35 x 4.6mm, 5µm,200Å	SCX	Analytical	SS	10-32 coned	35	4.6	5	200	510
3.54SE0503	PolySULFOETHYL A® Column, SS, 35 x 4.6mm, 5µm,300Å	SCX	Analytical	SS	10-32 coned	35	4.6	5	300	510
3.54SE0510	PolySULFOETHYL A® Column, SS, 35 x 4.6mm, 5µm,1000Å	SCX	Analytical	SS	10-32 coned	35	4.6	5	1000	510
054SE0502	PolySULFOETHYL A® Column, SS, 50 x 4.6mm, 5µm,200Å	SCX	Analytical	SS	10-32 coned	50	4.6	5	200	540

PolySULFOETHYL A® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes. Custom combinations available upon request info@polylc.com)

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
054SE0502	PolySULFOETHYL A® Column, SS, 50 x 4.6mm, 5µm,200Å	SCX	Analytical	SS	10-32 coned	50	4.6	5	200	540
054SE0503	PolySULFOETHYL A® Column, SS, 50 x 4.6mm, 5µm,300Å	SCX	Analytical	SS	10-32 coned	50	4.6	5	300	540
054SE0510	PolySULFOETHYL A® Column, SS, 50 x 4.6mm, 5µm,1000Å	SCX	Analytical	SS	10-32 coned	50	4.6	5	1000	540
104SE0502	PolySULFOETHYL A® Column, SS, 100 x 4.6mm, 5µm,200Å	SCX	Analytical	SS	10-32 coned	100	4.6	5	200	650
104SE0503	PolySULFOETHYL A® Column, SS, 100 x 4.6mm, 5µm,300Å	SCX	Analytical	SS	10-32 coned	100	4.6	5	300	650
104SE0510	PolySULFOETHYL A® Column, SS, 100 x 4.6mm, 5µm,1000Å	SCX	Analytical	SS	10-32 coned	100	4.6	5	1000	650
154SE0502	PolySULFOETHYL A® Column, SS, 150 x 4.6mm, 5µm,200Å	SCX	Analytical	SS	10-32 coned	150	4.6	5	200	725
154SE0503	PolySULFOETHYL A® Column, SS, 150 x 4.6mm, 5µm,300Å	SCX	Analytical	SS	10-32 coned	150	4.6	5	300	725
154SE0510	PolySULFOETHYL A® Column, SS, 150 x 4.6mm, 5µm,1000Å	SCX	Analytical	SS	10-32 coned	150	4.6	5	1000	725
204SE0502	PolySULFOETHYL A® Column, SS, 200 x 4.6mm, 5µm,200Å	SCX	Analytical	SS	10-32 coned	200	4.6	5	200	750
204SE0503	PolySULFOETHYL A® Column, SS, 200 x 4.6mm, 5µm,300Å	SCX	Analytical	SS	10-32 coned	200	4.6	5	300	750
204SE0510	PolySULFOETHYL A® Column, SS, 200 x 4.6mm, 5µm,1000Å	SCX	Analytical	SS	10-32 coned	200	4.6	5	1000	750
254SE0502	PolySULFOETHYL A® Column, SS, 250 x 4.6mm, 5µm,200Å	SCX	Analytical	SS	10-32 coned	250	4.6	5	200	775
254SE0503	PolySULFOETHYL A® Column, SS, 250 x 4.6mm, 5µm,300Å	SCX	Analytical	SS	10-32 coned	250	4.6	5	300	775
254SE0510	PolySULFOETHYL A® Column, SS, 250 x 4.6mm, 5µm,1000Å	SCX	Analytical	SS	10-32 coned	250	4.6	5	1000	775
109SE0502	PolySULFOETHYL A® Column, SS, 100 x 9.4mm, 5µm,200Å	SCX	Semi-Prep	SS	10-32 coned	100	9.4	5	200	1250
109SE0503	PolySULFOETHYL A® Column, SS, 100 x 9.4mm, 5µm,300Å	SCX	Semi-Prep	SS	10-32 coned	100	9.4	5	300	1250
109SE0510	PolySULFOETHYL A® Column, SS, 100 x 9.4mm, 5µm,1000Å	SCX	Semi-Prep	SS	10-32 coned	100	9.4	5	1000	1250
209SE0502	PolySULFOETHYL A® Column, SS, 200 x 9.4mm, 5µm,200Å	SCX	Semi-Prep	SS	10-32 coned	200	9.4	5	200	1500
209SE0503	PolySULFOETHYL A® Column, SS, 200 x 9.4mm, 5µm,300Å	SCX	Semi-Prep	SS	10-32 coned	200	9.4	5	300	1500
209SE0510	PolySULFOETHYL A® Column, SS, 200 x 9.4mm, 5µm,1000Å	SCX	Semi-Prep	SS	10-32 coned	200	9.4	5	1000	1500
259SE0502	PolySULFOETHYL A® Column, SS, 250 x 9.4mm, 5µm,200Å	SCX	Semi-Prep	SS	10-32 coned	250	9.4	5	200	1570
259SE0503	PolySULFOETHYL A® Column, SS, 250 x 9.4mm, 5µm,300Å	SCX	Semi-Prep	SS	10-32 coned	250	9.4	5	300	1570
259SE0510	PolySULFOETHYL A® Column, SS, 250 x 9.4mm, 5µm,1000Å	SCX	Semi-Prep	SS	10-32 coned	250	9.4	5	1000	1570
2521SE0502	PolySULFOETHYL A® Column, SS, 250 x 21mm, 5µm,200Å	SCX	Semi-Prep	SS	10-32 coned	250	21.0	5	200	4700
2521SE0503	PolySULFOETHYL A® Column, SS, 250 x 21mm, 5µm,300Å	SCX	Semi-Prep	SS	10-32 coned	250	21.0	5	300	4700
2521SE0510	PolySULFOETHYL A® Column, SS, 250 x 21mm, 5µm,1000Å	SCX	Semi-Prep	SS	10-32 coned	250	21.0	5	1000	4700
204SE1203	PolySULFOETHYL A® Column, SS, 200 x 4.6mm, 12µm,300Å	SCX	Analytical	SS	10-32 coned	200	4.6	12	300	430
204SE1215	PolySULFOETHYL A® Column, SS, 200 x 4.6mm, 12µm,1500Å	SCX	Analytical	SS	10-32 coned	200	4.6	12	1500	430
254SE1203	PolySULFOETHYL A® Column, SS, 250 x 4.6mm, 12µm,300Å	SCX	Analytical	SS	10-32 coned	250	4.6	12	300	445
254SE1215	PolySULFOETHYL A® Column, SS, 250 x 4.6mm, 12µm,1500Å	SCX	Analytical	SS	10-32 coned	250	4.6	12	1500	445
209SE1203	PolySULFOETHYL A® Column, SS, 200 x 9.4mm, 12µm,300Å	SCX	Semi-Prep	SS	10-32 coned	200	9.4	12	300	1060
209SE1215	PolySULFOETHYL A® Column, SS, 200 x 9.4mm, 12µm,1500Å	SCX	Semi-Prep	SS	10-32 coned	200	9.4	12	1500	1060
259SE1203	PolySULFOETHYL A® Column, SS, 250 x 9.4mm, 12µm,300Å	SCX	Semi-Prep	SS	10-32 coned	250	9.4	12	300	1160
259SE1215	PolySULFOETHYL A® Column, SS, 250 x 9.4mm, 12µm,1500Å	SCX	Semi-Prep	SS	10-32 coned	250	9.4	12	1500	1160
2521SE1203	PolySULFOETHYL A® Column, SS, 250 x 21mm, 12µm,300Å	SCX	Semi-Prep	SS	10-32 coned	250	21.0	12	300	4000
2521SE1215	PolySULFOETHYL A® Column, SS, 250 x 21mm, 12µm,1500Å	SCX	Semi-Prep	SS	10-32 coned	250	21.0	12	1500	4000
2550SE1203	PolySULFOETHYL A® Column, SS, 250 x 50.8mm, 12µm,300Å	SCX	Prep	SS	10-32 coned	250	50.8	12	300	14150
2550SE1215	PolySULFOETHYL A® Column, SS, 250 x 50.8mm, 12µm,1500Å	SCX	Prep	SS	10-32 coned	250	50.8	12	1500	14150

HPLC Capillaries



PolySULFOETHYL A® Capillaries (Available in 2, 3, 5, and 5µm particles with multiple pore sizes.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
050.075SE0203	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 2µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.075	2	300	520
050.075SE0210	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 2µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.075	2	1000	520
050.075SE0215	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 2µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.075	2	1500	520
100.075SE0203	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 2µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.075	2	300	580
100.075SE0210	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 2µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.075	2	1000	580
100.075SE0215	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 2µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.075	2	1500	580
150.075SE0203	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 2µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.075	2	300	640
150.075SE0210	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 2µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.075	2	1000	640
150.075SE0215	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 2µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.075	2	1500	640
200.075SE0203	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 2µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.075	2	300	700
200.075SE0210	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 2µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.075	2	1000	700
200.075SE0215	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 2µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.075	2	1500	700
050.15SE0203	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 2µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.150	2	300	520
050.15SE0210	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 2µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.150	2	1000	520
050.15SE0215	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 2µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.150	2	1500	520
100.15SE0203	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 2µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.150	2	300	580
100.15SE0210	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 2µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.150	2	1000	580
100.15SE0215	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 2µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.150	2	1500	580
150.15SE0203	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 2µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.150	2	300	640
150.15SE0210	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 2µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.150	2	1000	640
150.15SE0215	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 2µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.150	2	1500	640
200.15SE0203	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 2µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.150	2	300	700
200.15SE0210	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 2µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.150	2	1000	700
200.15SE0215	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 2µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.150	2	1500	700
050.30SE0203	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 2µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.300	2	300	520
050.30SE0210	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 2µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.300	2	1000	520
050.30SE0215	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 2µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.300	2	1500	520
100.30SE0203	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 2µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.300	2	300	580
100.30SE0210	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 2µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.300	2	1000	580
100.30SE0215	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 2µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.300	2	1500	580
150.30SE0203	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 2µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.300	2	300	640
150.30SE0210	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 2µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.300	2	1000	640
150.30SE0215	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 2µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.300	2	1500	640
200.30SE0203	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 2µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.300	2	300	700
200.30SE0210	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 2µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.300	2	1000	700
200.30SE0215	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 2µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.300	2	1500	700

PolySULFOETHYL A® Capillaries (Available in 2, 3, 5, and 5µm particles with multiple pore sizes.) Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
050.075SE0302	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.075	3	200	520
050.075SE0303	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.075	3	300	520
050.075SE0310	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.075	3	1000	520
050.075SE0315	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.075	3	1500	520
100.075SE0302	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.075	3	200	580
100.075SE0303	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.075	3	300	580
100.075SE0310	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.075	3	1000	580
100.075SE0315	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.075	3	1500	580
150.075SE0302	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.075	3	200	640
150.075SE0303	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.075	3	300	640
150.075SE0310	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.075	3	1000	640
150.075SE0315	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.075	3	1500	640
200.075SE0302	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.075	3	200	700
200.075SE0303	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.075	3	300	700
200.075SE0310	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.075	3	1000	700
200.075SE0315	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.075	3	1500	700
050.15SE0302	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.150	3	200	520
050.15SE0303	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.150	3	300	520
050.15SE0310	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.150	3	1000	520
050.15SE0315	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.150	3	1500	520
100.15SE0302	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.150	3	200	580
100.15SE0303	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.150	3	300	580
100.15SE0310	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.150	3	1000	580
100.15SE0315	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.150	3	1500	580
150.15SE0302	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.150	3	200	640
150.15SE0303	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.150	3	300	640
150.15SE0310	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.150	3	1000	640
150.15SE0315	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.150	3	1500	640
200.15SE0302	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.150	3	200	700
200.15SE0303	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.150	3	300	700
200.15SE0310	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.150	3	1000	700
200.15SE0315	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.150	3	1500	700
050.30SE0302	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.300	3	200	520
050.30SE0303	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.300	3	300	520
050.30SE0310	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.300	3	1000	520
050.30SE0315	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.300	3	1500	520
100.30SE0302	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.300	3	200	580
100.30SE0303	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.300	3	300	580
100.30SE0310	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.300	3	1000	580
100.30SE0315	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.300	3	1500	580
150.30SE0302	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.300	3	200	640
150.30SE0303	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.300	3	300	640
150.30SE0310	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.300	3	1000	640
150.30SE0315	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.300	3	1500	640
200.30SE0302	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.300	3	200	700
200.30SE0303	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.300	3	300	700

PolySULFOETHYL A® Capillaries (Available in 2, 3, 5, and 5µm particles with multiple pore sizes.) Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
200.30SE0310	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.300	3	1000	700
200.30SE0315	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,1500Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.300	3	1500	700
050.075SE0502	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 5µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.075	5	200	520
050.075SE0503	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 5µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.075	5	300	520
050.075SE0510	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 5µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.075	5	1000	520
100.075SE0502	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 5µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.075	5	200	580
100.075SE0503	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 5µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.075	5	300	580
100.075SE0510	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 5µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.075	5	1000	580
150.075SE0502	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 5µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.075	5	200	640
150.075SE0503	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 5µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.075	5	300	640
150.075SE0510	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 5µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.075	5	1000	640
200.075SE0502	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 5µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.075	5	200	700
200.075SE0503	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 5µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.075	5	300	700
200.075SE0510	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 5µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.075	5	1000	700
050.15SE0502	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 5µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.150	5	200	520
050.15SE0503	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 5µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.150	5	300	520
050.15SE0510	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 5µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.150	5	1000	520
100.15SE0502	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 5µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.150	5	200	580
100.15SE0503	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 5µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.150	5	300	580
100.15SE0510	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 5µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.150	5	1000	580
150.15SE0502	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 5µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.150	5	200	640
150.15SE0503	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 5µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.150	5	300	640
150.15SE0510	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 5µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.150	5	1000	640
200.15SE0502	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 5µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.150	5	200	700
200.15SE0503	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 5µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.150	5	300	700
200.15SE0510	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 5µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.150	5	1000	700
050.30SE0502	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 5µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.300	5	200	520
050.30SE0503	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 5µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.300	5	300	520
050.30SE0510	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 5µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	50	0.300	5	1000	520
100.30SE0502	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 5µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.300	5	200	580
100.30SE0503	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 5µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.300	5	300	580
100.30SE0510	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 5µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	100	0.300	5	1000	580
150.30SE0502	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 5µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.300	5	200	640
150.30SE0503	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 5µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.300	5	300	640
150.30SE0510	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 5µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	150	0.300	5	1000	640
200.30SE0502	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 5µm,200Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.300	5	200	700
200.30SE0503	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 5µm,300Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.300	5	300	700
200.30SE0510	PolySULFOETHYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 5µm,1000Å	SCX	Analytical	PEEK SIL	10-32 coned	200	0.300	5	1000	700

Guard Cartridges



PolySULFOETHYL A® Guard Cartridges (Available in 2, 3, and 5µm particles with pore sizes.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Conn	SE	Ion L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
GC1-SE0203	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 2µm, 300Å	SCX	Analytical	SS/PEEK	10-32	coned	10	1.0	2	300	120
GC1-SE0210	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 2µm, 1000Å	SCX	Analytical	SS/PEEK	10-32	coned	10	1.0	2	1000	120
GC1-SE0215	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 2µm, 1500Å	SCX	Analytical	SS/PEEK	10-32	coned	10	1.0	2	1500	120
GC2-SE0203	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 2µm, 300Å	SCX	Analytical	SS/PEEK	10-32	coned	10	2.0	2	300	120
GC2-SE0210	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 2µm, 1000Å	SCX	Analytical	SS/PEEK	10-32	coned	10	2.0	2	1000	120
GC2-SE0215	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 2µm, 1500Å	SCX	Analytical	SS/PEEK	10-32	coned	10	2.0	2	1500	120
GC3-SE0203	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 2µm, 300Å	SCX	Analytical	SS/PEEK	10-32	coned	10	3.0	2	300	120
GC3-SE0210	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 2µm, 1000Å	SCX	Analytical	SS/PEEK	10-32	coned	10	3.0	2	1000	120
GC3-SE0215	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 2µm, 1500Å	SCX	Analytical	SS/PEEK	10-32	coned	10	3.0	2	1500	120
GC4-SE0203	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 2µm, 300Å	SCX	Analytical	SS/PEEK	10-32	coned	10	4.0	2	300	120
GC4-SE0210	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 2µm, 1000Å	SCX	Analytical	SS/PEEK	10-32	coned	10	4.0	2	1000	120
GC4-SE0215	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 2µm, 1500Å	SCX	Analytical	SS/PEEK	10-32	coned	10	4.0	2	1500	120
GC1-SE0302	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 200Å	SCX	Analytical	SS/PEEK	10-32	coned	10	1.0	3	200	100
GC1-SE0303	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 300Å	SCX	Analytical	SS/PEEK	10-32	coned	10	1.0	3	300	100
GC1-SE0310	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 1000Å	SCX	Analytical	SS/PEEK	10-32	coned	10	1.0	3	1000	100
GC1-SE0315	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 1500Å	SCX	Analytical	SS/PEEK	10-32	coned	10	1.0	3	1500	100
GC2-SE0302	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 3µm, 200Å	SCX	Analytical	SS/PEEK	10-32	coned	10	2.0	3	200	100
GC2-SE0303	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 3µm, 300Å	SCX	Analytical	SS/PEEK	10-32	coned	10	2.0	3	300	100
GC2-SE0310	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 3µm, 1000Å	SCX	Analytical	SS/PEEK	10-32	coned	10	2.0	3	1000	100
GC2-SE0315	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 3µm, 1500Å	SCX	Analytical	SS/PEEK	10-32	coned	10	2.0	3	1500	100
GC3-SE0302	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 3µm, 200Å	SCX	Analytical	SS/PEEK	10-32	coned	10	3.0	3	200	100
GC3-SE0303	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 3µm, 300Å	SCX	Analytical	SS/PEEK	10-32	coned	10	3.0	3	300	100
GC3-SE0310	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 3µm, 1000Å	SCX	Analytical	SS/PEEK	10-32	coned	10	3.0	3	1000	100
GC3-SE0315	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 3µm, 1500Å	SCX	Analytical	SS/PEEK	10-32	coned	10	3.0	3	1500	100
GC4-SE0302	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 3µm, 200Å	SCX	Analytical	SS/PEEK	10-32	coned	10	4.0	3	200	100
GC4-SE0303	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 3µm, 300Å	SCX	Analytical	SS/PEEK	10-32	coned	10	4.0	3	300	100
GC4-SE0310	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 3µm, 1000Å	SCX	Analytical	SS/PEEK	10-32	coned	10	4.0	3	1000	100
GC4-SE0315	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 3µm, 1500Å	SCX	Analytical	SS/PEEK	10-32	coned	10	4.0	3	1500	100
GC1-SE0502	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 5µm, 200Å	SCX	Analytical	SS/PEEK	10-32	coned	10	1.0	5	200	90
GC1-SE0503	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 5µm, 300Å	SCX	Analytical	SS/PEEK	10-32	coned	10	1.0	5	300	90
GC1-SE0510	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 5µm, 1000Å	SCX	Analytical	SS/PEEK	10-32	coned	10	1.0	5	1000	90
GC2-SE0502	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 5µm, 200Å	SCX	Analytical	SS/PEEK	10-32	coned	10	2.0	5	200	90
GC2-SE0503	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 5µm, 300Å	SCX	Analytical	SS/PEEK	10-32	coned	10	2.0	5	300	90
GC2-SE0510	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 5µm, 1000Å	SCX	Analytical	SS/PEEK	10-32	coned	10	2.0	5	1000	90
GC3-SE0502	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 5µm, 200Å	SCX	Analytical	SS/PEEK	10-32	coned	10	3.0	5	200	90

PolySULFOETHYL A® Guard Cartridges (Available in 2, 3, and 5µm particles with pore sizes.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	ConneSEion	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
GC3-SE0503	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 5µm, 300Å	SCX	Analytical	SS/PEEK	10-32 coned	10	3.0	5	300	90
GC3-SE0510	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 5µm, 1000Å	SCX	Analytical	SS/PEEK	10-32 coned	10	3.0	5	1000	90
GC4-SE0502	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 5µm, 200Å	SCX	Analytical	SS/PEEK	10-32 coned	10	4.0	5	200	90
GC4-SE0503	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 5µm, 300Å	SCX	Analytical	SS/PEEK	10-32 coned	10	4.0	5	300	90
GC4-SE0510	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 5µm, 1000Å	SCX	Analytical	SS/PEEK	10-32 coned	10	4.0	5	1000	90
DGCH10	Direct Guard Cartridge Holder, 10mm (For all analytical guards, REUSABLE)	N/A	Analytical	SS/PEEK	10-32 coned	10	N/A	N/A	N/A	75
IGCH10	Inline Guard Cartridge Holder, 10mm (For all analytical guards, REUSABLE)	N/A	Analytical	SS/PEEK	10-32 coned	10	N/A	N/A	N/A	75
*** NOTE: DGCH10 for Direct Connection to Column (female inlet/male outlet) IGCH10 for Connection to Line (two female ports)										
GC10-SE0302	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 200Å	SCX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	200	225
GC10-SE0303	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 300Å	SCX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	300	225
GC10-SE0310	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 1000Å	SCX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	1000	225
GC10-SE0315	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 1500Å	SCX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	1500	225
GC10-SE0502	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 5µm, 200Å	SCX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	5	200	185
GC10-SE0503	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 5µm, 300Å	SCX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	5	300	185
GC10-SE0510	PolySULFOETHYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 5µm, 1000Å	SCX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	5	1000	185
CH10	Semi-Prep Guard Cartridge Holder, 10mm (For 9.4mm columns, REUSABLE)	SCX	Semi-Prep	SS/PEEK	10-32 coned	10	N/A	N/A	N/A	175

*** NOTE: CH10 Holder Required for Semi-Prep Guard Cartridges (10 x 10mm) Connection to Line (two female ports)

Solid Phase Extraction Cartridges


PolySULFOETHYL A® SPE Cartridges

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	ConneSEion	V (ml)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
SPESE1203	PolySULFOETHYL A® SPE Cartridge 1mL, PE, 12µm, 300Å (Pack of 10 cartridges)	SCX	SPE	PE	Lure	1	N/A	12	300	130
SPESE1215	PolySULFOETHYL A® SPE Cartridge 1mL, PE, 12µm, 1500Å (Pack of 10 cartridges)	SCX	SPE	PE	Lure	1	N/A	12	1500	130
SPESE2003	PolySULFOETHYL A® SPE Cartridge 1mL, PE, 20µm, 300Å (Pack of 10 cartridges)	SCX	SPE	PE	Lure	1	N/A	20	300	130
SPESE2015	PolySULFOETHYL A® SPE Cartridge 1mL, PE, 20µm, 1500Å (Pack of 10 cartridges)	SCX	SPE	PE	Lure	1	N/A	20	1500	130



Spin Columns



PolySULFOETHYL A® BioPureSPN Columns

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Bed (mg)	Pore (Å)	SCX Capacity	HILIC Capacity	2026 Price (USD)
HMM HIL-SCX	PolySULFOETHYL A® BioPureSPN MACRO for SCX or HILIC 50-200µl (Box of 100 each)	SCX	Spin	100	300	0.2 -1mg	35 - 350µg	940
HMM HIL-SCX.20	PolySULFOETHYL A® BioPureSPN MACRO for SCX or HILIC 50-200µl (Pack of 20 each)	SCX	Spin	100	300	0.2 -1mg	35 - 350µg	265
HMM HIL-SCX-L	PolySULFOETHYL A® BioPureSPN MACRO PLATE 96 Well for SCX or HILIC 50-200µl	SCX	Spin	100	300	0.2 -1mg	35 - 350µg	800
HEM HIL-SCX	PolySULFOETHYL A® BioPureSPN MIDI for SCX/HILIC 50-200µl (Box of 100 each)	SCX	Spin	50	300	100-500µg	17 - 170µg	820
HEM HIL-SCX.20	PolySULFOETHYL A® BioPureSPN MIDI for SCX/HILIC 50-200µl (Pack of 20 each)	SCX	Spin	50	300	100-500µg	17 - 170µg	240
HNS HIL-SCX-M	PolySULFOETHYL A® BioPureSPN MIDI PLATE 96 Well for SCX/HILIC 50-200µl	SCX	Spin	50	300	100-500µg	17 - 170µg	575
HUM HIL-SCX	PolySULFOETHYL A® BioPureSPN MINI for SCX/HILIC 20-100µl (Box of 100 each)	SCX	Spin	20	300	40-200µg	7 - 70µg	735
HUM HIL-SCX.20	PolySULFOETHYL A® BioPureSPN MINI for SCX/HILIC 20-100µl (Pack of 20 each)	SCX	Spin	20	300	40-200µg	7 - 70µg	225
HNS HIL-SCX	PolySULFOETHYL A® BioPureSPN MINI PLATE 96 Well for SCX/HILIC 20-100µl	SCX	Spin	20	300	40-200µg	7 - 70µg	575

Bulk Material



PolySULFOETHYL A® Bulk Material (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore size options.)

Special request? Contact info@polylc.com

PN	Product Description	Phase	Type	Hardware	QTY (g)	Particle (µm)	Pore (Å)	2026 Price (USD)
BMSE0203	PolySULFOETHYL A® Bulk Material, 2µm, 300Å	SCX	BM	N/A	0.5	2	300	110
BMSE0210	PolySULFOETHYL A® Bulk Material, 2µm, 1000Å	SCX	BM	N/A	0.5	2	1000	110
BMSE0215	PolySULFOETHYL A® Bulk Material, 2µm, 1500Å	SCX	BM	N/A	0.5	2	1500	110
BMSE0302	PolySULFOETHYL A® Bulk Material, 3µm, 200Å	SCX	BM	N/A	1	3	200	145
BMSE0303	PolySULFOETHYL A® Bulk Material, 3µm, 300Å	SCX	BM	N/A	1	3	300	145
BMSE0310	PolySULFOETHYL A® Bulk Material, 3µm, 1000Å	SCX	BM	N/A	1	3	1000	145
BMSE0315	PolySULFOETHYL A® Bulk Material, 3µm, 1500Å	SCX	BM	N/A	1	3	1500	145
BMSE0502	PolySULFOETHYL A® Bulk Material, 5µm, 200Å	SCX	BM	N/A	1	5	200	95
BMSE0503	PolySULFOETHYL A® Bulk Material, 5µm, 300Å	SCX	BM	N/A	1	5	300	95
BMSE0510	PolySULFOETHYL A® Bulk Material, 5µm, 1000Å	SCX	BM	N/A	1	5	1000	95
BMSE1203	PolySULFOETHYL A® Bulk Material, 12µm, 300Å	SCX	BM	N/A	1	12	300	42
BMSE1215	PolySULFOETHYL A® Bulk Material, 12µm, 1500Å	SCX	BM	N/A	1	12	1500	42
BMSE2003	PolySULFOETHYL A® Bulk Material, 20µm, 300Å	SCX	BM	N/A	1	20	300	42
BMSE2015	PolySULFOETHYL A® Bulk Material, 20µm, 1500Å	SCX	BM	N/A	1	20	1500	42
BMSE3010	PolySULFOETHYL A® Bulk Material, 30µm, 1000Å	SCX	BM	N/A	1	30	1000	42

Please contact info@polylc.com with any questions or for information on any product combination not listed in the table above

PolySULFOETHYL A[®]

(Strong-Cation Exchange)

2026 PolyLC Inc. Retail Pricelist

PolyLC Inc., 9151 Rumsey Rd., Suite 180, Columbia, MD 21045 410-992-5400 www.polylc.com



2026 PolyLC Inc. Retail Pricelist

PolyWAX LP®

(Weak-Anion Exchange)



Offering high-performance liquid chromatography media, columns, capillaries, cartridges, and instrumentation. We specialize in advanced alternatives to reversed-phase for proteins, peptides, and other biomolecules, supported by deep application expertise and customized solutions for the most demanding separations.

HPLC Columns

PolyWAX LP® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
051WX0210	PolyWAX LP® Column, Glass Lined SS, 50 x 1.0mm, 2µm,1000Å	WAX	Analytical	GLT	10-32 coned	50	1.0	2	1000	615
101WX0210	PolyWAX LP® Column, Glass Lined SS, 100 x 1.0mm, 2µm,1000Å	WAX	Analytical	GLT	10-32 coned	100	1.0	2	1000	750
151WX0210	PolyWAX LP® Column, Glass Lined SS, 150 x 1.0mm, 2µm,1000Å	WAX	Analytical	GLT	10-32 coned	150	1.0	2	1000	800
3.52WX0210	PolyWAX LP® Column, SS, 35 x 2.1mm, 2µm,1000Å	WAX	Analytical	SS	10-32 coned	35	2.1	2	1000	570
052WX0210	PolyWAX LP® Column, SS, 50 x 2.1mm, 2µm,1000Å	WAX	Analytical	SS	10-32 coned	50	2.1	2	1000	580
102WX0210	PolyWAX LP® Column, SS, 100 x 2.1mm, 2µm,1000Å	WAX	Analytical	SS	10-32 coned	100	2.1	2	1000	750
3.54WX0210	PolyWAX LP® Column, SS, 35 x 4.6mm, 2µm,1000Å	WAX	Analytical	SS	10-32 coned	35	4.6	2	1000	570
054WX0210	PolyWAX LP® Column, SS, 50 x 4.6mm, 2µm,1000Å	WAX	Analytical	SS	10-32 coned	50	4.6	2	1000	580
104WX0210	PolyWAX LP® Column, SS, 100 x 4.6mm, 2µm,1000Å	WAX	Analytical	SS	10-32 coned	100	4.6	2	1000	750
154WX0210	PolyWAX LP® Column, SS, 150 x 4.6mm, 2µm,1000Å	WAX	Analytical	SS	10-32 coned	150	4.6	2	1000	790
204WX0210	PolyWAX LP® Column, SS, 200 x 4.6mm, 2µm,1000Å	WAX	Analytical	SS	10-32 coned	200	4.6	2	1000	810
109WX0210	PolyWAX LP® Column, SS, 100 x 9.4mm, 2µm,1000Å	WAX	Semi-Prep	SS	10-32 coned	100	9.4	2	1000	1525
209WX0210	PolyWAX LP® Column, SS, 200 x 9.4mm, 2µm,1000Å	WAX	Semi-Prep	SS	10-32 coned	200	9.4	2	1000	2450
051WX0301	PolyWAX LP® Column, Glass Lined SS, 50 x 1.0mm, 3µm,100Å	WAX	Analytical	GLT	10-32 coned	50	1.0	3	100	570
051WX0303	PolyWAX LP® Column, Glass Lined SS, 50 x 1.0mm, 3µm,300Å	WAX	Analytical	GLT	10-32 coned	50	1.0	3	300	570
051WX0310	PolyWAX LP® Column, Glass Lined SS, 50 x 1.0mm, 3µm,1000Å	WAX	Analytical	GLT	10-32 coned	50	1.0	3	1000	570
051WX0315	PolyWAX LP® Column, Glass Lined SS, 50 x 1.0mm, 3µm,1500Å	WAX	Analytical	GLT	10-32 coned	50	1.0	3	1500	570
101WX0301	PolyWAX LP® Column, Glass Lined SS, 100 x 1.0mm, 3µm,100Å	WAX	Analytical	GLT	10-32 coned	100	1.0	3	100	660
101WX0303	PolyWAX LP® Column, Glass Lined SS, 100 x 1.0mm, 3µm,300Å	WAX	Analytical	GLT	10-32 coned	100	1.0	3	300	660
101WX0310	PolyWAX LP® Column, Glass Lined SS, 100 x 1.0mm, 3µm,1000Å	WAX	Analytical	GLT	10-32 coned	100	1.0	3	1000	660
101WX0315	PolyWAX LP® Column, Glass Lined SS, 100 x 1.0mm, 3µm,1500Å	WAX	Analytical	GLT	10-32 coned	100	1.0	3	1500	660
151WX0301	PolyWAX LP® Column, Glass Lined SS, 150 x 1.0mm, 3µm,100Å	WAX	Analytical	GLT	10-32 coned	150	1.0	3	100	735
151WX0303	PolyWAX LP® Column, Glass Lined SS, 150 x 1.0mm, 3µm,300Å	WAX	Analytical	GLT	10-32 coned	150	1.0	3	300	735
151WX0310	PolyWAX LP® Column, Glass Lined SS, 150 x 1.0mm, 3µm,1000Å	WAX	Analytical	GLT	10-32 coned	150	1.0	3	1000	735
151WX0315	PolyWAX LP® Column, Glass Lined SS, 150 x 1.0mm, 3µm,1500Å	WAX	Analytical	GLT	10-32 coned	150	1.0	3	1500	735
3.52WX0301	PolyWAX LP® Column, SS, 35 x 2.1mm, 3µm,100Å	WAX	Analytical	SS	10-32 coned	35	2.1	3	100	550
3.52WX0303	PolyWAX LP® Column, SS, 35 x 2.1mm, 3µm,300Å	WAX	Analytical	SS	10-32 coned	35	2.1	3	300	550
3.52WX0310	PolyWAX LP® Column, SS, 35 x 2.1mm, 3µm,1000Å	WAX	Analytical	SS	10-32 coned	35	2.1	3	1000	550
3.52WX0315	PolyWAX LP® Column, SS, 35 x 2.1mm, 3µm,1500Å	WAX	Analytical	SS	10-32 coned	35	2.1	3	1500	550
052WX0301	PolyWAX LP® Column, SS, 50 x 2.1mm, 3µm,100Å	WAX	Analytical	SS	10-32 coned	50	2.1	3	100	560
052WX0303	PolyWAX LP® Column, SS, 50 x 2.1mm, 3µm,300Å	WAX	Analytical	SS	10-32 coned	50	2.1	3	300	560

PolyWAX LP® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
052WX0310	PolyWAX LP® Column, SS, 50 x 2.1mm, 3µm,1000Å	WAX	Analytical	SS	10-32 coned	50	2.1	3	1000	560
052WX0315	PolyWAX LP® Column, SS, 50 x 2.1mm, 3µm,1500Å	WAX	Analytical	SS	10-32 coned	50	2.1	3	1500	560
102WX0301	PolyWAX LP® Column, SS, 100 x 2.1mm, 3µm,100Å	WAX	Analytical	SS	10-32 coned	100	2.1	3	100	710
102WX0303	PolyWAX LP® Column, SS, 100 x 2.1mm, 3µm,300Å	WAX	Analytical	SS	10-32 coned	100	2.1	3	300	710
102WX0310	PolyWAX LP® Column, SS, 100 x 2.1mm, 3µm,1000Å	WAX	Analytical	SS	10-32 coned	100	2.1	3	1000	710
102WX0315	PolyWAX LP® Column, SS, 100 x 2.1mm, 3µm,1500Å	WAX	Analytical	SS	10-32 coned	100	2.1	3	1500	710
3.54WX0301	PolyWAX LP® Column, SS, 35 x 4.6mm, 3µm,100Å	WAX	Analytical	SS	10-32 coned	35	4.6	3	100	545
3.54WX0303	PolyWAX LP® Column, SS, 35 x 4.6mm, 3µm,300Å	WAX	Analytical	SS	10-32 coned	35	4.6	3	300	545
3.54WX0310	PolyWAX LP® Column, SS, 35 x 4.6mm, 3µm,1000Å	WAX	Analytical	SS	10-32 coned	35	4.6	3	1000	545
3.54WX0315	PolyWAX LP® Column, SS, 35 x 4.6mm, 3µm,1500Å	WAX	Analytical	SS	10-32 coned	35	4.6	3	1500	545
054WX0301	PolyWAX LP® Column, SS, 50 x 4.6mm, 3µm,100Å	WAX	Analytical	SS	10-32 coned	50	4.6	3	100	560
054WX0303	PolyWAX LP® Column, SS, 50 x 4.6mm, 3µm,300Å	WAX	Analytical	SS	10-32 coned	50	4.6	3	300	560
054WX0310	PolyWAX LP® Column, SS, 50 x 4.6mm, 3µm,1000Å	WAX	Analytical	SS	10-32 coned	50	4.6	3	1000	560
054WX0315	PolyWAX LP® Column, SS, 50 x 4.6mm, 3µm,1500Å	WAX	Analytical	SS	10-32 coned	50	4.6	3	1500	560
104WX0301	PolyWAX LP® Column, SS, 100 x 4.6mm, 3µm,100Å	WAX	Analytical	SS	10-32 coned	100	4.6	3	100	710
104WX0303	PolyWAX LP® Column, SS, 100 x 4.6mm, 3µm,300Å	WAX	Analytical	SS	10-32 coned	100	4.6	3	300	710
104WX0310	PolyWAX LP® Column, SS, 100 x 4.6mm, 3µm,1000Å	WAX	Analytical	SS	10-32 coned	100	4.6	3	1000	710
104WX0315	PolyWAX LP® Column, SS, 100 x 4.6mm, 3µm,1500Å	WAX	Analytical	SS	10-32 coned	100	4.6	3	1500	710
154WX0301	PolyWAX LP® Column, SS, 150 x 4.6mm, 3µm,100Å	WAX	Analytical	SS	10-32 coned	150	4.6	3	100	740
154WX0303	PolyWAX LP® Column, SS, 150 x 4.6mm, 3µm,300Å	WAX	Analytical	SS	10-32 coned	150	4.6	3	300	740
154WX0310	PolyWAX LP® Column, SS, 150 x 4.6mm, 3µm,1000Å	WAX	Analytical	SS	10-32 coned	150	4.6	3	1000	740
154WX0315	PolyWAX LP® Column, SS, 150 x 4.6mm, 3µm,1500Å	WAX	Analytical	SS	10-32 coned	150	4.6	3	1500	740
204WX0301	PolyWAX LP® Column, SS, 200 x 4.6mm, 3µm,100Å	WAX	Analytical	SS	10-32 coned	200	4.6	3	100	760
204WX0303	PolyWAX LP® Column, SS, 200 x 4.6mm, 3µm,300Å	WAX	Analytical	SS	10-32 coned	200	4.6	3	300	760
204WX0310	PolyWAX LP® Column, SS, 200 x 4.6mm, 3µm,1000Å	WAX	Analytical	SS	10-32 coned	200	4.6	3	1000	760
204WX0315	PolyWAX LP® Column, SS, 200 x 4.6mm, 3µm,1500Å	WAX	Analytical	SS	10-32 coned	200	4.6	3	1500	760
254WX0301	PolyWAX LP® Column, SS, 250 x 4.6mm, 3µm,100Å	WAX	Analytical	SS	10-32 coned	250	4.6	3	100	780
254WX0303	PolyWAX LP® Column, SS, 250 x 4.6mm, 3µm,300Å	WAX	Analytical	SS	10-32 coned	250	4.6	3	300	780
254WX0310	PolyWAX LP® Column, SS, 250 x 4.6mm, 3µm,1000Å	WAX	Analytical	SS	10-32 coned	250	4.6	3	1000	780
254WX0315	PolyWAX LP® Column, SS, 250 x 4.6mm, 3µm,1500Å	WAX	Analytical	SS	10-32 coned	250	4.6	3	1500	780
109WX0301	PolyWAX LP® Column, SS, 100 x 9.4mm, 3µm,100Å	WAX	Analytical	SS	10-32 coned	100	9.4	3	100	1400
109WX0303	PolyWAX LP® Column, SS, 100 x 9.4mm, 3µm,300Å	WAX	Analytical	SS	10-32 coned	100	9.4	3	300	1400
109WX0310	PolyWAX LP® Column, SS, 100 x 9.4mm, 3µm,1000Å	WAX	Analytical	SS	10-32 coned	100	9.4	3	1000	1400
109WX0315	PolyWAX LP® Column, SS, 100 x 9.4mm, 3µm,1500Å	WAX	Analytical	SS	10-32 coned	100	9.4	3	1500	1400
209WX0301	PolyWAX LP® Column, SS, 200 x 9.4mm, 3µm,100Å	WAX	Analytical	SS	10-32 coned	200	9.4	3	100	1620
209WX0303	PolyWAX LP® Column, SS, 200 x 9.4mm, 3µm,300Å	WAX	Analytical	SS	10-32 coned	200	9.4	3	300	1620
209WX0310	PolyWAX LP® Column, SS, 200 x 9.4mm, 3µm,1000Å	WAX	Analytical	SS	10-32 coned	200	9.4	3	1000	1620
209WX0315	PolyWAX LP® Column, SS, 200 x 9.4mm, 3µm,1500Å	WAX	Analytical	SS	10-32 coned	200	9.4	3	1500	1620
259WX0301	PolyWAX LP® Column, SS, 250 x 9.4mm, 3µm,100Å	WAX	Analytical	SS	10-32 coned	250	9.4	3	100	1750
259WX0303	PolyWAX LP® Column, SS, 250 x 9.4mm, 3µm,300Å	WAX	Analytical	SS	10-32 coned	250	9.4	3	300	1750
259WX0310	PolyWAX LP® Column, SS, 250 x 9.4mm, 3µm,1000Å	WAX	Analytical	SS	10-32 coned	250	9.4	3	1000	1750
259WX0315	PolyWAX LP® Column, SS, 250 x 9.4mm, 3µm,1500Å	WAX	Analytical	SS	10-32 coned	250	9.4	3	1500	1750
051WX0501	PolyWAX LP® Column, Glass Lined SS, 50 x 1.0mm, 5µm,100Å	WAX	Analytical	GLT	10-32 coned	50	1.0	5	100	570
051WX0503	PolyWAX LP® Column, Glass Lined SS, 50 x 1.0mm, 5µm,300Å	WAX	Analytical	GLT	10-32 coned	50	1.0	5	300	570
051WX0510	PolyWAX LP® Column, Glass Lined SS, 50 x 1.0mm, 5µm,1000Å	WAX	Analytical	GLT	10-32 coned	50	1.0	5	1000	570

PolyWAX LP® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
101WX0501	PolyWAX LP® Column, Glass Lined SS, 100 x 1.0mm, 5µm,100Å	WAX	Analytical	GLT	10-32 coned	100	1.0	5	100	660
101WX0503	PolyWAX LP® Column, Glass Lined SS, 100 x 1.0mm, 5µm,300Å	WAX	Analytical	GLT	10-32 coned	100	1.0	5	300	660
101WX0510	PolyWAX LP® Column, Glass Lined SS, 100 x 1.0mm, 5µm,1000Å	WAX	Analytical	GLT	10-32 coned	100	1.0	5	1000	660
151WX0501	PolyWAX LP® Column, Glass Lined SS, 150 x 1.0mm, 5µm,100Å	WAX	Analytical	GLT	10-32 coned	150	1.0	5	100	740
151WX0503	PolyWAX LP® Column, Glass Lined SS, 150 x 1.0mm, 5µm,300Å	WAX	Analytical	GLT	10-32 coned	150	1.0	5	300	740
151WX0510	PolyWAX LP® Column, Glass Lined SS, 150 x 1.0mm, 5µm,1000Å	WAX	Analytical	GLT	10-32 coned	150	1.0	5	1000	740
3.52WX0501	PolyWAX LP® Column, SS, 35 x 2.1mm, 5µm,100Å	WAX	Analytical	SS	10-32 coned	35	2.1	5	100	500
3.52WX0503	PolyWAX LP® Column, SS, 35 x 2.1mm, 5µm,300Å	WAX	Analytical	SS	10-32 coned	35	2.1	5	300	500
3.52WX0510	PolyWAX LP® Column, SS, 35 x 2.1mm, 5µm,1000Å	WAX	Analytical	SS	10-32 coned	35	2.1	5	1000	500
052WX0501	PolyWAX LP® Column, SS, 50 x 2.1mm, 5µm,100Å	WAX	Analytical	SS	10-32 coned	50	2.1	5	100	520
052WX0503	PolyWAX LP® Column, SS, 50 x 2.1mm, 5µm,300Å	WAX	Analytical	SS	10-32 coned	50	2.1	5	300	520
052WX0510	PolyWAX LP® Column, SS, 50 x 2.1mm, 5µm,1000Å	WAX	Analytical	SS	10-32 coned	50	2.1	5	1000	520
102WX0501	PolyWAX LP® Column, SS, 100 x 2.1mm, 5µm,100Å	WAX	Analytical	SS	10-32 coned	100	2.1	5	100	600
102WX0503	PolyWAX LP® Column, SS, 100 x 2.1mm, 5µm,300Å	WAX	Analytical	SS	10-32 coned	100	2.1	5	300	600
102WX0510	PolyWAX LP® Column, SS, 100 x 2.1mm, 5µm,1000Å	WAX	Analytical	SS	10-32 coned	100	2.1	5	1000	600
152WX0501	PolyWAX LP® Column, SS, 150 x 2.1mm, 5µm,100Å	WAX	Analytical	SS	10-32 coned	150	2.1	5	100	650
152WX0503	PolyWAX LP® Column, SS, 150 x 2.1mm, 5µm,300Å	WAX	Analytical	SS	10-32 coned	150	2.1	5	300	650
152WX0510	PolyWAX LP® Column, SS, 150 x 2.1mm, 5µm,1000Å	WAX	Analytical	SS	10-32 coned	150	2.1	5	1000	650
202WX0501	PolyWAX LP® Column, SS, 200 x 2.1mm, 5µm,100Å	WAX	Analytical	SS	10-32 coned	200	2.1	5	100	710
202WX0503	PolyWAX LP® Column, SS, 200 x 2.1mm, 5µm,300Å	WAX	Analytical	SS	10-32 coned	200	2.1	5	300	710
202WX0510	PolyWAX LP® Column, SS, 200 x 2.1mm, 5µm,1000Å	WAX	Analytical	SS	10-32 coned	200	2.1	5	1000	710
3.54WX0501	PolyWAX LP® Column, SS, 35 x 4.6mm, 5µm,100Å	WAX	Analytical	SS	10-32 coned	35	4.6	5	100	500
3.54WX0503	PolyWAX LP® Column, SS, 35 x 4.6mm, 5µm,300Å	WAX	Analytical	SS	10-32 coned	35	4.6	5	300	500
3.54WX0510	PolyWAX LP® Column, SS, 35 x 4.6mm, 5µm,1000Å	WAX	Analytical	SS	10-32 coned	35	4.6	5	1000	500
054WX0501	PolyWAX LP® Column, SS, 50 x 4.6mm, 5µm,100Å	WAX	Analytical	SS	10-32 coned	50	4.6	5	100	520
054WX0503	PolyWAX LP® Column, SS, 50 x 4.6mm, 5µm,300Å	WAX	Analytical	SS	10-32 coned	50	4.6	5	300	520
054WX0510	PolyWAX LP® Column, SS, 50 x 4.6mm, 5µm,1000Å	WAX	Analytical	SS	10-32 coned	50	4.6	5	1000	520
104WX0501	PolyWAX LP® Column, SS, 100 x 4.6mm, 5µm,100Å	WAX	Analytical	SS	10-32 coned	100	4.6	5	100	600
104WX0503	PolyWAX LP® Column, SS, 100 x 4.6mm, 5µm,300Å	WAX	Analytical	SS	10-32 coned	100	4.6	5	300	600
104WX0510	PolyWAX LP® Column, SS, 100 x 4.6mm, 5µm,1000Å	WAX	Analytical	SS	10-32 coned	100	4.6	5	1000	600
154WX0501	PolyWAX LP® Column, SS, 150 x 4.6mm, 5µm,100Å	WAX	Analytical	SS	10-32 coned	150	4.6	5	100	650
154WX0503	PolyWAX LP® Column, SS, 150 x 4.6mm, 5µm,300Å	WAX	Analytical	SS	10-32 coned	150	4.6	5	300	650
154WX0510	PolyWAX LP® Column, SS, 150 x 4.6mm, 5µm,1000Å	WAX	Analytical	SS	10-32 coned	150	4.6	5	1000	650
204WX0501	PolyWAX LP® Column, SS, 200 x 4.6mm, 5µm,100Å	WAX	Analytical	SS	10-32 coned	200	4.6	5	100	710
204WX0503	PolyWAX LP® Column, SS, 200 x 4.6mm, 5µm,300Å	WAX	Analytical	SS	10-32 coned	200	4.6	5	300	710
204WX0510	PolyWAX LP® Column, SS, 200 x 4.6mm, 5µm,1000Å	WAX	Analytical	SS	10-32 coned	200	4.6	5	1000	710
254WX0501	PolyWAX LP® Column, SS, 250 x 4.6mm, 5µm,100Å	WAX	Analytical	SS	10-32 coned	250	4.6	5	100	740
254WX0503	PolyWAX LP® Column, SS, 250 x 4.6mm, 5µm,300Å	WAX	Analytical	SS	10-32 coned	250	4.6	5	300	740
254WX0510	PolyWAX LP® Column, SS, 250 x 4.6mm, 5µm,1000Å	WAX	Analytical	SS	10-32 coned	250	4.6	5	1000	740
109WX0501	PolyWAX LP® Column, SS, 100 x 9.4mm, 5µm,100Å	WAX	Semi-Prep	SS	10-32 coned	100	9.4	5	100	1110
109WX0503	PolyWAX LP® Column, SS, 100 x 9.4mm, 5µm,300Å	WAX	Semi-Prep	SS	10-32 coned	100	9.4	5	300	1110
109WX0510	PolyWAX LP® Column, SS, 100 x 9.4mm, 5µm,1000Å	WAX	Semi-Prep	SS	10-32 coned	100	9.4	5	1000	1110
209WX0501	PolyWAX LP® Column, SS, 200 x 9.4mm, 5µm,100Å	WAX	Semi-Prep	SS	10-32 coned	200	9.4	5	100	1400
209WX0503	PolyWAX LP® Column, SS, 200 x 9.4mm, 5µm,300Å	WAX	Semi-Prep	SS	10-32 coned	200	9.4	5	300	1400
209WX0510	PolyWAX LP® Column, SS, 200 x 9.4mm, 5µm,1000Å	WAX	Semi-Prep	SS	10-32 coned	200	9.4	5	1000	1400

PolyWAX LP® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
259WX0501	PolyWAX LP® Column, SS, 250 x 9.4mm, 5µm,100Å	WAX	Semi-Prep	SS	10-32 coned	250	9.4	5	100	1500
259WX0503	PolyWAX LP® Column, SS, 250 x 9.4mm, 5µm,300Å	WAX	Semi-Prep	SS	10-32 coned	250	9.4	5	300	1500
259WX0510	PolyWAX LP® Column, SS, 250 x 9.4mm, 5µm,1000Å	WAX	Semi-Prep	SS	10-32 coned	250	9.4	5	1000	1500
2521WX0501	PolyWAX LP® Column, SS, 250 x 21mm, 5µm,100Å	WAX	Semi-Prep	SS	10-32 coned	250	21.0	5	100	4500
2521WX0503	PolyWAX LP® Column, SS, 250 x 21mm, 5µm,300Å	WAX	Semi-Prep	SS	10-32 coned	250	21.0	5	300	4500
2521WX0510	PolyWAX LP® Column, SS, 250 x 21mm, 5µm,1000Å	WAX	Semi-Prep	SS	10-32 coned	250	21.0	5	1000	4500
204WX1201	PolyWAX LP® Column, SS, 200 x 4.6mm, 12µm,100Å	WAX	Analytical	SS	10-32 coned	200	4.6	12	100	420
204WX1203	PolyWAX LP® Column, SS, 200 x 4.6mm, 12µm,300Å	WAX	Analytical	SS	10-32 coned	200	4.6	12	300	420
204WX1215	PolyWAX LP® Column, SS, 200 x 4.6mm, 12µm,1500Å	WAX	Analytical	SS	10-32 coned	200	4.6	12	1500	420
254WX1201	PolyWAX LP® Column, SS, 250 x 4.6mm, 12µm,100Å	WAX	Analytical	SS	10-32 coned	250	4.6	12	100	500
254WX1203	PolyWAX LP® Column, SS, 250 x 4.6mm, 12µm,300Å	WAX	Analytical	SS	10-32 coned	250	4.6	12	300	500
254WX1215	PolyWAX LP® Column, SS, 250 x 4.6mm, 12µm,1500Å	WAX	Analytical	SS	10-32 coned	250	4.6	12	1500	500
209WX1201	PolyWAX LP® Column, SS, 200 x 9.4mm, 12µm,100Å	WAX	Semi-Prep	SS	10-32 coned	200	9.4	12	100	1000
209WX1203	PolyWAX LP® Column, SS, 200 x 9.4mm, 12µm,300Å	WAX	Semi-Prep	SS	10-32 coned	200	9.4	12	300	1000
209WX1215	PolyWAX LP® Column, SS, 200 x 9.4mm, 12µm,1500Å	WAX	Semi-Prep	SS	10-32 coned	200	9.4	12	1500	1000
259WX1201	PolyWAX LP® Column, SS, 250 x 9.4mm, 12µm,100Å	WAX	Semi-Prep	SS	10-32 coned	250	9.4	12	100	1100
259WX1203	PolyWAX LP® Column, SS, 250 x 9.4mm, 12µm,300Å	WAX	Semi-Prep	SS	10-32 coned	250	9.4	12	300	1100
259WX1215	PolyWAX LP® Column, SS, 250 x 9.4mm, 12µm,1500Å	WAX	Semi-Prep	SS	10-32 coned	250	9.4	12	1500	1100
2521WX1201	PolyWAX LP® Column, SS, 250 x 21mm, 12µm,100Å	WAX	Semi-Prep	SS	10-32 coned	250	21.0	12	100	3200
2521WX1203	PolyWAX LP® Column, SS, 250 x 21mm, 12µm,300Å	WAX	Semi-Prep	SS	10-32 coned	250	21.0	12	300	3200
2521WX1215	PolyWAX LP® Column, SS, 250 x 21mm, 12µm,1500Å	WAX	Semi-Prep	SS	10-32 coned	250	21.0	12	1500	3200
2550WX1201	PolyWAX LP® Column, SS, 250 x 50.8mm, 12µm,100Å	WAX	Prep	SS	10-32 coned	250	50.8	12	100	12250
2550WX1203	PolyWAX LP® Column, SS, 250 x 50.8mm, 12µm,300Å	WAX	Prep	SS	10-32 coned	250	50.8	12	300	12250
2550WX1215	PolyWAX LP® Column, SS, 250 x 50.8mm, 12µm,1500Å	WAX	Prep	SS	10-32 coned	250	50.8	12	1500	12250

HPLC Capillaries

PolyWAX LP® Capillaries (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
050.075WX0210	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 2µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.075	2	1000	520
100.075WX0210	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 2µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.075	2	1000	580
150.075WX0210	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 2µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.075	2	1000	640
200.075WX0210	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 2µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.075	2	1000	700
050.15WX0210	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 2µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.150	2	1000	520
100.15WX0210	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 2µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.150	2	1000	580
150.15WX0210	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 2µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.150	2	1000	640
200.15WX0210	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 2µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.150	2	1000	700



PolyWAX LP® Capillaries (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
050.30WX0210	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 2µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.300	2	1000	520
100.30WX0210	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 2µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.300	2	1000	580
150.30WX0210	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 2µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.300	2	1000	640
200.30WX0210	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 2µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.300	2	1000	700
050.075WX0301	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.075	3	100	520
050.075WX0303	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.075	3	300	520
050.075WX0310	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.075	3	1000	520
050.075WX0315	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,1500Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.075	3	1500	520
100.075WX0301	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.075	3	100	580
100.075WX0303	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.075	3	300	580
100.075WX0310	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.075	3	1000	580
100.075WX0315	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,1500Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.075	3	1500	580
150.075WX0301	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.075	3	100	640
150.075WX0303	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.075	3	300	640
150.075WX0310	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.075	3	1000	640
150.075WX0315	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,1500Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.075	3	1500	640
200.075WX0301	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.075	3	100	700
200.075WX0303	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.075	3	300	700
200.075WX0310	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.075	3	1000	700
200.075WX0315	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,1500Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.075	3	1500	700
050.15WX0301	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.150	3	100	520
050.15WX0303	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.150	3	300	520
050.15WX0310	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.150	3	1000	520
050.15WX0315	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,1500Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.150	3	1500	520
100.15WX0301	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.150	3	100	580
100.15WX0303	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.150	3	300	580
100.15WX0310	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.150	3	1000	580
100.15WX0315	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,1500Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.150	3	1500	580
150.15WX0301	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.150	3	100	640
150.15WX0303	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.150	3	300	640
150.15WX0310	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.150	3	1000	640
150.15WX0315	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,1500Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.150	3	1500	640
200.15WX0301	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.150	3	100	700
200.15WX0303	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.150	3	300	700
200.15WX0310	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.150	3	1000	700
200.15WX0315	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,1500Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.150	3	1500	700
050.30WX0301	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.300	3	100	520
050.30WX0303	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.300	3	300	520
050.30WX0310	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.300	3	1000	520
050.30WX0315	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,1500Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.300	3	1500	520
100.30WX0301	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.300	3	100	580
100.30WX0303	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.300	3	300	580
100.30WX0310	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.300	3	1000	580
100.30WX0315	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,1500Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.300	3	1500	580
150.30WX0301	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.300	3	100	640
150.30WX0303	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.300	3	300	640

PolyWAX LP® Capillaries (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
150.30WX0310	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.300	3	1000	640
150.30WX0315	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,1500Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.300	3	1500	640
200.30WX0301	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.300	3	100	700
200.30WX0303	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.300	3	300	700
200.30WX0310	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.300	3	1000	700
200.30WX0315	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,1500Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.300	3	1500	700
050.075WX0501	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 5µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.075	5	100	520
050.075WX0503	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 5µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.075	5	300	520
050.075WX0510	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 5µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.075	5	1000	520
100.075WX0501	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 5µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.075	5	100	580
100.075WX0503	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 5µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.075	5	300	580
100.075WX0510	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 5µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.075	5	1000	580
150.075WX0501	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 5µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.075	5	100	640
150.075WX0503	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 5µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.075	5	300	640
150.075WX0510	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 5µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.075	5	1000	640
200.075WX0501	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 5µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.075	5	100	700
200.075WX0503	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 5µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.075	5	300	700
200.075WX0510	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 5µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.075	5	1000	700
050.15WX0501	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 5µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.150	5	100	520
050.15WX0503	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 5µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.150	5	300	520
050.15WX0510	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 5µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.150	5	1000	520
100.15WX0501	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 5µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.150	5	100	580
100.15WX0503	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 5µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.150	5	300	580
100.15WX0510	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 5µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.150	5	1000	580
150.15WX0501	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 5µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.150	5	100	640
150.15WX0503	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 5µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.150	5	300	640
150.15WX0510	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 5µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.150	5	1000	640
200.15WX0501	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 5µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.150	5	100	700
200.15WX0503	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 5µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.150	5	300	700
200.15WX0510	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 5µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.150	5	1000	700
050.30WX0501	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 5µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.300	5	100	520
050.30WX0503	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 5µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.300	5	300	520
050.30WX0510	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 5µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	50	0.300	5	1000	520
100.30WX0501	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 5µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.300	5	100	580
100.30WX0503	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 5µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.300	5	300	580
100.30WX0510	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 5µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	100	0.300	5	1000	580
150.30WX0501	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 5µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.300	5	100	640
150.30WX0503	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 5µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.300	5	300	640
150.30WX0510	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 5µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	150	0.300	5	1000	640
200.30WX0501	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 5µm,100Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.300	5	100	700
200.30WX0503	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 5µm,300Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.300	5	300	700
200.30WX0510	PolyWAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 5µm,1000Å	WAX	Analytical	PEEK SIL	10-32 coned	200	0.300	5	1000	700



Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	Wax Inlet L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
GC1-WX0210	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 1mm, 2µm, 1000Å	WAX	Analytical	SS/PEEK	10-32 coned	10	1.0	2	1000	120
GC2-WX0210	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 2µm, 1000Å	WAX	Analytical	SS/PEEK	10-32 coned	10	2.0	2	1000	120
GC3-WX0210	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 2µm, 1000Å	WAX	Analytical	SS/PEEK	10-32 coned	10	3.0	2	1000	120
GC4-WX0210	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 2µm, 1000Å	WAX	Analytical	SS/PEEK	10-32 coned	10	4.0	2	1000	120
GC1-WX0301	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 100Å	WAX	Analytical	SS/PEEK	10-32 coned	10	1.0	3	100	100
GC1-WX0303	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 300Å	WAX	Analytical	SS/PEEK	10-32 coned	10	1.0	3	300	100
GC1-WX0310	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 1000Å	WAX	Analytical	SS/PEEK	10-32 coned	10	1.0	3	1000	100
GC1-WX0315	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 1500Å	WAX	Analytical	SS/PEEK	10-32 coned	10	1.0	3	1500	100
GC2-WX0301	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 3µm, 100Å	WAX	Analytical	SS/PEEK	10-32 coned	10	2.0	3	100	100
GC2-WX0303	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 3µm, 300Å	WAX	Analytical	SS/PEEK	10-32 coned	10	2.0	3	300	100
GC2-WX0310	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 3µm, 1000Å	WAX	Analytical	SS/PEEK	10-32 coned	10	2.0	3	1000	100
GC2-WX0315	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 3µm, 1500Å	WAX	Analytical	SS/PEEK	10-32 coned	10	2.0	3	1500	100
GC3-WX0301	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 3µm, 100Å	WAX	Analytical	SS/PEEK	10-32 coned	10	3.0	3	100	100
GC3-WX0303	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 3µm, 300Å	WAX	Analytical	SS/PEEK	10-32 coned	10	3.0	3	300	100
GC3-WX0310	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 3µm, 1000Å	WAX	Analytical	SS/PEEK	10-32 coned	10	3.0	3	1000	100
GC3-WX0315	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 3µm, 1500Å	WAX	Analytical	SS/PEEK	10-32 coned	10	3.0	3	1500	100
GC4-WX0301	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 3µm, 100Å	WAX	Analytical	SS/PEEK	10-32 coned	10	4.0	3	100	100
GC4-WX0303	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 3µm, 300Å	WAX	Analytical	SS/PEEK	10-32 coned	10	4.0	3	300	100
GC4-WX0310	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 3µm, 1000Å	WAX	Analytical	SS/PEEK	10-32 coned	10	4.0	3	1000	100
GC4-WX0315	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 3µm, 1500Å	WAX	Analytical	SS/PEEK	10-32 coned	10	4.0	3	1500	100
GC1-WX0501	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 1mm, 5µm, 100Å	WAX	Analytical	SS/PEEK	10-32 coned	10	1.0	5	100	90
GC1-WX0503	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 1mm, 5µm, 300Å	WAX	Analytical	SS/PEEK	10-32 coned	10	1.0	5	300	90
GC1-WX0510	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 1mm, 5µm, 1000Å	WAX	Analytical	SS/PEEK	10-32 coned	10	1.0	5	1000	90
GC2-WX0501	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 5µm, 100Å	WAX	Analytical	SS/PEEK	10-32 coned	10	2.0	5	100	90
GC2-WX0503	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 5µm, 300Å	WAX	Analytical	SS/PEEK	10-32 coned	10	2.0	5	300	90
GC2-WX0510	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 5µm, 1000Å	WAX	Analytical	SS/PEEK	10-32 coned	10	2.0	5	1000	90
GC3-WX0501	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 5µm, 100Å	WAX	Analytical	SS/PEEK	10-32 coned	10	3.0	5	100	90
GC3-WX0503	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 5µm, 300Å	WAX	Analytical	SS/PEEK	10-32 coned	10	3.0	5	300	90
GC3-WX0510	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 5µm, 1000Å	WAX	Analytical	SS/PEEK	10-32 coned	10	3.0	5	1000	90
GC4-WX0501	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 5µm, 100Å	WAX	Analytical	SS/PEEK	10-32 coned	10	4.0	5	100	90
GC4-WX0503	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 5µm, 300Å	WAX	Analytical	SS/PEEK	10-32 coned	10	4.0	5	300	90
GC4-WX0510	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 5µm, 1000Å	WAX	Analytical	SS/PEEK	10-32 coned	10	4.0	5	1000	90
DGCH10	DireWX Guard Cartridge Holder, 10mm (For all analytical guards, REUSABLE)	N/A	Analytical	SS/PEEK	10-32 coned	10	N/A	N/A	N/A	75
IGCH10	Inline Guard Cartridge Holder, 10mm (For all analytical guards, REUSABLE)	N/A	Analytical	SS/PEEK	10-32 coned	10	N/A	N/A	N/A	75
*** NOTE: DGCH10 for DirectConnection to Column (female inlet/male outlet) IGCH10 for Connection to Line (two female ports)										

GC10-WX0301	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 100Å	WAX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	100	225
GC10-WX0303	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 300Å	WAX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	300	225
GC10-WX0310	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 1000Å	WAX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	1000	225
GC10-WX0315	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 1500Å	WAX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	1500	225
GC10-WX0501	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 10mm, 5µm, 100Å	WAX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	5	100	185
GC10-WX0503	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 10mm, 5µm, 300Å	WAX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	5	300	185
GC10-WX0510	PolyWAX LP® Guard Cartridge, SS/PEEK, 10mm x 10mm, 5µm, 1000Å	WAX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	5	1000	185
CH10	Semi-Prep Guard Cartridge Holder, 10mm (For 9.4mm columns, REUSABLE)	WAX	Semi-Prep	SS/PEEK	10-32 coned	10	N/A	N/A	N/A	175

*** NOTE: CH10 Holder Required for Semi-Prep Guard Cartridges (10 x 10mm) Connection to Line (two female ports)

Solid Phase Extraction Cartridges

PolyWAX LP® SPE Cartridges

PN	Product Description	Phase	Type	Hardware	ConneXion	V (ml)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
SPEWX1201	PolyWAX LP® SPE Cartridge 1mL, PE, 12µm, 100Å (Pack of 10 cartridges)	WAX	SPE	PE	Lure	1	N/A	12	100	130
SPEWX1203	PolyWAX LP® SPE Cartridge 1mL, PE, 12µm, 300Å (Pack of 10 cartridges)	WAX	SPE	PE	Lure	1	N/A	12	300	130
SPEWX1215	PolyWAX LP® SPE Cartridge 1mL, PE, 12µm, 1500Å (Pack of 10 cartridges)	WAX	SPE	PE	Lure	1	N/A	12	1500	130
SPEWX2003	PolyWAX LP® SPE Cartridge 1mL, PE, 20µm, 300Å (Pack of 10 cartridges)	WAX	SPE	PE	Lure	1	N/A	20	300	130
SPEWX2015	PolyWAX LP® SPE Cartridge 1mL, PE, 20µm, 1500Å (Pack of 10 cartridges)	WAX	SPE	PE	Lure	1	N/A	20	1500	130

Custom combinations available upon request info@polylc.com

Spin Columns

PolyWAX LP® BioPureSPN Columns

PN	Product Description	Phase	Type	Bed (mg)	Pore (Å)	WAX Capacity	HILIC Capacity	2026 Price (USD)
HMM HIL-WAX	PolyWAX LP® BioPureSPN MACRO for WAX or HILIC 50-200µl (Box of 100 each)	WAX	Spin	100	300	0.2 -1mg	35 - 350µg	940
HMM HIL-WAX.20	PolyWAX LP® BioPureSPN MACRO for WAX or HILIC 50-200µl (Pack of 20 each)	WAX	Spin	100	300	0.2 -1mg	35 - 350µg	265
HMM HIL-WAX-L	PolyWAX LP® BioPureSPN MACRO PLATE 96 Well for WAX or HILIC 50-200µl	WAX	Spin	100	300	0.2 -1mg	35 - 350µg	800
HEM HIL-WAX	PolyWAX LP® BioPureSPN MIDI for WAX/HILIC 50-200µl (Box of 100 each)	WAX	Spin	50	300	100-500µg	17 - 170µg	820
HEM HIL-WAX.20	PolyWAX LP® BioPureSPN MIDI for WAX/HILIC 50-200µl (Pack of 20 each)	WAX	Spin	50	300	100-500µg	17 - 170µg	240
HNS HIL-WAX-M	PolyWAX LP® BioPureSPN MIDI PLATE 96 Well for WAX/HILIC 50-200µl	WAX	Spin	50	300	100-500µg	17 - 170µg	575
HUM HIL-WAX	PolyWAX LP® BioPureSPN MINI for WAX/HILIC 20-100µl (Box of 100 each)	WAX	Spin	20	300	40-200µg	7 - 70µg	735
HUM HIL-WAX.20	PolyWAX LP® BioPureSPN MINI for WAX/HILIC 20-100µl (Pack of 20 each)	WAX	Spin	20	300	40-200µg	7 - 70µg	225
HNS HIL-WAX	PolyWAX LP® BioPureSPN MINI PLATE 96 Well for WAX/HILIC 20-100µl	WAX	Spin	20	300	40-200µg	7 - 70µg	735

Custom combinations available upon request info@polylc.com

Bulk Material



PolyWAX LP® Bulk Material (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.) Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	QTY (g)	Particle (µm)	Pore (Å)	2026 Price (USD)
BMW0210	PolyWAX LP® Bulk Material, 2µm, 1000Å	WAX	BM	N/A	0.5	2	1000	100
BMW0301	PolyWAX LP® Bulk Material, 3µm, 100Å	WAX	BM	N/A	1	3	100	130
BMW0303	PolyWAX LP® Bulk Material, 3µm, 300Å	WAX	BM	N/A	1	3	300	130
BMW0310	PolyWAX LP® Bulk Material, 3µm, 1000Å	WAX	BM	N/A	1	3	1000	130
BMW0315	PolyWAX LP® Bulk Material, 3µm, 1500Å	WAX	BM	N/A	1	3	1500	130
BMW0501	PolyWAX LP® Bulk Material, 5µm, 100Å	WAX	BM	N/A	1	5	100	85
BMW0503	PolyWAX LP® Bulk Material, 5µm, 300Å	WAX	BM	N/A	1	5	300	85
BMW0510	PolyWAX LP® Bulk Material, 5µm, 1000Å	WAX	BM	N/A	1	5	1000	85
BMW1203	PolyWAX LP® Bulk Material, 12µm, 300Å	WAX	BM	N/A	1	12	300	55
BMW1215	PolyWAX LP® Bulk Material, 12µm, 1000Å	WAX	BM	N/A	1	12	1000	55
BMW2001	PolyWAX LP® Bulk Material, 20µm, 100Å	WAX	BM	N/A	1	20	100	55
BMW2003	PolyWAX LP® Bulk Material, 20µm, 300Å	WAX	BM	N/A	1	20	300	55
BMW2015	PolyWAX LP® Bulk Material, 20µm, 1500Å	WAX	BM	N/A	1	20	1500	55
BMW3003	PolyWAX LP® Bulk Material, 30µm, 300Å	WAX	BM	N/A	1	30	300	55
BMW3010	PolyWAX LP® Bulk Material, 30µm, 1000Å	WAX	BM	N/A	1	30	1000	55

Please contact info@polylc.com with any questions or for information on any product combination not listed in the table above

PolyWAX LP®

(Weak-Anion Exchange)

2026 PolyLC Inc. Retail Pricelist

PolyLC Inc., 9151 Rumsey Rd., Suite 180, Columbia, MD 21045 410-992-5400 www.polylc.com



2026 PolyLC Inc. Retail Pricelist



Offering high-performance liquid chromatography media, columns, capillaries, cartridges, and instrumentation. We specialize in advanced alternatives to reversed-phase for proteins, peptides, and other biomolecules, supported by deep application expertise and customized solutions for the most demanding separations.

PolySAX LP®

(Strong-Anion Exchange)



HPLC Columns

PolySAX LP® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
051SX0210	PolySAX LP® Column, Glass Lined SS, 50 x 1.0mm, 2µm,1000Å	SAX	Analytical	GLT	10-32 coned	50	1.0	2	1000	615
101SX0210	PolySAX LP® Column, Glass Lined SS, 100 x 1.0mm, 2µm,1000Å	SAX	Analytical	GLT	10-32 coned	100	1.0	2	1000	750
151SX0210	PolySAX LP® Column, Glass Lined SS, 150 x 1.0mm, 2µm,1000Å	SAX	Analytical	GLT	10-32 coned	150	1.0	2	1000	800
3.52SX0210	PolySAX LP® Column, SS, 35 x 2.1mm, 2µm,1000Å	SAX	Analytical	SS	10-32 coned	35	2.1	2	1000	570
052SX0210	PolySAX LP® Column, SS, 50 x 2.1mm, 2µm,1000Å	SAX	Analytical	SS	10-32 coned	50	2.1	2	1000	580
102SX0210	PolySAX LP® Column, SS, 100 x 2.1mm, 2µm,1000Å	SAX	Analytical	SS	10-32 coned	100	2.1	2	1000	750
3.54SX0210	PolySAX LP® Column, SS, 35 x 4.6mm, 2µm,1000Å	SAX	Analytical	SS	10-32 coned	35	4.6	2	1000	570
054SX0210	PolySAX LP® Column, SS, 50 x 4.6mm, 2µm,1000Å	SAX	Analytical	SS	10-32 coned	50	4.6	2	1000	580
104SX0210	PolySAX LP® Column, SS, 100 x 4.6mm, 2µm,1000Å	SAX	Analytical	SS	10-32 coned	100	4.6	2	1000	750
154SX0210	PolySAX LP® Column, SS, 150 x 4.6mm, 2µm,1000Å	SAX	Analytical	SS	10-32 coned	150	4.6	2	1000	790
204SX0210	PolySAX LP® Column, SS, 200 x 4.6mm, 2µm,1000Å	SAX	Analytical	SS	10-32 coned	200	4.6	2	1000	810
109SX0210	PolySAX LP® Column, SS, 100 x 9.4mm, 2µm,1000Å	SAX	Semi-Prep	SS	10-32 coned	100	9.4	2	1000	1525
209SX0210	PolySAX LP® Column, SS, 200 x 9.4mm, 2µm,1000Å	SAX	Semi-Prep	SS	10-32 coned	200	9.4	2	1000	2450
051SX0301	PolySAX LP® Column, Glass Lined SS, 50 x 1.0mm, 3µm,100Å	SAX	Analytical	GLT	10-32 coned	50	1.0	3	100	570
051SX0303	PolySAX LP® Column, Glass Lined SS, 50 x 1.0mm, 3µm,300Å	SAX	Analytical	GLT	10-32 coned	50	1.0	3	300	570
051SX0310	PolySAX LP® Column, Glass Lined SS, 50 x 1.0mm, 3µm,1000Å	SAX	Analytical	GLT	10-32 coned	50	1.0	3	1000	570
051SX0315	PolySAX LP® Column, Glass Lined SS, 50 x 1.0mm, 3µm,1500Å	SAX	Analytical	GLT	10-32 coned	50	1.0	3	1500	570
101SX0301	PolySAX LP® Column, Glass Lined SS, 100 x 1.0mm, 3µm,100Å	SAX	Analytical	GLT	10-32 coned	100	1.0	3	100	660
101SX0303	PolySAX LP® Column, Glass Lined SS, 100 x 1.0mm, 3µm,300Å	SAX	Analytical	GLT	10-32 coned	100	1.0	3	300	660
101SX0310	PolySAX LP® Column, Glass Lined SS, 100 x 1.0mm, 3µm,1000Å	SAX	Analytical	GLT	10-32 coned	100	1.0	3	1000	660
101SX0315	PolySAX LP® Column, Glass Lined SS, 100 x 1.0mm, 3µm,1500Å	SAX	Analytical	GLT	10-32 coned	100	1.0	3	1500	660
151SX0301	PolySAX LP® Column, Glass Lined SS, 150 x 1.0mm, 3µm,100Å	SAX	Analytical	GLT	10-32 coned	150	1.0	3	100	735
151SX0303	PolySAX LP® Column, Glass Lined SS, 150 x 1.0mm, 3µm,300Å	SAX	Analytical	GLT	10-32 coned	150	1.0	3	300	735
151SX0310	PolySAX LP® Column, Glass Lined SS, 150 x 1.0mm, 3µm,1000Å	SAX	Analytical	GLT	10-32 coned	150	1.0	3	1000	735
151SX0315	PolySAX LP® Column, Glass Lined SS, 150 x 1.0mm, 3µm,1500Å	SAX	Analytical	GLT	10-32 coned	150	1.0	3	1500	735
3.52SX0301	PolySAX LP® Column, SS, 35 x 2.1mm, 3µm,100Å	SAX	Analytical	SS	10-32 coned	35	2.1	3	100	550
3.52SX0303	PolySAX LP® Column, SS, 35 x 2.1mm, 3µm,300Å	SAX	Analytical	SS	10-32 coned	35	2.1	3	300	550
3.52SX0310	PolySAX LP® Column, SS, 35 x 2.1mm, 3µm,1000Å	SAX	Analytical	SS	10-32 coned	35	2.1	3	1000	550
3.52SX0315	PolySAX LP® Column, SS, 35 x 2.1mm, 3µm,1500Å	SAX	Analytical	SS	10-32 coned	35	2.1	3	1500	550
052SX0301	PolySAX LP® Column, SS, 50 x 2.1mm, 3µm,100Å	SAX	Analytical	SS	10-32 coned	50	2.1	3	100	560
052SX0303	PolySAX LP® Column, SS, 50 x 2.1mm, 3µm,300Å	SAX	Analytical	SS	10-32 coned	50	2.1	3	300	560

PolySAX LP® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
052SX0310	PolySAX LP® Column, SS, 50 x 2.1mm, 3µm,1000Å	SAX	Analytical	SS	10-32 coned	50	2.1	3	1000	560
052SX0315	PolySAX LP® Column, SS, 50 x 2.1mm, 3µm,1500Å	SAX	Analytical	SS	10-32 coned	50	2.1	3	1500	560
102SX0301	PolySAX LP® Column, SS, 100 x 2.1mm, 3µm,100Å	SAX	Analytical	SS	10-32 coned	100	2.1	3	100	710
102SX0303	PolySAX LP® Column, SS, 100 x 2.1mm, 3µm,300Å	SAX	Analytical	SS	10-32 coned	100	2.1	3	300	710
102SX0310	PolySAX LP® Column, SS, 100 x 2.1mm, 3µm,1000Å	SAX	Analytical	SS	10-32 coned	100	2.1	3	1000	710
102SX0315	PolySAX LP® Column, SS, 100 x 2.1mm, 3µm,1500Å	SAX	Analytical	SS	10-32 coned	100	2.1	3	1500	710
3.54SX0301	PolySAX LP® Column, SS, 35 x 4.6mm, 3µm,100Å	SAX	Analytical	SS	10-32 coned	35	4.6	3	100	545
3.54SX0303	PolySAX LP® Column, SS, 35 x 4.6mm, 3µm,300Å	SAX	Analytical	SS	10-32 coned	35	4.6	3	300	545
3.54SX0310	PolySAX LP® Column, SS, 35 x 4.6mm, 3µm,1000Å	SAX	Analytical	SS	10-32 coned	35	4.6	3	1000	545
3.54SX0315	PolySAX LP® Column, SS, 35 x 4.6mm, 3µm,1500Å	SAX	Analytical	SS	10-32 coned	35	4.6	3	1500	545
054SX0301	PolySAX LP® Column, SS, 50 x 4.6mm, 3µm,100Å	SAX	Analytical	SS	10-32 coned	50	4.6	3	100	560
054SX0303	PolySAX LP® Column, SS, 50 x 4.6mm, 3µm,300Å	SAX	Analytical	SS	10-32 coned	50	4.6	3	300	560
054SX0310	PolySAX LP® Column, SS, 50 x 4.6mm, 3µm,1000Å	SAX	Analytical	SS	10-32 coned	50	4.6	3	1000	560
054SX0315	PolySAX LP® Column, SS, 50 x 4.6mm, 3µm,1500Å	SAX	Analytical	SS	10-32 coned	50	4.6	3	1500	560
104SX0301	PolySAX LP® Column, SS, 100 x 4.6mm, 3µm,100Å	SAX	Analytical	SS	10-32 coned	100	4.6	3	100	710
104SX0303	PolySAX LP® Column, SS, 100 x 4.6mm, 3µm,300Å	SAX	Analytical	SS	10-32 coned	100	4.6	3	300	710
104SX0310	PolySAX LP® Column, SS, 100 x 4.6mm, 3µm,1000Å	SAX	Analytical	SS	10-32 coned	100	4.6	3	1000	710
104SX0315	PolySAX LP® Column, SS, 100 x 4.6mm, 3µm,1500Å	SAX	Analytical	SS	10-32 coned	100	4.6	3	1500	710
154SX0301	PolySAX LP® Column, SS, 150 x 4.6mm, 3µm,100Å	SAX	Analytical	SS	10-32 coned	150	4.6	3	100	740
154SX0303	PolySAX LP® Column, SS, 150 x 4.6mm, 3µm,300Å	SAX	Analytical	SS	10-32 coned	150	4.6	3	300	740
154SX0310	PolySAX LP® Column, SS, 150 x 4.6mm, 3µm,1000Å	SAX	Analytical	SS	10-32 coned	150	4.6	3	1000	740
154SX0315	PolySAX LP® Column, SS, 150 x 4.6mm, 3µm,1500Å	SAX	Analytical	SS	10-32 coned	150	4.6	3	1500	740
204SX0301	PolySAX LP® Column, SS, 200 x 4.6mm, 3µm,100Å	SAX	Analytical	SS	10-32 coned	200	4.6	3	100	760
204SX0303	PolySAX LP® Column, SS, 200 x 4.6mm, 3µm,300Å	SAX	Analytical	SS	10-32 coned	200	4.6	3	300	760
204SX0310	PolySAX LP® Column, SS, 200 x 4.6mm, 3µm,1000Å	SAX	Analytical	SS	10-32 coned	200	4.6	3	1000	760
204SX0315	PolySAX LP® Column, SS, 200 x 4.6mm, 3µm,1500Å	SAX	Analytical	SS	10-32 coned	200	4.6	3	1500	760
254SX0301	PolySAX LP® Column, SS, 250 x 4.6mm, 3µm,100Å	SAX	Analytical	SS	10-32 coned	250	4.6	3	100	780
254SX0303	PolySAX LP® Column, SS, 250 x 4.6mm, 3µm,300Å	SAX	Analytical	SS	10-32 coned	250	4.6	3	300	780
254SX0310	PolySAX LP® Column, SS, 250 x 4.6mm, 3µm,1000Å	SAX	Analytical	SS	10-32 coned	250	4.6	3	1000	780
254SX0315	PolySAX LP® Column, SS, 250 x 4.6mm, 3µm,1500Å	SAX	Analytical	SS	10-32 coned	250	4.6	3	1500	780
109SX0301	PolySAX LP® Column, SS, 100 x 9.4mm, 3µm,100Å	SAX	Analytical	SS	10-32 coned	100	9.4	3	100	1400
109SX0303	PolySAX LP® Column, SS, 100 x 9.4mm, 3µm,300Å	SAX	Analytical	SS	10-32 coned	100	9.4	3	300	1400
109SX0310	PolySAX LP® Column, SS, 100 x 9.4mm, 3µm,1000Å	SAX	Analytical	SS	10-32 coned	100	9.4	3	1000	1400
109SX0315	PolySAX LP® Column, SS, 100 x 9.4mm, 3µm,1500Å	SAX	Analytical	SS	10-32 coned	100	9.4	3	1500	1400
209SX0301	PolySAX LP® Column, SS, 200 x 9.4mm, 3µm,100Å	SAX	Analytical	SS	10-32 coned	200	9.4	3	100	1620
209SX0303	PolySAX LP® Column, SS, 200 x 9.4mm, 3µm,300Å	SAX	Analytical	SS	10-32 coned	200	9.4	3	300	1620
209SX0310	PolySAX LP® Column, SS, 200 x 9.4mm, 3µm,1000Å	SAX	Analytical	SS	10-32 coned	200	9.4	3	1000	1620
209SX0315	PolySAX LP® Column, SS, 200 x 9.4mm, 3µm,1500Å	SAX	Analytical	SS	10-32 coned	200	9.4	3	1500	1620
259SX0301	PolySAX LP® Column, SS, 250 x 9.4mm, 3µm,100Å	SAX	Analytical	SS	10-32 coned	250	9.4	3	100	1750
259SX0303	PolySAX LP® Column, SS, 250 x 9.4mm, 3µm,300Å	SAX	Analytical	SS	10-32 coned	250	9.4	3	300	1750
259SX0310	PolySAX LP® Column, SS, 250 x 9.4mm, 3µm,1000Å	SAX	Analytical	SS	10-32 coned	250	9.4	3	1000	1750
259SX0315	PolySAX LP® Column, SS, 250 x 9.4mm, 3µm,1500Å	SAX	Analytical	SS	10-32 coned	250	9.4	3	1500	1750
051SX0501	PolySAX LP® Column, Glass Lined SS, 50 x 1.0mm, 5µm,100Å	SAX	Analytical	GLT	10-32 coned	50	1.0	5	100	570
051SX0503	PolySAX LP® Column, Glass Lined SS, 50 x 1.0mm, 5µm,300Å	SAX	Analytical	GLT	10-32 coned	50	1.0	5	300	570
051SX0510	PolySAX LP® Column, Glass Lined SS, 50 x 1.0mm, 5µm,1000Å	SAX	Analytical	GLT	10-32 coned	50	1.0	5	1000	570

PolySAX LP® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
101SX0501	PolySAX LP® Column, Glass Lined SS, 100 x 1.0mm, 5µm,100Å	SAX	Analytical	GLT	10-32 coned	100	1.0	5	100	660
101SX0503	PolySAX LP® Column, Glass Lined SS, 100 x 1.0mm, 5µm,300Å	SAX	Analytical	GLT	10-32 coned	100	1.0	5	300	660
101SX0510	PolySAX LP® Column, Glass Lined SS, 100 x 1.0mm, 5µm,1000Å	SAX	Analytical	GLT	10-32 coned	100	1.0	5	1000	660
151SX0501	PolySAX LP® Column, Glass Lined SS, 150 x 1.0mm, 5µm,100Å	SAX	Analytical	GLT	10-32 coned	150	1.0	5	100	740
151SX0503	PolySAX LP® Column, Glass Lined SS, 150 x 1.0mm, 5µm,300Å	SAX	Analytical	GLT	10-32 coned	150	1.0	5	300	740
151SX0510	PolySAX LP® Column, Glass Lined SS, 150 x 1.0mm, 5µm,1000Å	SAX	Analytical	GLT	10-32 coned	150	1.0	5	1000	740
3.52SX0501	PolySAX LP® Column, SS, 35 x 2.1mm, 5µm,100Å	SAX	Analytical	SS	10-32 coned	35	2.1	5	100	500
3.52SX0503	PolySAX LP® Column, SS, 35 x 2.1mm, 5µm,300Å	SAX	Analytical	SS	10-32 coned	35	2.1	5	300	500
3.52SX0510	PolySAX LP® Column, SS, 35 x 2.1mm, 5µm,1000Å	SAX	Analytical	SS	10-32 coned	35	2.1	5	1000	500
052SX0501	PolySAX LP® Column, SS, 50 x 2.1mm, 5µm,100Å	SAX	Analytical	SS	10-32 coned	50	2.1	5	100	520
052SX0503	PolySAX LP® Column, SS, 50 x 2.1mm, 5µm,300Å	SAX	Analytical	SS	10-32 coned	50	2.1	5	300	520
052SX0510	PolySAX LP® Column, SS, 50 x 2.1mm, 5µm,1000Å	SAX	Analytical	SS	10-32 coned	50	2.1	5	1000	520
102SX0501	PolySAX LP® Column, SS, 100 x 2.1mm, 5µm,100Å	SAX	Analytical	SS	10-32 coned	100	2.1	5	100	600
102SX0503	PolySAX LP® Column, SS, 100 x 2.1mm, 5µm,300Å	SAX	Analytical	SS	10-32 coned	100	2.1	5	300	600
102SX0510	PolySAX LP® Column, SS, 100 x 2.1mm, 5µm,1000Å	SAX	Analytical	SS	10-32 coned	100	2.1	5	1000	600
152SX0501	PolySAX LP® Column, SS, 150 x 2.1mm, 5µm,100Å	SAX	Analytical	SS	10-32 coned	150	2.1	5	100	650
152SX0503	PolySAX LP® Column, SS, 150 x 2.1mm, 5µm,300Å	SAX	Analytical	SS	10-32 coned	150	2.1	5	300	650
152SX0510	PolySAX LP® Column, SS, 150 x 2.1mm, 5µm,1000Å	SAX	Analytical	SS	10-32 coned	150	2.1	5	1000	650
202SX0501	PolySAX LP® Column, SS, 200 x 2.1mm, 5µm,100Å	SAX	Analytical	SS	10-32 coned	200	2.1	5	100	710
202SX0503	PolySAX LP® Column, SS, 200 x 2.1mm, 5µm,300Å	SAX	Analytical	SS	10-32 coned	200	2.1	5	300	710
202SX0510	PolySAX LP® Column, SS, 200 x 2.1mm, 5µm,1000Å	SAX	Analytical	SS	10-32 coned	200	2.1	5	1000	710
3.54SX0501	PolySAX LP® Column, SS, 35 x 4.6mm, 5µm,100Å	SAX	Analytical	SS	10-32 coned	35	4.6	5	100	500
3.54SX0503	PolySAX LP® Column, SS, 35 x 4.6mm, 5µm,300Å	SAX	Analytical	SS	10-32 coned	35	4.6	5	300	500
3.54SX0510	PolySAX LP® Column, SS, 35 x 4.6mm, 5µm,1000Å	SAX	Analytical	SS	10-32 coned	35	4.6	5	1000	500
054SX0501	PolySAX LP® Column, SS, 50 x 4.6mm, 5µm,100Å	SAX	Analytical	SS	10-32 coned	50	4.6	5	100	520
054SX0503	PolySAX LP® Column, SS, 50 x 4.6mm, 5µm,300Å	SAX	Analytical	SS	10-32 coned	50	4.6	5	300	520
054SX0510	PolySAX LP® Column, SS, 50 x 4.6mm, 5µm,1000Å	SAX	Analytical	SS	10-32 coned	50	4.6	5	1000	520
104SX0501	PolySAX LP® Column, SS, 100 x 4.6mm, 5µm,100Å	SAX	Analytical	SS	10-32 coned	100	4.6	5	100	600
104SX0503	PolySAX LP® Column, SS, 100 x 4.6mm, 5µm,300Å	SAX	Analytical	SS	10-32 coned	100	4.6	5	300	600
104SX0510	PolySAX LP® Column, SS, 100 x 4.6mm, 5µm,1000Å	SAX	Analytical	SS	10-32 coned	100	4.6	5	1000	600
154SX0501	PolySAX LP® Column, SS, 150 x 4.6mm, 5µm,100Å	SAX	Analytical	SS	10-32 coned	150	4.6	5	100	650
154SX0503	PolySAX LP® Column, SS, 150 x 4.6mm, 5µm,300Å	SAX	Analytical	SS	10-32 coned	150	4.6	5	300	650
154SX0510	PolySAX LP® Column, SS, 150 x 4.6mm, 5µm,1000Å	SAX	Analytical	SS	10-32 coned	150	4.6	5	1000	650
204SX0501	PolySAX LP® Column, SS, 200 x 4.6mm, 5µm,100Å	SAX	Analytical	SS	10-32 coned	200	4.6	5	100	710
204SX0503	PolySAX LP® Column, SS, 200 x 4.6mm, 5µm,300Å	SAX	Analytical	SS	10-32 coned	200	4.6	5	300	710
204SX0510	PolySAX LP® Column, SS, 200 x 4.6mm, 5µm,1000Å	SAX	Analytical	SS	10-32 coned	200	4.6	5	1000	710
254SX0501	PolySAX LP® Column, SS, 250 x 4.6mm, 5µm,100Å	SAX	Analytical	SS	10-32 coned	250	4.6	5	100	740
254SX0503	PolySAX LP® Column, SS, 250 x 4.6mm, 5µm,300Å	SAX	Analytical	SS	10-32 coned	250	4.6	5	300	740
254SX0510	PolySAX LP® Column, SS, 250 x 4.6mm, 5µm,1000Å	SAX	Analytical	SS	10-32 coned	250	4.6	5	1000	740
109SX0501	PolySAX LP® Column, SS, 100 x 9.4mm, 5µm,100Å	SAX	Semi-Prep	SS	10-32 coned	100	9.4	5	100	1110
109SX0503	PolySAX LP® Column, SS, 100 x 9.4mm, 5µm,300Å	SAX	Semi-Prep	SS	10-32 coned	100	9.4	5	300	1110
109SX0510	PolySAX LP® Column, SS, 100 x 9.4mm, 5µm,1000Å	SAX	Semi-Prep	SS	10-32 coned	100	9.4	5	1000	1110
209SX0501	PolySAX LP® Column, SS, 200 x 9.4mm, 5µm,100Å	SAX	Semi-Prep	SS	10-32 coned	200	9.4	5	100	1400
209SX0503	PolySAX LP® Column, SS, 200 x 9.4mm, 5µm,300Å	SAX	Semi-Prep	SS	10-32 coned	200	9.4	5	300	1400
209SX0510	PolySAX LP® Column, SS, 200 x 9.4mm, 5µm,1000Å	SAX	Semi-Prep	SS	10-32 coned	200	9.4	5	1000	1400

PolySAX LP® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
259SX0501	PolySAX LP® Column, SS, 250 x 9.4mm, 5µm,100Å	SAX	Semi-Prep	SS	10-32 coned	250	9.4	5	100	1500
259SX0503	PolySAX LP® Column, SS, 250 x 9.4mm, 5µm,300Å	SAX	Semi-Prep	SS	10-32 coned	250	9.4	5	300	1500
259SX0510	PolySAX LP® Column, SS, 250 x 9.4mm, 5µm,1000Å	SAX	Semi-Prep	SS	10-32 coned	250	9.4	5	1000	1500
2521SX0501	PolySAX LP® Column, SS, 250 x 21mm, 5µm,100Å	SAX	Semi-Prep	SS	10-32 coned	250	21.0	5	100	4500
2521SX0503	PolySAX LP® Column, SS, 250 x 21mm, 5µm,300Å	SAX	Semi-Prep	SS	10-32 coned	250	21.0	5	300	4500
2521SX0510	PolySAX LP® Column, SS, 250 x 21mm, 5µm,1000Å	SAX	Semi-Prep	SS	10-32 coned	250	21.0	5	1000	4500
204SX1201	PolySAX LP® Column, SS, 200 x 4.6mm, 12µm,100Å	SAX	Analytical	SS	10-32 coned	200	4.6	12	100	420
204SX1203	PolySAX LP® Column, SS, 200 x 4.6mm, 12µm,300Å	SAX	Analytical	SS	10-32 coned	200	4.6	12	300	420
204SX1215	PolySAX LP® Column, SS, 200 x 4.6mm, 12µm,1500Å	SAX	Analytical	SS	10-32 coned	200	4.6	12	1500	420
254SX1201	PolySAX LP® Column, SS, 250 x 4.6mm, 12µm,100Å	SAX	Analytical	SS	10-32 coned	250	4.6	12	100	500
254SX1203	PolySAX LP® Column, SS, 250 x 4.6mm, 12µm,300Å	SAX	Analytical	SS	10-32 coned	250	4.6	12	300	500
254SX1215	PolySAX LP® Column, SS, 250 x 4.6mm, 12µm,1500Å	SAX	Analytical	SS	10-32 coned	250	4.6	12	1500	500
209SX1201	PolySAX LP® Column, SS, 200 x 9.4mm, 12µm,100Å	SAX	Semi-Prep	SS	10-32 coned	200	9.4	12	100	1000
209SX1203	PolySAX LP® Column, SS, 200 x 9.4mm, 12µm,300Å	SAX	Semi-Prep	SS	10-32 coned	200	9.4	12	300	1000
209SX1215	PolySAX LP® Column, SS, 200 x 9.4mm, 12µm,1500Å	SAX	Semi-Prep	SS	10-32 coned	200	9.4	12	1500	1000
259SX1201	PolySAX LP® Column, SS, 250 x 9.4mm, 12µm,100Å	SAX	Semi-Prep	SS	10-32 coned	250	9.4	12	100	1100
259SX1203	PolySAX LP® Column, SS, 250 x 9.4mm, 12µm,300Å	SAX	Semi-Prep	SS	10-32 coned	250	9.4	12	300	1100
259SX1215	PolySAX LP® Column, SS, 250 x 9.4mm, 12µm,1500Å	SAX	Semi-Prep	SS	10-32 coned	250	9.4	12	1500	1100
2521SX1201	PolySAX LP® Column, SS, 250 x 21mm, 12µm,100Å	SAX	Semi-Prep	SS	10-32 coned	250	21.0	12	100	3200
2521SX1203	PolySAX LP® Column, SS, 250 x 21mm, 12µm,300Å	SAX	Semi-Prep	SS	10-32 coned	250	21.0	12	300	3200
2521SX1215	PolySAX LP® Column, SS, 250 x 21mm, 12µm,1500Å	SAX	Semi-Prep	SS	10-32 coned	250	21.0	12	1500	3200
2550SX1201	PolySAX LP® Column, SS, 250 x 50.8mm, 12µm,100Å	SAX	Prep	SS	10-32 coned	250	50.8	12	100	12250
2550SX1203	PolySAX LP® Column, SS, 250 x 50.8mm, 12µm,300Å	SAX	Prep	SS	10-32 coned	250	50.8	12	300	12250
2550SX1215	PolySAX LP® Column, SS, 250 x 50.8mm, 12µm,1500Å	SAX	Prep	SS	10-32 coned	250	50.8	12	1500	12250



HPLC Capillaries

PolySAX LP® Capillaries (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
050.075SX0210	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 2µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.075	2	1000	520
100.075SX0210	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 2µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.075	2	1000	580
150.075SX0210	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 2µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.075	2	1000	640
200.075SX0210	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 2µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.075	2	1000	700
050.15SX0210	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 2µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.150	2	1000	520
100.15SX0210	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 2µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.150	2	1000	580
150.15SX0210	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 2µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.150	2	1000	640
200.15SX0210	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 2µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.150	2	1000	700

PolySAX LP® Capillaries (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
050.30SX0210	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 2µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.300	2	1000	520
100.30SX0210	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 2µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.300	2	1000	580
150.30SX0210	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 2µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.300	2	1000	640
200.30SX0210	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 2µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.300	2	1000	700
050.075SX0301	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.075	3	100	520
050.075SX0303	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.075	3	300	520
050.075SX0310	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.075	3	1000	520
050.075SX0315	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,1500Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.075	3	1500	520
100.075SX0301	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.075	3	100	580
100.075SX0303	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.075	3	300	580
100.075SX0310	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.075	3	1000	580
100.075SX0315	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,1500Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.075	3	1500	580
150.075SX0301	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.075	3	100	640
150.075SX0303	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.075	3	300	640
150.075SX0310	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.075	3	1000	640
150.075SX0315	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,1500Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.075	3	1500	640
200.075SX0301	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.075	3	100	700
200.075SX0303	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.075	3	300	700
200.075SX0310	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.075	3	1000	700
200.075SX0315	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,1500Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.075	3	1500	700
050.15SX0301	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.150	3	100	520
050.15SX0303	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.150	3	300	520
050.15SX0310	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.150	3	1000	520
050.15SX0315	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,1500Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.150	3	1500	520
100.15SX0301	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.150	3	100	580
100.15SX0303	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.150	3	300	580
100.15SX0310	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.150	3	1000	580
100.15SX0315	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,1500Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.150	3	1500	580
150.15SX0301	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.150	3	100	640
150.15SX0303	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.150	3	300	640
150.15SX0310	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.150	3	1000	640
150.15SX0315	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,1500Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.150	3	1500	640
200.15SX0301	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.150	3	100	700
200.15SX0303	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.150	3	300	700
200.15SX0310	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.150	3	1000	700
200.15SX0315	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,1500Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.150	3	1500	700
050.30SX0301	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.300	3	100	520
050.30SX0303	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.300	3	300	520
050.30SX0310	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.300	3	1000	520
050.30SX0315	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,1500Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.300	3	1500	520
100.30SX0301	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.300	3	100	580
100.30SX0303	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.300	3	300	580
100.30SX0310	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.300	3	1000	580
100.30SX0315	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,1500Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.300	3	1500	580
150.30SX0301	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.300	3	100	640
150.30SX0303	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.300	3	300	640

PolySAX LP® Capillaries (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.) Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
150.30SX0310	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.300	3	1000	640
150.30SX0315	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,1500Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.300	3	1500	640
200.30SX0301	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.300	3	100	700
200.30SX0303	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.300	3	300	700
200.30SX0310	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.300	3	1000	700
200.30SX0315	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,1500Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.300	3	1500	700
050.075SX0501	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 5µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.075	5	100	520
050.075SX0503	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 5µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.075	5	300	520
050.075SX0510	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 5µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.075	5	1000	520
100.075SX0501	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 5µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.075	5	100	580
100.075SX0503	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 5µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.075	5	300	580
100.075SX0510	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 5µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.075	5	1000	580
150.075SX0501	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 5µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.075	5	100	640
150.075SX0503	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 5µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.075	5	300	640
150.075SX0510	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 5µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.075	5	1000	640
200.075SX0501	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 5µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.075	5	100	700
200.075SX0503	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 5µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.075	5	300	700
200.075SX0510	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 5µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.075	5	1000	700
050.15SX0501	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 5µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.150	5	100	520
050.15SX0503	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 5µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.150	5	300	520
050.15SX0510	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 5µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.150	5	1000	520
100.15SX0501	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 5µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.150	5	100	580
100.15SX0503	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 5µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.150	5	300	580
100.15SX0510	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 5µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.150	5	1000	580
150.15SX0501	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 5µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.150	5	100	640
150.15SX0503	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 5µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.150	5	300	640
150.15SX0510	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 5µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.150	5	1000	640
200.15SX0501	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 5µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.150	5	100	700
200.15SX0503	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 5µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.150	5	300	700
200.15SX0510	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 5µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.150	5	1000	700
050.30SX0501	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 5µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.300	5	100	520
050.30SX0503	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 5µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.300	5	300	520
050.30SX0510	PolySAX LP® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 5µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	50	0.300	5	1000	520
100.30SX0501	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 5µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.300	5	100	580
100.30SX0503	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 5µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.300	5	300	580
100.30SX0510	PolySAX LP® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 5µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	100	0.300	5	1000	580
150.30SX0501	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 5µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.300	5	100	640
150.30SX0503	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 5µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.300	5	300	640
150.30SX0510	PolySAX LP® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 5µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	150	0.300	5	1000	640
200.30SX0501	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 5µm,100Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.300	5	100	700
200.30SX0503	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 5µm,300Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.300	5	300	700
200.30SX0510	PolySAX LP® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 5µm,1000Å	SAX	Analytical	PEEK SIL	10-32 coned	200	0.300	5	1000	700

Guard Cartridges



PolySAX LP® Guard Cartridges (Available in a 2, 3, and 5µm particles with pore size options.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	ConneSXion	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
GC1-SX0210	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 1mm, 2µm, 1000Å	SAX	Analytical	SS/PEEK	10-32 coned	10	1.0	2	1000	120
GC2-SX0210	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 2µm, 1000Å	SAX	Analytical	SS/PEEK	10-32 coned	10	2.0	2	1000	120
GC3-SX0210	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 2µm, 1000Å	SAX	Analytical	SS/PEEK	10-32 coned	10	3.0	2	1000	120
GC4-SX0210	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 2µm, 1000Å	SAX	Analytical	SS/PEEK	10-32 coned	10	4.0	2	1000	120
GC1-SX0301	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 100Å	SAX	Analytical	SS/PEEK	10-32 coned	10	1.0	3	100	100
GC1-SX0303	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 300Å	SAX	Analytical	SS/PEEK	10-32 coned	10	1.0	3	300	100
GC1-SX0310	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 1000Å	SAX	Analytical	SS/PEEK	10-32 coned	10	1.0	3	1000	100
GC1-SX0315	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 1500Å	SAX	Analytical	SS/PEEK	10-32 coned	10	1.0	3	1500	100
GC2-SX0301	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 3µm, 100Å	SAX	Analytical	SS/PEEK	10-32 coned	10	2.0	3	100	100
GC2-SX0303	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 3µm, 300Å	SAX	Analytical	SS/PEEK	10-32 coned	10	2.0	3	300	100
GC2-SX0310	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 3µm, 1000Å	SAX	Analytical	SS/PEEK	10-32 coned	10	2.0	3	1000	100
GC2-SX0315	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 3µm, 1500Å	SAX	Analytical	SS/PEEK	10-32 coned	10	2.0	3	1500	100
GC3-SX0301	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 3µm, 100Å	SAX	Analytical	SS/PEEK	10-32 coned	10	3.0	3	100	100
GC3-SX0303	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 3µm, 300Å	SAX	Analytical	SS/PEEK	10-32 coned	10	3.0	3	300	100
GC3-SX0310	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 3µm, 1000Å	SAX	Analytical	SS/PEEK	10-32 coned	10	3.0	3	1000	100
GC3-SX0315	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 3µm, 1500Å	SAX	Analytical	SS/PEEK	10-32 coned	10	3.0	3	1500	100
GC4-SX0301	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 3µm, 100Å	SAX	Analytical	SS/PEEK	10-32 coned	10	4.0	3	100	100
GC4-SX0303	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 3µm, 300Å	SAX	Analytical	SS/PEEK	10-32 coned	10	4.0	3	300	100
GC4-SX0310	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 3µm, 1000Å	SAX	Analytical	SS/PEEK	10-32 coned	10	4.0	3	1000	100
GC4-SX0315	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 3µm, 1500Å	SAX	Analytical	SS/PEEK	10-32 coned	10	4.0	3	1500	100
GC1-SX0501	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 1mm, 5µm, 100Å	SAX	Analytical	SS/PEEK	10-32 coned	10	1.0	5	100	90
GC1-SX0503	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 1mm, 5µm, 300Å	SAX	Analytical	SS/PEEK	10-32 coned	10	1.0	5	300	90
GC1-SX0510	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 1mm, 5µm, 1000Å	SAX	Analytical	SS/PEEK	10-32 coned	10	1.0	5	1000	90
GC2-SX0501	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 5µm, 100Å	SAX	Analytical	SS/PEEK	10-32 coned	10	2.0	5	100	90
GC2-SX0503	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 5µm, 300Å	SAX	Analytical	SS/PEEK	10-32 coned	10	2.0	5	300	90
GC2-SX0510	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 5µm, 1000Å	SAX	Analytical	SS/PEEK	10-32 coned	10	2.0	5	1000	90
GC3-SX0501	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 5µm, 100Å	SAX	Analytical	SS/PEEK	10-32 coned	10	3.0	5	100	90
GC3-SX0503	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 5µm, 300Å	SAX	Analytical	SS/PEEK	10-32 coned	10	3.0	5	300	90
GC3-SX0510	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 5µm, 1000Å	SAX	Analytical	SS/PEEK	10-32 coned	10	3.0	5	1000	90
GC4-SX0501	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 5µm, 100Å	SAX	Analytical	SS/PEEK	10-32 coned	10	4.0	5	100	90
GC4-SX0503	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 5µm, 300Å	SAX	Analytical	SS/PEEK	10-32 coned	10	4.0	5	300	90
GC4-SX0510	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 5µm, 1000Å	SAX	Analytical	SS/PEEK	10-32 coned	10	4.0	5	1000	90
DGCH10	DireSX Guard Cartridge Holder, 10mm (For all analytical guards, REUSABLE)	N/A	Analytical	SS/PEEK	10-32 coned	10	N/A	N/A	N/A	75
IGCH10	Inline Guard Cartridge Holder, 10mm (For all analytical guards, REUSABLE)	N/A	Analytical	SS/PEEK	10-32 coned	10	N/A	N/A	N/A	75

*** NOTE: DGCH10 for DirectConnection to Column (female inlet/male outlet) IGCH10 for Connection to Line (two female ports)

GC10-SX0301	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 100Å	SAX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	100	225
GC10-SX0303	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 300Å	SAX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	300	225
GC10-SX0310	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 1000Å	SAX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	1000	225
GC10-SX0315	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 1500Å	SAX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	1500	225
GC10-SX0501	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 10mm, 5µm, 100Å	SAX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	5	100	185
GC10-SX0503	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 10mm, 5µm, 300Å	SAX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	5	300	185
GC10-SX0510	PolySAX LP® Guard Cartridge, SS/PEEK, 10mm x 10mm, 5µm, 1000Å	SAX	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	5	1000	185
CH10	Semi-Prep Guard Cartridge Holder, 10mm (For 9.4mm columns, REUSABLE)	SAX	Semi-Prep	SS/PEEK	10-32 coned	10	N/A	N/A	N/A	175

*** NOTE: CH10 Holder Required for Semi-Prep Guard Cartridges (10 x 10mm) Connection to Line (two female ports)

Solid Phase Extraction Cartridges

PolySAX LP® SPE Cartridges

PN	Product Description	Phase	Type	Hardware	Connexion	V (ml)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
SPESX1201	PolySAX LP® SPE Cartridge 1mL, PE, 12µm, 100Å (Pack of 10 cartridges)	SAX	SPE	PE	Lure	1	N/A	12	100	130
SPESX1203	PolySAX LP® SPE Cartridge 1mL, PE, 12µm, 300Å (Pack of 10 cartridges)	SAX	SPE	PE	Lure	1	N/A	12	300	130
SPESX1215	PolySAX LP® SPE Cartridge 1mL, PE, 12µm, 1500Å (Pack of 10 cartridges)	SAX	SPE	PE	Lure	1	N/A	12	1500	130
SPESX2003	PolySAX LP® SPE Cartridge 1mL, PE, 20µm, 300Å (Pack of 10 cartridges)	SAX	SPE	PE	Lure	1	N/A	20	300	130
SPESX2015	PolySAX LP® SPE Cartridge 1mL, PE, 20µm, 1500Å (Pack of 10 cartridges)	SAX	SPE	PE	Lure	1	N/A	20	1500	130

Custom combinations available upon request info@polylc.com

Spin Columns

PolySAX LP® BioPureSPN Columns

PN	Product Description	Phase	Type	Bed (mg)	Pore (Å)	SAX Capacity	HILIC Capacity	2026 Price (USD)
HMM HIL-SAX	PolySAX LP® BioPureSPN MACRO for SAX or HILIC 50-200µl (Box of 100 each)	SAX	Spin	100	300	0.2 -1mg	35 - 350µg	940
HMM HIL-SAX.20	PolySAX LP® BioPureSPN MACRO for SAX or HILIC 50-200µl (Pack of 20 each)	SAX	Spin	100	300	0.2 -1mg	35 - 350µg	265
HMM HIL-SAX-L	PolySAX LP® BioPureSPN MACRO PLATE 96 Well for SAX or HILIC 50-200µl	SAX	Spin	100	300	0.2 -1mg	35 - 350µg	800
HEM HIL-SAX	PolySAX LP® BioPureSPN MIDI for SAX/HILIC 50-200µl (Box of 100 each)	SAX	Spin	50	300	100-500µg	17 - 170µg	820
HEM HIL-SAX.20	PolySAX LP® BioPureSPN MIDI for SAX/HILIC 50-200µl (Pack of 20 each)	SAX	Spin	50	300	100-500µg	17 - 170µg	240
HNS HIL-SAX-M	PolySAX LP® BioPureSPN MIDI PLATE 96 Well for SAX/HILIC 50-200µl	SAX	Spin	50	300	100-500µg	17 - 170µg	575
HUM HIL-SAX	PolySAX LP® BioPureSPN MINI for SAX/HILIC 20-100µl (Box of 100 each)	SAX	Spin	20	300	40-200µg	7 - 70µg	735
HUM HIL-SAX.20	PolySAX LP® BioPureSPN MINI for SAX/HILIC 20-100µl (Pack of 20 each)	SAX	Spin	20	300	40-200µg	7 - 70µg	225
HNS HIL-SAX	PolySAX LP® BioPureSPN MINI PLATE 96 Well for SAX/HILIC 20-100µl	SAX	Spin	20	300	40-200µg	7 - 70µg	735

Custom combinations available upon request info@polylc.com

Bulk Material



PolySAX LP® Bulk Material (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.) Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	QTY (g)	Particle (µm)	Pore (Å)	2026 Price (USD)
BMSX0210	PolySAX LP® Bulk Material, 2µm, 1000Å	SAX	BM	N/A	0.5	2	1000	100
BMSX0301	PolySAX LP® Bulk Material, 3µm, 100Å	SAX	BM	N/A	1	3	100	130
BMSX0303	PolySAX LP® Bulk Material, 3µm, 300Å	SAX	BM	N/A	1	3	300	130
BMSX0310	PolySAX LP® Bulk Material, 3µm, 1000Å	SAX	BM	N/A	1	3	1000	130
BMSX0315	PolySAX LP® Bulk Material, 3µm, 1500Å	SAX	BM	N/A	1	3	1500	130
BMSX0501	PolySAX LP® Bulk Material, 5µm, 100Å	SAX	BM	N/A	1	5	100	85
BMSX0503	PolySAX LP® Bulk Material, 5µm, 300Å	SAX	BM	N/A	1	5	300	85
BMSX0510	PolySAX LP® Bulk Material, 5µm, 1000Å	SAX	BM	N/A	1	5	1000	85
BMSX1203	PolySAX LP® Bulk Material, 12µm, 300Å	SAX	BM	N/A	1	12	300	55
BMSX1215	PolySAX LP® Bulk Material, 12µm, 1000Å	SAX	BM	N/A	1	12	1000	55
BMSX2001	PolySAX LP® Bulk Material, 20µm, 100Å	SAX	BM	N/A	1	20	100	55
BMSX2003	PolySAX LP® Bulk Material, 20µm, 300Å	SAX	BM	N/A	1	20	300	55
BMSX2015	PolySAX LP® Bulk Material, 20µm, 1500Å	SAX	BM	N/A	1	20	1500	55
BMSX3003	PolySAX LP® Bulk Material, 30µm, 300Å	SAX	BM	N/A	1	30	300	55
BMSX3010	PolySAX LP® Bulk Material, 30µm, 1000Å	SAX	BM	N/A	1	30	1000	55

Please contact info@polylc.com with any questions or for information on any product combination not listed in the table above

PolySAX LP®

(Strong-Anion Exchange)

2026 PolyLC Inc. Retail Pricelist

PolyLC Inc., 9151 Rumsey Rd., Suite 180, Columbia, MD 21045 410-992-5400 www.polylc.com



2026 PolyLC Inc. Retail Pricelist

Hydrophobic Interaction Chromatography (HIC) Resins



Offering high-performance liquid chromatography media, columns, capillaries, cartridges, and instrumentation. We specialize in advanced alternatives to reversed-phase for proteins, peptides, and other biomolecules, supported by deep application expertise and customized solutions for the most demanding separations.

We offer a broad range of HIC stationary phases with ligand lengths spanning one to ten carbons. All HIC materials share the same dimensional-based pricing structure.

The products listed below are shown using our standard PolyPROPYL A® resin. To select a different ligand length, simply replace the material identifier within the existing part number with the appropriate ligand code.

Example: To convert listed PN 051PR0203 to a ten-carbon ligand, substitute DEC to obtain 051DEC0203 (PolyDECYL A®).

Please contact info@polylc.com if you have questions or need assistance selecting the appropriate material.

Length	HIC Material	Material Identifier
C1	PolyMETHYL A	ME
C2	PolyETHYL A	ET
C3	PolyPROPYL A	PR
C4	PolyBUTYL A	BU
C5	PolyPENTYL A	PE
C6	PolyHEXYL A	HX
C7	PolyHEPTYL A	HP
C8	PolyOCTYL A	OCT
C9	PolyNONYL A	NON
C10	PolyDECYL A	DEC

HPLC Columns

PolyPROPYL A® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
051PR0203	PolyPROPYL A® Column, Glass Lined SS, 50 x 1.0mm, 2µm,300Å	HIC	Analytical	GLT	10-32 coned	50	1.0	2	300	615
051PR0210	PolyPROPYL A® Column, Glass Lined SS, 50 x 1.0mm, 2µm,1000Å	HIC	Analytical	GLT	10-32 coned	50	1.0	2	1000	615
051PR0215	PolyPROPYL A® Column, Glass Lined SS, 50 x 1.0mm, 2µm,1500Å	HIC	Analytical	GLT	10-32 coned	50	1.0	2	1500	615
101PR0203	PolyPROPYL A® Column, Glass Lined SS, 100 x 1.0mm, 2µm,300Å	HIC	Analytical	GLT	10-32 coned	100	1.0	2	300	750
101PR0210	PolyPROPYL A® Column, Glass Lined SS, 100 x 1.0mm, 2µm,1000Å	HIC	Analytical	GLT	10-32 coned	100	1.0	2	1000	750
101PR0215	PolyPROPYL A® Column, Glass Lined SS, 100 x 1.0mm, 2µm,1500Å	HIC	Analytical	GLT	10-32 coned	100	1.0	2	1500	750
151PR0203	PolyPROPYL A® Column, Glass Lined SS, 150 x 1.0mm, 2µm,300Å	HIC	Analytical	GLT	10-32 coned	150	1.0	2	300	800
151PR0210	PolyPROPYL A® Column, Glass Lined SS, 150 x 1.0mm, 2µm,1000Å	HIC	Analytical	GLT	10-32 coned	150	1.0	2	1000	800
151PR0215	PolyPROPYL A® Column, Glass Lined SS, 150 x 1.0mm, 2µm,1500Å	HIC	Analytical	GLT	10-32 coned	150	1.0	2	1500	800
3.52PR0203	PolyPROPYL A® Column, SS, 35 x 2.1mm, 2µm,300Å	HIC	Analytical	SS	10-32 coned	35	2.1	2	300	570
3.52PR0210	PolyPROPYL A® Column, SS, 35 x 2.1mm, 2µm,1000Å	HIC	Analytical	SS	10-32 coned	35	2.1	2	1000	570
3.52PR0215	PolyPROPYL A® Column, SS, 35 x 2.1mm, 2µm,1500Å	HIC	Analytical	SS	10-32 coned	35	2.1	2	1500	570
052PR0203	PolyPROPYL A® Column, SS, 50 x 2.1mm, 2µm,300Å	HIC	Analytical	SS	10-32 coned	50	2.1	2	300	580
052PR0210	PolyPROPYL A® Column, SS, 50 x 2.1mm, 2µm,1000Å	HIC	Analytical	SS	10-32 coned	50	2.1	2	1000	580
052PR0215	PolyPROPYL A® Column, SS, 50 x 2.1mm, 2µm,1500Å	HIC	Analytical	SS	10-32 coned	50	2.1	2	1500	580
102PR0203	PolyPROPYL A® Column, SS, 100 x 2.1mm, 2µm,300Å	HIC	Analytical	SS	10-32 coned	100	2.1	2	300	750
102PR0210	PolyPROPYL A® Column, SS, 100 x 2.1mm, 2µm,1000Å	HIC	Analytical	SS	10-32 coned	100	2.1	2	1000	750
102PR0215	PolyPROPYL A® Column, SS, 100 x 2.1mm, 2µm,1500Å	HIC	Analytical	SS	10-32 coned	100	2.1	2	1500	750

PolyPROPYL A® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.)

 Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
3.54PR0203	PolyPROPYL A® Column, SS, 35 x 4.6mm, 2µm,300Å	HIC	Analytical	SS	10-32 coned	35	4.6	2	300	570
3.54PR0210	PolyPROPYL A® Column, SS, 35 x 4.6mm, 2µm,1000Å	HIC	Analytical	SS	10-32 coned	35	4.6	2	1000	570
3.54PR0215	PolyPROPYL A® Column, SS, 35 x 4.6mm, 2µm,1500Å	HIC	Analytical	SS	10-32 coned	35	4.6	2	1500	570
054PR0203	PolyPROPYL A® Column, SS, 50 x 4.6mm, 2µm,300Å	HIC	Analytical	SS	10-32 coned	50	4.6	2	300	580
054PR0210	PolyPROPYL A® Column, SS, 50 x 4.6mm, 2µm,1000Å	HIC	Analytical	SS	10-32 coned	50	4.6	2	1000	580
054PR0215	PolyPROPYL A® Column, SS, 50 x 4.6mm, 2µm,1500Å	HIC	Analytical	SS	10-32 coned	50	4.6	2	1500	580
104PR0203	PolyPROPYL A® Column, SS, 100 x 4.6mm, 2µm,300Å	HIC	Analytical	SS	10-32 coned	100	4.6	2	300	750
104PR0210	PolyPROPYL A® Column, SS, 100 x 4.6mm, 2µm,1000Å	HIC	Analytical	SS	10-32 coned	100	4.6	2	1000	750
104PR0215	PolyPROPYL A® Column, SS, 100 x 4.6mm, 2µm,1500Å	HIC	Analytical	SS	10-32 coned	100	4.6	2	1500	750
154PR0203	PolyPROPYL A® Column, SS, 150 x 4.6mm, 2µm,300Å	HIC	Analytical	SS	10-32 coned	150	4.6	2	300	790
154PR0210	PolyPROPYL A® Column, SS, 150 x 4.6mm, 2µm,1000Å	HIC	Analytical	SS	10-32 coned	150	4.6	2	1000	790
154PR0215	PolyPROPYL A® Column, SS, 150 x 4.6mm, 2µm,1500Å	HIC	Analytical	SS	10-32 coned	150	4.6	2	1500	790
204PR0203	PolyPROPYL A® Column, SS, 200 x 4.6mm, 2µm,300Å	HIC	Analytical	SS	10-32 coned	200	4.6	2	300	810
204PR0210	PolyPROPYL A® Column, SS, 200 x 4.6mm, 2µm,1000Å	HIC	Analytical	SS	10-32 coned	200	4.6	2	1000	810
204PR0215	PolyPROPYL A® Column, SS, 250 x 4.6mm, 2µm,1500Å	HIC	Analytical	SS	10-32 coned	200	4.6	2	1500	810
109PR0203	PolyPROPYL A® Column, SS, 100 x 9.4mm, 2µm,300Å	HIC	Semi-Prep	SS	10-32 coned	100	9.4	2	300	1525
109PR0210	PolyPROPYL A® Column, SS, 100 x 9.4mm, 2µm,1000Å	HIC	Semi-Prep	SS	10-32 coned	100	9.4	2	1000	1525
109PR0215	PolyPROPYL A® Column, SS, 100 x 9.4mm, 2µm,1500Å	HIC	Semi-Prep	SS	10-32 coned	100	9.4	2	1500	1525
209PR0203	PolyPROPYL A® Column, SS, 200 x 9.4mm, 2µm,300Å	HIC	Semi-Prep	SS	10-32 coned	200	9.4	2	300	2450
209PR0210	PolyPROPYL A® Column, SS, 200 x 9.4mm, 2µm,1000Å	HIC	Semi-Prep	SS	10-32 coned	200	9.4	2	1000	2450
209PR0215	PolyPROPYL A® Column, SS, 200 x 9.4mm, 2µm,1500Å	HIC	Semi-Prep	SS	10-32 coned	200	9.4	2	1500	2450
051PR0303	PolyPROPYL A® Column, Glass Lined SS, 50 x 1.0mm, 3µm,300Å	HIC	Analytical	GLT	10-32 coned	50	1.0	3	300	570
051PR0310	PolyPROPYL A® Column, Glass Lined SS, 50 x 1.0mm, 3µm,1000Å	HIC	Analytical	GLT	10-32 coned	50	1.0	3	1000	570
051PR0315	PolyPROPYL A® Column, Glass Lined SS, 50 x 1.0mm, 3µm,1500Å	HIC	Analytical	GLT	10-32 coned	50	1.0	3	1500	570
101PR0303	PolyPROPYL A® Column, Glass Lined SS, 100 x 1.0mm, 3µm,300Å	HIC	Analytical	GLT	10-32 coned	100	1.0	3	300	660
101PR0310	PolyPROPYL A® Column, Glass Lined SS, 100 x 1.0mm, 3µm,1000Å	HIC	Analytical	GLT	10-32 coned	100	1.0	3	1000	660
101PR0315	PolyPROPYL A® Column, Glass Lined SS, 100 x 1.0mm, 3µm,1500Å	HIC	Analytical	GLT	10-32 coned	100	1.0	3	1500	660
151PR0303	PolyPROPYL A® Column, Glass Lined SS, 150 x 1.0mm, 3µm,300Å	HIC	Analytical	GLT	10-32 coned	150	1.0	3	300	740
151PR0310	PolyPROPYL A® Column, Glass Lined SS, 150 x 1.0mm, 3µm,1000Å	HIC	Analytical	GLT	10-32 coned	150	1.0	3	1000	740
151PR0315	PolyPROPYL A® Column, Glass Lined SS, 150 x 1.0mm, 3µm,1500Å	HIC	Analytical	GLT	10-32 coned	150	1.0	3	1500	740
3.52PR0303	PolyPROPYL A® Column, SS, 35 x 2.1mm, 3µm,300Å	HIC	Analytical	SS	10-32 coned	35	2.1	3	300	545
3.52PR0310	PolyPROPYL A® Column, SS, 35 x 2.1mm, 3µm,1000Å	HIC	Analytical	SS	10-32 coned	35	2.1	3	1000	545
3.52PR0315	PolyPROPYL A® Column, SS, 35 x 2.1mm, 3µm,1500Å	HIC	Analytical	SS	10-32 coned	35	2.1	3	1500	545
052PR0303	PolyPROPYL A® Column, SS, 50 x 2.1mm, 3µm,300Å	HIC	Analytical	SS	10-32 coned	50	2.1	3	300	550
052PR0310	PolyPROPYL A® Column, SS, 50 x 2.1mm, 3µm,1000Å	HIC	Analytical	SS	10-32 coned	50	2.1	3	1000	550
052PR0315	PolyPROPYL A® Column, SS, 50 x 2.1mm, 3µm,1500Å	HIC	Analytical	SS	10-32 coned	50	2.1	3	1500	550
102PR0303	PolyPROPYL A® Column, SS, 100 x 2.1mm, 3µm,300Å	HIC	Analytical	SS	10-32 coned	100	2.1	3	300	710
102PR0310	PolyPROPYL A® Column, SS, 100 x 2.1mm, 3µm,1000Å	HIC	Analytical	SS	10-32 coned	100	2.1	3	1000	710
102PR0315	PolyPROPYL A® Column, SS, 100 x 2.1mm, 3µm,1500Å	HIC	Analytical	SS	10-32 coned	100	2.1	3	1500	710
3.54PR0303	PolyPROPYL A® Column, SS, 35 x 4.6mm, 3µm,300Å	HIC	Analytical	SS	10-32 coned	35	4.6	3	300	545
3.54PR0310	PolyPROPYL A® Column, SS, 35 x 4.6mm, 3µm,1000Å	HIC	Analytical	SS	10-32 coned	35	4.6	3	1000	545
3.54PR0315	PolyPROPYL A® Column, SS, 35 x 4.6mm, 3µm,1500Å	HIC	Analytical	SS	10-32 coned	35	4.6	3	1500	545
054PR0303	PolyPROPYL A® Column, SS, 50 x 4.6mm, 3µm,300Å	HIC	Analytical	SS	10-32 coned	50	4.6	3	300	550
054PR0310	PolyPROPYL A® Column, SS, 50 x 4.6mm, 3µm,1000Å	HIC	Analytical	SS	10-32 coned	50	4.6	3	1000	550
054PR0315	PolyPROPYL A® Column, SS, 50 x 4.6mm, 3µm,1500Å	HIC	Analytical	SS	10-32 coned	50	4.6	3	1500	550
104PR0303	PolyPROPYL A® Column, SS, 100 x 4.6mm, 3µm,300Å	HIC	Analytical	SS	10-32 coned	100	4.6	3	300	710
104PR0310	PolyPROPYL A® Column, SS, 100 x 4.6mm, 3µm,1000Å	HIC	Analytical	SS	10-32 coned	100	4.6	3	1000	710

PolyPROPYL A® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.)

 Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
104PR0315	PolyPROPYL A® Column, SS, 100 x 4.6mm, 3µm,1500Å	HIC	Analytical	SS	10-32 coned	100	4.6	3	1500	710
154PR0303	PolyPROPYL A® Column, SS, 150 x 4.6mm, 3µm,300Å	HIC	Analytical	SS	10-32 coned	150	4.6	3	300	740
154PR0310	PolyPROPYL A® Column, SS, 150 x 4.6mm, 3µm,1000Å	HIC	Analytical	SS	10-32 coned	150	4.6	3	1000	740
154PR0315	PolyPROPYL A® Column, SS, 150 x 4.6mm, 3µm,1500Å	HIC	Analytical	SS	10-32 coned	150	4.6	3	1500	740
204PR0303	PolyPROPYL A® Column, SS, 200 x 4.6mm, 3µm,300Å	HIC	Analytical	SS	10-32 coned	200	4.6	3	300	750
204PR0310	PolyPROPYL A® Column, SS, 200 x 4.6mm, 3µm,1000Å	HIC	Analytical	SS	10-32 coned	200	4.6	3	1000	750
204PR0315	PolyPROPYL A® Column, SS, 200 x 4.6mm, 3µm,1500Å	HIC	Analytical	SS	10-32 coned	200	4.6	3	1500	750
254PR0303	PolyPROPYL A® Column, SS, 250 x 4.6mm, 3µm,300Å	HIC	Analytical	SS	10-32 coned	250	4.6	3	300	800
254PR0310	PolyPROPYL A® Column, SS, 250 x 4.6mm, 3µm,1000Å	HIC	Analytical	SS	10-32 coned	250	4.6	3	1000	800
254PR0315	PolyPROPYL A® Column, SS, 250 x 4.6mm, 3µm,1500Å	HIC	Analytical	SS	10-32 coned	250	4.6	3	1500	800
109PR0303	PolyPROPYL A® Column, SS, 100 x 9.4mm, 3µm,300Å	HIC	Analytical	SS	10-32 coned	100	9.4	3	300	1400
109PR0310	PolyPROPYL A® Column, SS, 100 x 9.4mm, 3µm,1000Å	HIC	Analytical	SS	10-32 coned	100	9.4	3	1000	1400
109PR0315	PolyPROPYL A® Column, SS, 100 x 9.4mm, 3µm,1500Å	HIC	Analytical	SS	10-32 coned	100	9.4	3	1500	1400
209PR0303	PolyPROPYL A® Column, SS, 200 x 9.4mm, 3µm,300Å	HIC	Analytical	SS	10-32 coned	200	9.4	3	300	1620
209PR0310	PolyPROPYL A® Column, SS, 200 x 9.4mm, 3µm,1000Å	HIC	Analytical	SS	10-32 coned	200	9.4	3	1000	1620
209PR0315	PolyPROPYL A® Column, SS, 200 x 9.4mm, 3µm,1500Å	HIC	Analytical	SS	10-32 coned	200	9.4	3	1500	1620
259PR0303	PolyPROPYL A® Column, SS, 250 x 9.4mm, 3µm,300Å	HIC	Analytical	SS	10-32 coned	250	9.4	3	300	1750
259PR0310	PolyPROPYL A® Column, SS, 250 x 9.4mm, 3µm,1000Å	HIC	Analytical	SS	10-32 coned	250	9.4	3	1000	1750
259PR0315	PolyPROPYL A® Column, SS, 250 x 9.4mm, 3µm,1500Å	HIC	Analytical	SS	10-32 coned	250	9.4	3	1500	1750
051PR0503	PolyPROPYL A® Column, Glass Lined SS, 50 x 1.0mm, 5µm,300Å	HIC	Analytical	GLT	10-32 coned	50	1.0	5	300	570
051PR0510	PolyPROPYL A® Column, Glass Lined SS, 50 x 1.0mm, 5µm,1000Å	HIC	Analytical	GLT	10-32 coned	50	1.0	5	1000	570
101PR0503	PolyPROPYL A® Column, Glass Lined SS, 100 x 1.0mm, 5µm,300Å	HIC	Analytical	GLT	10-32 coned	100	1.0	5	300	660
101PR0510	PolyPROPYL A® Column, Glass Lined SS, 100 x 1.0mm, 5µm,1000Å	HIC	Analytical	GLT	10-32 coned	100	1.0	5	1000	660
151PR0503	PolyPROPYL A® Column, Glass Lined SS, 150 x 1.0mm, 5µm,300Å	HIC	Analytical	GLT	10-32 coned	150	1.0	5	300	740
151PR0510	PolyPROPYL A® Column, Glass Lined SS, 150 x 1.0mm, 5µm,1000Å	HIC	Analytical	GLT	10-32 coned	150	1.0	5	1000	740
3.52PR0503	PolyPROPYL A® Column, SS, 35 x 2.1mm, 5µm,300Å	HIC	Analytical	SS	10-32 coned	35	2.1	5	300	500
3.52PR0510	PolyPROPYL A® Column, SS, 35 x 2.1mm, 5µm,1000Å	HIC	Analytical	SS	10-32 coned	35	2.1	5	1000	500
052PR0503	PolyPROPYL A® Column, SS, 50 x 2.1mm, 5µm,300Å	HIC	Analytical	SS	10-32 coned	50	2.1	5	300	520
052PR0510	PolyPROPYL A® Column, SS, 50 x 2.1mm, 5µm,1000Å	HIC	Analytical	SS	10-32 coned	50	2.1	5	1000	520
102PR0503	PolyPROPYL A® Column, SS, 100 x 2.1mm, 5µm,300Å	HIC	Analytical	SS	10-32 coned	100	2.1	5	300	600
102PR0510	PolyPROPYL A® Column, SS, 100 x 2.1mm, 5µm,1000Å	HIC	Analytical	SS	10-32 coned	100	2.1	5	1000	600
152PR0503	PolyPROPYL A® Column, SS, 150 x 2.1mm, 5µm,300Å	HIC	Analytical	SS	10-32 coned	150	2.1	5	300	650
152PR0510	PolyPROPYL A® Column, SS, 150 x 2.1mm, 5µm,1000Å	HIC	Analytical	SS	10-32 coned	150	2.1	5	1000	650
202PR0503	PolyPROPYL A® Column, SS, 200 x 2.1mm, 5µm,300Å	HIC	Analytical	SS	10-32 coned	200	2.1	5	300	710
202PR0510	PolyPROPYL A® Column, SS, 200 x 2.1mm, 5µm,1000Å	HIC	Analytical	SS	10-32 coned	200	2.1	5	1000	710
3.54PR0503	PolyPROPYL A® Column, SS, 35 x 4.6mm, 5µm,300Å	HIC	Analytical	SS	10-32 coned	35	4.6	5	300	500
3.54PR0510	PolyPROPYL A® Column, SS, 35 x 4.6mm, 5µm,1000Å	HIC	Analytical	SS	10-32 coned	35	4.6	5	1000	500
054PR0503	PolyPROPYL A® Column, SS, 50 x 4.6mm, 5µm,300Å	HIC	Analytical	SS	10-32 coned	50	4.6	5	300	520
054PR0510	PolyPROPYL A® Column, SS, 50 x 4.6mm, 5µm,1000Å	HIC	Analytical	SS	10-32 coned	50	4.6	5	1000	520
104PR0503	PolyPROPYL A® Column, SS, 100 x 4.6mm, 5µm,300Å	HIC	Analytical	SS	10-32 coned	100	4.6	5	300	600
104PR0510	PolyPROPYL A® Column, SS, 100 x 4.6mm, 5µm,1000Å	HIC	Analytical	SS	10-32 coned	100	4.6	5	1000	600
154PR0503	PolyPROPYL A® Column, SS, 150 x 4.6mm, 5µm,300Å	HIC	Analytical	SS	10-32 coned	150	4.6	5	300	650
154PR0510	PolyPROPYL A® Column, SS, 150 x 4.6mm, 5µm,1000Å	HIC	Analytical	SS	10-32 coned	150	4.6	5	1000	650
204PR0503	PolyPROPYL A® Column, SS, 200 x 4.6mm, 5µm,300Å	HIC	Analytical	SS	10-32 coned	200	4.6	5	300	710
204PR0510	PolyPROPYL A® Column, SS, 200 x 4.6mm, 5µm,1000Å	HIC	Analytical	SS	10-32 coned	200	4.6	5	1000	710
254PR0503	PolyPROPYL A® Column, SS, 250 x 4.6mm, 5µm,300Å	HIC	Analytical	SS	10-32 coned	250	4.6	5	300	740
254PR0510	PolyPROPYL A® Column, SS, 250 x 4.6mm, 5µm,1000Å	HIC	Analytical	SS	10-32 coned	250	4.6	5	1000	740

PolyPROPYL A® Columns (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
109PR0503	PolyPROPYL A® Column, SS, 100 x 9.4mm, 5µm,300Å	HIC	Semi-Prep	SS	10-32 coned	100	9.4	5	300	1110
109PR0510	PolyPROPYL A® Column, SS, 100 x 9.4mm, 5µm,1000Å	HIC	Semi-Prep	SS	10-32 coned	100	9.4	5	1000	1110
209PR0503	PolyPROPYL A® Column, SS, 200 x 9.4mm, 5µm,300Å	HIC	Semi-Prep	SS	10-32 coned	200	9.4	5	300	1400
209PR0510	PolyPROPYL A® Column, SS, 200 x 9.4mm, 5µm,1000Å	HIC	Semi-Prep	SS	10-32 coned	200	9.4	5	1000	1400
259PR0503	PolyPROPYL A® Column, SS, 250 x 9.4mm, 5µm,300Å	HIC	Semi-Prep	SS	10-32 coned	250	9.4	5	300	1500
259PR0510	PolyPROPYL A® Column, SS, 250 x 9.4mm, 5µm,1000Å	HIC	Semi-Prep	SS	10-32 coned	250	9.4	5	1000	1500
2521PR0503	PolyPROPYL A® Column, SS, 250 x 21mm, 5µm,300Å	HIC	Semi-Prep	SS	10-32 coned	250	21.0	5	300	4500
2521PR0510	PolyPROPYL A® Column, SS, 250 x 21mm, 5µm,1000Å	HIC	Semi-Prep	SS	10-32 coned	250	21.0	5	1000	4500
204PR1203	PolyPROPYL A® Column, SS, 200 x 4.6mm, 12µm,300Å	HIC	Analytical	SS	10-32 coned	200	4.6	12	300	420
204PR1215	PolyPROPYL A® Column, SS, 200 x 4.6mm, 12µm,1500Å	HIC	Analytical	SS	10-32 coned	200	4.6	12	1500	420
254PR1203	PolyPROPYL A® Column, SS, 250 x 4.6mm, 12µm,300Å	HIC	Analytical	SS	10-32 coned	250	4.6	12	300	500
254PR1215	PolyPROPYL A® Column, SS, 250 x 4.6mm, 12µm,1500Å	HIC	Analytical	SS	10-32 coned	250	4.6	12	1500	500
209PR1203	PolyPROPYL A® Column, SS, 200 x 9.4mm, 12µm,300Å	HIC	Semi-Prep	SS	10-32 coned	200	9.4	12	300	1000
209PR1215	PolyPROPYL A® Column, SS, 200 x 9.4mm, 12µm,1500Å	HIC	Semi-Prep	SS	10-32 coned	200	9.4	12	1500	1000
259PR1203	PolyPROPYL A® Column, SS, 250 x 9.4mm, 12µm,300Å	HIC	Semi-Prep	SS	10-32 coned	250	9.4	12	300	1100
259PR1215	PolyPROPYL A® Column, SS, 250 x 9.4mm, 12µm,1500Å	HIC	Semi-Prep	SS	10-32 coned	250	9.4	12	1500	1100
2521PR1203	PolyPROPYL A® Column, SS, 250 x 21mm, 12µm,300Å	HIC	Semi-Prep	SS	10-32 coned	250	21.0	12	300	3200
2521PR1215	PolyPROPYL A® Column, SS, 250 x 21mm, 12µm,1500Å	HIC	Semi-Prep	SS	10-32 coned	250	21.0	12	1500	3200
2550PR1203	PolyPROPYL A® Column, SS, 250 x 50.8mm, 12µm,300Å	HIC	Prep	SS	10-32 coned	250	50.8	12	300	12250
2550PR1215	PolyPROPYL A® Column, SS, 250 x 50.8mm, 12µm,1500Å	HIC	Prep	SS	10-32 coned	250	50.8	12	1500	12250



HPLC Capillaries

PolyPROPYL A® Capillaries (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
050.075PR0203	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 2µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.075	2	300	520
050.075PR0210	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 2µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.075	2	1000	520
050.075PR0215	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 2µm,1500Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.075	2	1500	520
100.075PR0203	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 2µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.075	2	300	580
100.075PR0210	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 2µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.075	2	1000	580
100.075PR0215	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 2µm,1500Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.075	2	1500	580
150.075PR0203	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 2µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.075	2	300	640
150.075PR0210	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 2µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.075	2	1000	640
150.075PR0215	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 2µm,1500Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.075	2	1500	640
200.075PR0203	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 2µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.075	2	300	700
200.075PR0210	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 2µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.075	2	1000	700
200.075PR0215	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 2µm,1500Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.075	2	1500	700
050.15PR0203	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 2µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.150	2	300	520
050.15PR0210	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 2µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.150	2	1000	520
050.15PR0215	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 2µm,1500Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.150	2	1500	520

PolyPROPYL A® Capillaries (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
100.15PR0203	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 2µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.150	2	300	580
100.15PR0210	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 2µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.150	2	1000	580
100.15PR0215	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 2µm,1500Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.150	2	1500	580
150.15PR0203	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 2µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.150	2	300	640
150.15PR0210	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 2µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.150	2	1000	640
150.15PR0215	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 2µm,1500Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.150	2	1500	640
200.15PR0203	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 2µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.150	2	300	700
200.15PR0210	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 2µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.150	2	1000	700
200.15PR0215	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 2µm,1500Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.150	2	1500	700
050.30PR0203	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 2µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.300	2	300	520
050.30PR0210	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 2µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.300	2	1000	520
050.30PR0215	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 2µm,1500Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.300	2	1500	520
100.30PR0203	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 2µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.300	2	300	580
100.30PR0210	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 2µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.300	2	1000	580
100.30PR0215	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 2µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.300	2	1500	580
150.30PR0203	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 2µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.300	2	300	640
150.30PR0210	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 2µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.300	2	1000	640
150.30PR0215	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 2µm,1500Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.300	2	1500	640
200.30PR0203	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 2µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.300	2	300	700
200.30PR0210	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 2µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.300	2	1000	700
200.30PR0215	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 2µm,1500Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.300	2	1500	700
050.075PR0303	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.075	3	300	520
050.075PR0310	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.075	3	1000	520
050.075PR0315	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,1500Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.075	3	1500	520
100.075PR0303	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.075	3	300	580
100.075PR0310	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.075	3	1000	580
100.075PR0315	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,1500Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.075	3	1500	580
150.075PR0303	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.075	3	300	640
150.075PR0310	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.075	3	1000	640
150.075PR0315	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,1500Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.075	3	1500	640
200.075PR0303	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.075	3	300	700
200.075PR0310	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.075	3	1000	700
200.075PR0315	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,1500Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.075	3	1500	700
050.15PR0303	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.150	3	300	520
050.15PR0310	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.150	3	1000	520
050.15PR0315	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,1500Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.150	3	1500	520
100.15PR0303	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.150	3	300	580
100.15PR0310	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.150	3	1000	580
100.15PR0315	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,1500Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.150	3	1500	580
150.15PR0303	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.150	3	300	640
150.15PR0310	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.150	3	1000	640
150.15PR0315	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,1500Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.150	3	1500	640
200.15PR0303	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.150	3	300	700
200.15PR0310	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.150	3	1000	700
200.15PR0315	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,1500Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.150	3	1500	700
050.30PR0303	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.300	3	300	520
050.30PR0310	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.300	3	1000	520

PolyPROPYL A® Capillaries (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
050.30PR0315	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,1500Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.300	3	1500	520
100.30PR0303	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.300	3	300	580
100.30PR0310	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.300	3	1000	580
100.30PR0315	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,1500Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.300	3	1500	580
150.30PR0303	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.300	3	300	640
150.30PR0310	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.300	3	1000	640
150.30PR0315	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,1500Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.300	3	1500	640
200.30PR0303	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.300	3	300	700
200.30PR0310	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.300	3	1000	700
200.30PR0315	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,1500Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.300	3	1500	700
050.075PR0503	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 5µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.075	5	300	520
050.075PR0510	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 5µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.075	5	1000	520
100.075PR0503	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 5µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.075	5	300	580
100.075PR0510	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 5µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.075	5	1000	580
150.075PR0503	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 5µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.075	5	300	640
150.075PR0510	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 5µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.075	5	1000	640
200.075PR0503	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 5µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.075	5	300	700
200.075PR0510	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 5µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.075	5	1000	700
050.15PR0503	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 5µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.150	5	300	520
050.15PR0510	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 5µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.150	5	1000	520
100.15PR0503	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 5µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.150	5	300	580
100.15PR0510	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 5µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.150	5	1000	580
150.15PR0503	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 5µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.150	5	300	640
150.15PR0510	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 5µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.150	5	1000	640
200.15PR0503	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 5µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.150	5	300	700
200.15PR0510	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 5µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.150	5	1000	700
050.30PR0503	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 5µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.300	5	300	520
050.30PR0510	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 5µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	50	0.300	5	1000	520
100.30PR0503	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 5µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.300	5	300	580
100.30PR0510	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 5µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	100	0.300	5	1000	580
150.30PR0503	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 5µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.300	5	300	640
150.30PR0510	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 5µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	150	0.300	5	1000	640
200.30PR0503	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 5µm,300Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.300	5	300	700
200.30PR0510	PolyPROPYL A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 5µm,1000Å	HIC	Analytical	PEEK SIL	10-32 coned	200	0.300	5	1000	700

Guard Cartridges



PolyPROPYL A® Guard Cartridges (Available in a 2, 3, and 5µm particles with pore size options.)

Custom combinations available upon request

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
GC1-PR0203	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 2µm, 300Å	HIC	Analytical	SS/PEEK	10-32 coned	10	1.0	2	300	120
GC1-PR0210	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 2µm, 1000Å	HIC	Analytical	SS/PEEK	10-32 coned	10	1.0	2	1000	120
GC1-PR0215	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 2µm, 1500Å	HIC	Analytical	SS/PEEK	10-32 coned	10	1.0	2	1500	120
GC2-PR0203	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 2µm, 300Å	HIC	Analytical	SS/PEEK	10-32 coned	10	2.0	2	300	120
GC2-PR0210	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 2µm, 1000Å	HIC	Analytical	SS/PEEK	10-32 coned	10	2.0	2	1000	120
GC2-PR0215	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 2µm, 1500Å	HIC	Analytical	SS/PEEK	10-32 coned	10	2.0	2	1500	120
GC3-PR0203	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 2µm, 300Å	HIC	Analytical	SS/PEEK	10-32 coned	10	3.0	2	300	120
GC3-PR0210	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 2µm, 1000Å	HIC	Analytical	SS/PEEK	10-32 coned	10	3.0	2	1000	120
GC3-PR0215	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 2µm, 1500Å	HIC	Analytical	SS/PEEK	10-32 coned	10	3.0	2	1500	120
GC4-PR0203	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 2µm, 300Å	HIC	Analytical	SS/PEEK	10-32 coned	10	4.0	2	300	120
GC4-PR0210	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 2µm, 1000Å	HIC	Analytical	SS/PEEK	10-32 coned	10	4.0	2	1000	120
GC4-PR0215	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 2µm, 1500Å	HIC	Analytical	SS/PEEK	10-32 coned	10	4.0	2	1500	120
GC1-PR0303	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 300Å	HIC	Analytical	SS/PEEK	10-32 coned	10	1.0	3	300	100
GC1-PR0310	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 1000Å	HIC	Analytical	SS/PEEK	10-32 coned	10	1.0	3	1000	100
GC1-PR0315	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 3µm, 1500Å	HIC	Analytical	SS/PEEK	10-32 coned	10	1.0	3	1500	100
GC2-PR0303	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 3µm, 300Å	HIC	Analytical	SS/PEEK	10-32 coned	10	2.0	3	300	100
GC2-PR0310	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 3µm, 1000Å	HIC	Analytical	SS/PEEK	10-32 coned	10	2.0	3	1000	100
GC2-PR0315	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 3µm, 1500Å	HIC	Analytical	SS/PEEK	10-32 coned	10	2.0	3	1500	100
GC3-PR0303	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 3µm, 300Å	HIC	Analytical	SS/PEEK	10-32 coned	10	3.0	3	300	100
GC3-PR0310	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 3µm, 1000Å	HIC	Analytical	SS/PEEK	10-32 coned	10	3.0	3	1000	100
GC3-PR0315	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 3µm, 1500Å	HIC	Analytical	SS/PEEK	10-32 coned	10	3.0	3	1500	100
GC4-PR0301	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 3µm, 100Å	HIC	Analytical	SS/PEEK	10-32 coned	10	4.0	3	100	100
GC4-PR0303	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 3µm, 300Å	HIC	Analytical	SS/PEEK	10-32 coned	10	4.0	3	300	100
GC4-PR0310	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 3µm, 1000Å	HIC	Analytical	SS/PEEK	10-32 coned	10	4.0	3	1000	100
GC4-PR0315	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 3µm, 1500Å	HIC	Analytical	SS/PEEK	10-32 coned	10	4.0	3	1500	100
GC1-PR0503	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 5µm, 300Å	HIC	Analytical	SS/PEEK	10-32 coned	10	1.0	5	300	90
GC1-PR0510	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 1mm, 5µm, 1000Å	HIC	Analytical	SS/PEEK	10-32 coned	10	1.0	5	1000	90
GC2-PR0503	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 5µm, 300Å	HIC	Analytical	SS/PEEK	10-32 coned	10	2.0	5	300	90
GC2-PR0510	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 2.0mm, 5µm, 1000Å	HIC	Analytical	SS/PEEK	10-32 coned	10	2.0	5	1000	90
GC3-PR0503	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 5µm, 300Å	HIC	Analytical	SS/PEEK	10-32 coned	10	3.0	5	300	90
GC3-PR0510	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 3.0mm, 5µm, 1000Å	HIC	Analytical	SS/PEEK	10-32 coned	10	3.0	5	1000	90
GC4-PR0503	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 5µm, 300Å	HIC	Analytical	SS/PEEK	10-32 coned	10	4.0	5	300	90
GC4-PR0510	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 4.0mm, 5µm, 1000Å	HIC	Analytical	SS/PEEK	10-32 coned	10	4.0	5	1000	90
DGCH10	DirePR Guard Cartridge Holder, 10mm (For all analytical guards, REUSABLE)	N/A	Analytical	SS/PEEK	10-32 coned	10	N/A	N/A	N/A	75
IGCH10	Inline Guard Cartridge Holder, 10mm (For all analytical guards, REUSABLE)	N/A	Analytical	SS/PEEK	10-32 coned	10	N/A	N/A	N/A	75

*** NOTE: DGCH10 for DirectConnection to Column (female inlet/male outlet) IGCH10 for Connection to Line (two female ports)

GC10-PR0303	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 300Å	HIC	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	300	225
GC10-PR0310	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 1000Å	HIC	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	1000	225
GC10-PR0315	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 3µm, 1500Å	HIC	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	3	1500	225
GC10-PR0503	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 5µm, 300Å	HIC	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	5	300	185
GC10-PR0510	PolyPROPYL A® Guard Cartridge, SS/PEEK, 10mm x 10mm, 5µm, 1000Å	HIC	Semi-Prep	SS/PEEK	10-32 coned	10	10.0	5	1000	185
CH10	Semi-Prep Guard Cartridge Holder, 10mm (For 9.4mm columns, REUSABLE)	HIC	Semi-Prep	SS/PEEK	10-32 coned	10	N/A	N/A	N/A	175

*** NOTE: CH10 Holder Required for Semi-Prep Guard Cartridges (10 x 10mm) Connection to Line (two female ports)

Solid Phase Extraction Cartridges

PolyPROPYL A® SPE Cartridges

PN	Product Description	Phase	Type	Hardware	Connection	V (ml)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
SPEPR1203	PolyPROPYL A® SPE Cartridge 1mL, PE, 12µm, 300Å (Pack of 10 cartridges)	HIC	SPE	PE	Lure	1	N/A	12	300	130
SPEPR1215	PolyPROPYL A® SPE Cartridge 1mL, PE, 12µm, 1500Å (Pack of 10 cartridges)	HIC	SPE	PE	Lure	1	N/A	12	1500	130
SPEPR2003	PolyPROPYL A® SPE Cartridge 1mL, PE, 20µm, 300Å (Pack of 10 cartridges)	HIC	SPE	PE	Lure	1	N/A	20	300	130
SPEPR2015	PolyPROPYL A® SPE Cartridge 1mL, PE, 20µm, 1500Å (Pack of 10 cartridges)	HIC	SPE	PE	Lure	1	N/A	20	1500	130


Spin Columns

PolyPROPYL A® BioPureSPN Columns

PN	Product Description	Phase	Type	Bed (mg)	Pore (Å)	HIC Capacity	HILIC Capacity	2026 Price (USD)
HMM HIC	PolyPROPYL A® BioPureSPN MACRO for HIC 50-200µl (Box of 100 each)	HIC	Spin	100	300	0.2 -1mg	N/A	940
HMM HIC.20	PolyPROPYL A® BioPureSPN MACRO for HIC 50-200µl (Pack of 20 each)	HIC	Spin	100	300	0.2 -1mg	N/A	265
HMM HIC-L	PolyPROPYL A® BioPureSPN MACRO PLATE 96 Well for HIC 50-200µl	HIC	Spin	100	300	0.2 -1mg	N/A	800
HEM HIC	PolyPROPYL A® BioPureSPN MIDI for HIC 50-200µl (Box of 100 each)	HIC	Spin	50	300	100-500µg	N/A	820
HEM HIC.20	PolyPROPYL A® BioPureSPN MIDI for HIC 50-200µl (Pack of 20 each)	HIC	Spin	50	300	100-500µg	N/A	240
HNS HIC-M	PolyPROPYL A® BioPureSPN MIDI PLATE 96 Well for HIC 50-200µl	HIC	Spin	50	300	100-500µg	N/A	575
HUM HIC	PolyPROPYL A® BioPureSPN MINI for HIC 20-100µl (Box of 100 each)	HIC	Spin	20	300	40-200µg	N/A	735
HUM HIC.20	PolyPROPYL A® BioPureSPN MINI for HIC 20-100µl (Pack of 20 each)	HIC	Spin	20	300	40-200µg	N/A	225
HNS HIC	PolyPROPYL A® BioPureSPN MINI PLATE 96 Well for HIC 20-100µl	HIC	Spin	20	300	40-200µg	N/A	735



Custom combinations available upon request info@polylc.com



Custom combinations available upon request info@polylc.com



Bulk Material

PolyPROPYL A® Bulk Material (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.) Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	QTY (g)	Particle (µm)	Pore (Å)	2026 Price (USD)
BMPR0203	PolyPROPYL A® Bulk Material, 2µm, 300Å	HIC	BM	N/A	0.5	2	300	100
BMPR0210	PolyPROPYL A® Bulk Material, 2µm, 1000Å	HIC	BM	N/A	0.5	2	1000	100
BMPR0215	PolyPROPYL A® Bulk Material, 2µm, 1500Å	HIC	BM	N/A	0.5	2	1500	100
BMPR0303	PolyPROPYL A® Bulk Material, 3µm, 300Å	HIC	BM	N/A	1	3	300	130
BMPR0310	PolyPROPYL A® Bulk Material, 3µm, 1000Å	HIC	BM	N/A	1	3	1000	130
BMPR0315	PolyPROPYL A® Bulk Material, 3µm, 1500Å	HIC	BM	N/A	1	3	1500	130
BMPR0503	PolyPROPYL A® Bulk Material, 5µm, 300Å	HIC	BM	N/A	1	5	300	85
BMPR0510	PolyPROPYL A® Bulk Material, 5µm, 1000Å	HIC	BM	N/A	1	5	1000	85
BMPR1203	PolyPROPYL A® Bulk Material, 12µm, 300Å	HIC	BM	N/A	1	12	300	32
BMPR1215	PolyPROPYL A® Bulk Material, 12µm, 1000Å	HIC	BM	N/A	1	12	1000	32
BMPR2003	PolyPROPYL A® Bulk Material, 20µm, 300Å	HIC	BM	N/A	1	20	300	32
BMPR2015	PolyPROPYL A® Bulk Material, 20µm, 1500Å	HIC	BM	N/A	1	20	1500	32
BMPR3003	PolyPROPYL A® Bulk Material, 30µm, 300Å	HIC	BM	N/A	1	30	300	32
BMPR3010	PolyPROPYL A® Bulk Material, 30µm, 1000Å	HIC	BM	N/A	1	30	1000	32

Please contact info@polylc.com with any questions or for information on any product combination not listed in the table above

Hydrophobic Interaction Chromatography

(HIC)

2026 PolyLC Inc. Retail Pricelist

PolyLC Inc., 9151 Rumsey Rd., Suite 180, Columbia, MD 21045 410-992-5400 www.polylc.com



2026 PolyLC Inc. Retail Pricelist

Mixed-Bed Chromatography IEX Resins




Offering high-performance liquid chromatography media, columns, capillaries, cartridges, and instrumentation. We specialize in advanced alternatives to reversed-phase for proteins, peptides, and other biomolecules, supported by deep application expertise and customized solutions for the most demanding separations.

Our resins can be blended within a column to permit retention/separation of negative and positive analytes in the same run.

Columns are in ratios to meet customers needs. Our most common combination is 1:1, but we are happy to make combinations like 2:1, 1:3 or anything you might be interested in. Both material identifiers are present in the part number and the ratio is added to the end.

For example, a 1 to 1 mixture of our weak cation exchange material and weak anion exchange material in a 200 x 4.6mm column with 5µm particles and 1000Å pores would be PN: 204CTWX0510_11. A combination containing twice as much WCX resin as WAX would be PN: 204CTWX0510_21



HPLC Columns

Mixed Bed IEX Columns (WCX-WAX) (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.)

Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
051CTWX0203_11	PolyCATWAX® Column, Glass Lined SS, 50 x 1.0mm, 2µm,300Å	MX-IEX	Analytical	GLT	10-32 coned	50	1.0	2	300	615
051CTWX0210_11	PolyCATWAX® Column, Glass Lined SS, 50 x 1.0mm, 2µm,1000Å	MX-IEX	Analytical	GLT	10-32 coned	50	1.0	2	1000	615
101CTWX0203_11	PolyCATWAX® Column, Glass Lined SS, 100 x 1.0mm, 2µm,300Å	MX-IEX	Analytical	GLT	10-32 coned	100	1.0	2	300	750
101CTWX0210_11	PolyCATWAX® Column, Glass Lined SS, 100 x 1.0mm, 2µm,1000Å	MX-IEX	Analytical	GLT	10-32 coned	100	1.0	2	1000	750
151CTWX0203_11	PolyCATWAX® Column, Glass Lined SS, 150 x 1.0mm, 2µm,300Å	MX-IEX	Analytical	GLT	10-32 coned	150	1.0	2	300	800
151CTWX0210_11	PolyCATWAX® Column, Glass Lined SS, 150 x 1.0mm, 2µm,1000Å	MX-IEX	Analytical	GLT	10-32 coned	150	1.0	2	1000	800
3.52CTWX0203_11	PolyCATWAX® Column , SS, 35 x 2.1mm, 2µm,300Å	MX-IEX	Analytical	SS	10-32 coned	35	2.1	2	300	570
3.52CTWX0210_11	PolyCATWAX® Column , SS, 35 x 2.1mm, 2µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	35	2.1	2	1000	570
052CTWX0203_11	PolyCATWAX® Column , SS, 50 x 2.1mm, 2µm,300Å	MX-IEX	Analytical	SS	10-32 coned	50	2.1	2	300	580
052CTWX0210_11	PolyCATWAX® Column , SS, 50 x 2.1mm, 2µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	50	2.1	2	1000	580
102CTWX0203_11	PolyCATWAX® Column , SS, 100 x 2.1mm, 2µm,300Å	MX-IEX	Analytical	SS	10-32 coned	100	2.1	2	300	750
102CTWX0210_11	PolyCATWAX® Column , SS, 100 x 2.1mm, 2µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	100	2.1	2	1000	750
3.54CTWX0203_11	PolyCATWAX® Column, SS, 35 x 4.6mm, 2µm,300Å	MX-IEX	Analytical	SS	10-32 coned	35	4.6	2	300	570
3.54CTWX0210_11	PolyCATWAX® Column, SS, 35 x 4.6mm, 2µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	35	4.6	2	1000	570
054CTWX0203_11	PolyCATWAX® Column, SS, 50 x 4.6mm, 2µm,300Å	MX-IEX	Analytical	SS	10-32 coned	50	4.6	2	300	580
054CTWX0210_11	PolyCATWAX® Column, SS, 50 x 4.6mm, 2µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	50	4.6	2	1000	580
104CTWX0203_11	PolyCATWAX® Column, SS, 100 x 4.6mm, 2µm,300Å	MX-IEX	Analytical	SS	10-32 coned	100	4.6	2	300	750
104CTWX0210_11	PolyCATWAX® Column, SS, 100 x 4.6mm, 2µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	100	4.6	2	1000	750

Mixed Bed IEX Columns (WCX-WAX) (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.) Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
154CTWX0203_11	PolyCATWAX® Column, SS, 150 x 4.6mm, 2µm,300Å	MX-IEX	Analytical	SS	10-32 coned	150	4.6	2	300	790
154CTWX0210_11	PolyCATWAX® Column, SS, 150 x 4.6mm, 2µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	150	4.6	2	1000	790
204CTWX0203_11	PolyCATWAX® Column, SS, 200 x 4.6mm, 2µm,300Å	MX-IEX	Analytical	SS	10-32 coned	200	4.6	2	300	810
204CTWX0210_11	PolyCATWAX® Column, SS, 200 x 4.6mm, 2µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	200	4.6	2	1000	810
051CTWX0303_11	PolyCATWAX® Column, Glass Lined SS, 50 x 1.0mm, 3µm,300Å	MX-IEX	Analytical	GLT	10-32 coned	50	1.0	3	300	570
051CTWX0310_11	PolyCATWAX® Column, Glass Lined SS, 50 x 1.0mm, 3µm,1000Å	MX-IEX	Analytical	GLT	10-32 coned	50	1.0	3	1000	570
101CTWX0303_11	PolyCATWAX® Column, Glass Lined SS, 100 x 1.0mm, 3µm,300Å	MX-IEX	Analytical	GLT	10-32 coned	100	1.0	3	300	660
101CTWX0310_11	PolyCATWAX® Column, Glass Lined SS, 100 x 1.0mm, 3µm,1000Å	MX-IEX	Analytical	GLT	10-32 coned	100	1.0	3	1000	660
151CTWX0303_11	PolyCATWAX® Column, Glass Lined SS, 150 x 1.0mm, 3µm,300Å	MX-IEX	Analytical	GLT	10-32 coned	150	1.0	3	300	735
151CTWX0310_11	PolyCATWAX® Column, Glass Lined SS, 150 x 1.0mm, 3µm,1000Å	MX-IEX	Analytical	GLT	10-32 coned	150	1.0	3	1000	735
3.52CTWX0303_11	PolyCATWAX® Column, SS, 35 x 2.1mm, 3µm,300Å	MX-IEX	Analytical	SS	10-32 coned	35	2.1	3	300	545
3.52CTWX0310_11	PolyCATWAX® Column, SS, 35 x 2.1mm, 3µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	35	2.1	3	1000	545
052CTWX0303_11	PolyCATWAX® Column, SS, 50 x 2.1mm, 3µm,300Å	MX-IEX	Analytical	SS	10-32 coned	50	2.1	3	300	550
052CTWX0310_11	PolyCATWAX® Column, SS, 50 x 2.1mm, 3µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	50	2.1	3	1000	550
102CTWX0303_11	PolyCATWAX® Column, SS, 100 x 2.1mm, 3µm,300Å	MX-IEX	Analytical	SS	10-32 coned	100	2.1	3	300	710
102CTWX0310_11	PolyCATWAX® Column, SS, 100 x 2.1mm, 3µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	100	2.1	3	1000	710
3.54CTWX0303_11	PolyCATWAX® Column, SS, 35 x 4.6mm, 3µm,300Å	MX-IEX	Analytical	SS	10-32 coned	35	4.6	3	300	545
3.54CTWX0310_11	PolyCATWAX® Column, SS, 35 x 4.6mm, 3µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	35	4.6	3	1000	545
054CTWX0303_11	PolyCATWAX® Column, SS, 50 x 4.6mm, 3µm,300Å	MX-IEX	Analytical	SS	10-32 coned	50	4.6	3	300	550
054CTWX0310_11	PolyCATWAX® Column, SS, 50 x 4.6mm, 3µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	50	4.6	3	1000	550
104CTWX0303_11	PolyCATWAX® Column, SS, 100 x 4.6mm, 3µm,300Å	MX-IEX	Analytical	SS	10-32 coned	100	4.6	3	300	710
104CTWX0310_11	PolyCATWAX® Column, SS, 100 x 4.6mm, 3µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	100	4.6	3	1000	710
154CTWX0303_11	PolyCATWAX® Column, SS, 150 x 4.6mm, 3µm,300Å	MX-IEX	Analytical	SS	10-32 coned	150	4.6	3	300	740
154CTWX0310_11	PolyCATWAX® Column, SS, 150 x 4.6mm, 3µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	150	4.6	3	1000	740
204CTWX0303_11	PolyCATWAX® Column, SS, 200 x 4.6mm, 3µm,300Å	MX-IEX	Analytical	SS	10-32 coned	200	4.6	3	300	750
204CTWX0310_11	PolyCATWAX® Column, SS, 200 x 4.6mm, 3µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	200	4.6	3	1000	750
254CTWX0303_11	PolyCATWAX® Column, SS, 250 x 4.6mm, 3µm,300Å	MX-IEX	Analytical	SS	10-32 coned	250	4.6	3	300	800
254CTWX0310_11	PolyCATWAX® Column, SS, 250 x 4.6mm, 3µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	250	4.6	3	1000	800
109CTWX0303_11	PolyCATWAX® Column, SS, 100 x 9.4mm, 3µm,300Å	MX-IEX	Analytical	SS	10-32 coned	100	9.4	3	300	1400
109CTWX0310_11	PolyCATWAX® Column, SS, 100 x 9.4mm, 3µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	100	9.4	3	1000	1400
209CTWX0303_11	PolyCATWAX® Column, SS, 200 x 9.4mm, 3µm,300Å	MX-IEX	Analytical	SS	10-32 coned	200	9.4	3	300	1620
209CTWX0310_11	PolyCATWAX® Column, SS, 200 x 9.4mm, 3µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	200	9.4	3	1000	1620
259CTWX0303_11	PolyCATWAX® Column, SS, 250 x 9.4mm, 3µm,300Å	MX-IEX	Analytical	SS	10-32 coned	250	9.4	3	300	1750
259CTWX0310_11	PolyCATWAX® Column, SS, 250 x 9.4mm, 3µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	250	9.4	3	1000	1750
051CTWX0303_11	PolyCATWAX® Column, Glass Lined SS, 50 x 1.0mm, 5µm,300Å	MX-IEX	Analytical	GLT	10-32 coned	50	1.0	5	300	570
051CTWX0310_11	PolyCATWAX® Column, Glass Lined SS, 50 x 1.0mm, 5µm,300Å	MX-IEX	Analytical	GLT	10-32 coned	50	1.0	5	1000	570
101CTWX0303_11	PolyCATWAX® Column, Glass Lined SS, 100 x 1.0mm, 5µm,300Å	MX-IEX	Analytical	GLT	10-32 coned	100	1.0	5	300	660
101CTWX0310_11	PolyCATWAX® Column, Glass Lined SS, 100 x 1.0mm, 5µm,300Å	MX-IEX	Analytical	GLT	10-32 coned	100	1.0	5	1000	660
151CTWX0303_11	PolyCATWAX® Column, Glass Lined SS, 150 x 1.0mm, 5µm,300Å	MX-IEX	Analytical	GLT	10-32 coned	150	1.0	5	300	740
151CTWX0310_11	PolyCATWAX® Column, Glass Lined SS, 150 x 1.0mm, 5µm,300Å	MX-IEX	Analytical	GLT	10-32 coned	150	1.0	5	1000	740
3.52CTWX0503_11	PolyCATWAX® Column, SS, 35 x 2.1mm, 5µm,300Å	MX-IEX	Analytical	SS	10-32 coned	35	2.1	5	300	500
3.52CTWX0510_11	PolyCATWAX® Column, SS, 35 x 2.1mm, 5µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	35	2.1	5	1000	500
052CTWX0503_11	PolyCATWAX® Column, SS, 50 x 2.1mm, 5µm,300Å	MX-IEX	Analytical	SS	10-32 coned	50	2.1	5	300	520
052CTWX0510_11	PolyCATWAX® Column, SS, 50 x 2.1mm, 5µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	50	2.1	5	1000	520
102CTWX0503_11	PolyCATWAX® Column, SS, 100 x 2.1mm, 5µm,300Å	MX-IEX	Analytical	SS	10-32 coned	100	2.1	5	300	600
102CTWX0510_11	PolyCATWAX® Column, SS, 100 x 2.1mm, 5µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	100	2.1	5	1000	600
152CTWX0503_11	PolyCATWAX® Column, SS, 150 x 2.1mm, 5µm,300Å	MX-IEX	Analytical	SS	10-32 coned	150	2.1	5	300	650

Mixed Bed IEX Columns (WCX-WAX) (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.) Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
152CTWX0510_11	PolyCATWAX® Column, SS, 150 x 2.1mm, 5µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	150	2.1	5	1000	650
202CTWX0503_11	PolyCATWAX® Column, SS, 200 x 2.1mm, 5µm,300Å	MX-IEX	Analytical	SS	10-32 coned	200	2.1	5	300	710
202CTWX0510_11	PolyCATWAX® Column, SS, 200 x 2.1mm, 5µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	200	2.1	5	1000	710
3.54CTWX0503_11	PolyCATWAX® Column, SS, 35 x 4.6mm, 5µm,300Å	MX-IEX	Analytical	SS	10-32 coned	35	4.6	5	300	500
3.54CTWX0510_11	PolyCATWAX® Column, SS, 35 x 4.6mm, 5µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	35	4.6	5	1000	500
054CTWX0503_11	PolyCATWAX® Column, SS, 50 x 4.6mm, 5µm,300Å	MX-IEX	Analytical	SS	10-32 coned	50	4.6	5	300	520
054CTWX0510_11	PolyCATWAX® Column, SS, 50 x 4.6mm, 5µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	50	4.6	5	1000	520
104CTWX0503_11	PolyCATWAX® Column, SS, 100 x 4.6mm, 5µm,300Å	MX-IEX	Analytical	SS	10-32 coned	100	4.6	5	300	600
104CTWX0510_11	PolyCATWAX® Column, SS, 100 x 4.6mm, 5µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	100	4.6	5	1000	600
154CTWX0503_11	PolyCATWAX® Column, SS, 150 x 4.6mm, 5µm,300Å	MX-IEX	Analytical	SS	10-32 coned	150	4.6	5	300	650
154CTWX0510_11	PolyCATWAX® Column, SS, 150 x 4.6mm, 5µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	150	4.6	5	1000	650
204CTWX0503_11	PolyCATWAX® Column, SS, 200 x 4.6mm, 5µm,300Å	MX-IEX	Analytical	SS	10-32 coned	200	4.6	5	300	710
204CTWX0510_11	PolyCATWAX® Column, SS, 200 x 4.6mm, 5µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	200	4.6	5	1000	710
254CTWX0503_11	PolyCATWAX® Column, SS, 250 x 4.6mm, 5µm,300Å	MX-IEX	Analytical	SS	10-32 coned	250	4.6	5	300	740
254CTWX0510_11	PolyCATWAX® Column, SS, 250 x 4.6mm, 5µm,1000Å	MX-IEX	Analytical	SS	10-32 coned	250	4.6	5	1000	740
109CTWX0503_11	PolyCATWAX® Column, SS, 100 x 9.4mm, 5µm,300Å	MX-IEX	Semi-Prep	SS	10-32 coned	100	9.4	5	300	1110
109CTWX0510_11	PolyCATWAX® Column, SS, 100 x 9.4mm, 5µm,1000Å	MX-IEX	Semi-Prep	SS	10-32 coned	100	9.4	5	1000	1110
209CTWX0503_11	PolyCATWAX® Column, SS, 200 x 9.4mm, 5µm,300Å	MX-IEX	Semi-Prep	SS	10-32 coned	200	9.4	5	300	1400
209CTWX0510_11	PolyCATWAX® Column, SS, 200 x 9.4mm, 5µm,1000Å	MX-IEX	Semi-Prep	SS	10-32 coned	200	9.4	5	1000	1400
259CTWX0503_11	PolyCATWAX® Column, SS, 250 x 9.4mm, 5µm,300Å	MX-IEX	Semi-Prep	SS	10-32 coned	250	9.4	5	300	1500
259CTWX0510_11	PolyCATWAX® Column, SS, 250 x 9.4mm, 5µm,1000Å	MX-IEX	Semi-Prep	SS	10-32 coned	250	9.4	5	1000	1500
2521CTWX0503_11	PolyCATWAX® Column, SS, 250 x 21mm, 5µm,300Å	MX-IEX	Semi-Prep	SS	10-32 coned	250	21.0	5	300	4500
2521CTWX0510_11	PolyCATWAX® Column, SS, 250 x 21mm, 5µm,1000Å	MX-IEX	Semi-Prep	SS	10-32 coned	250	21.0	5	1000	4500
204CTWX1203_11	PolyCATWAX® Column, SS, 200 x 4.6mm, 12µm,300Å	MX-IEX	Analytical	SS	10-32 coned	200	4.6	12	300	420
204CTWX1215_11	PolyCATWAX® Column, SS, 200 x 4.6mm, 12µm,1500Å	MX-IEX	Analytical	SS	10-32 coned	200	4.6	12	1500	420
254CTWX1203_11	PolyCATWAX® Column, SS, 250 x 4.6mm, 12µm,300Å	MX-IEX	Analytical	SS	10-32 coned	250	4.6	12	300	500
254CTWX1215_11	PolyCATWAX® Column, SS, 250 x 4.6mm, 12µm,1500Å	MX-IEX	Analytical	SS	10-32 coned	250	4.6	12	1500	500
209CTWX1203_11	PolyCATWAX® Column, SS, 200 x 9.4mm, 12µm,300Å	MX-IEX	Semi-Prep	SS	10-32 coned	200	9.4	12	300	1000
209CTWX1215_11	PolyCATWAX® Column, SS, 200 x 9.4mm, 12µm,1500Å	MX-IEX	Semi-Prep	SS	10-32 coned	200	9.4	12	1500	1000
259CTWX1203_11	PolyCATWAX® Column, SS, 250 x 9.4mm, 12µm,300Å	MX-IEX	Semi-Prep	SS	10-32 coned	250	9.4	12	300	1100
259CTWX1215_11	PolyCATWAX® Column, SS, 250 x 9.4mm, 12µm,1500Å	MX-IEX	Semi-Prep	SS	10-32 coned	250	9.4	12	1500	1100
2521CTWX1203_11	PolyCATWAX® Column, SS, 250 x 21mm, 12µm,300Å	MX-IEX	Semi-Prep	SS	10-32 coned	250	21.0	12	300	3200
2521CTWX1215_11	PolyCATWAX® Column, SS, 250 x 21mm, 12µm,1500Å	MX-IEX	Semi-Prep	SS	10-32 coned	250	21.0	12	1500	3200
2550CTWX1203_11	PolyCATWAX® Column, SS, 250 x 50.8mm, 12µm,300Å	MX-IEX	Prep	SS	10-32 coned	250	50.8	12	300	12250
2550CTWX1215_11	PolyCATWAX® Column, SS, 250 x 50.8mm, 12µm,1500Å	MX-IEX	Prep	SS	10-32 coned	250	50.8	12	1500	12250



HPLC Capillaries

Mixed Bed IEX Capillaries (WCX-WAX) (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.) Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
050.075CTWX0203_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 2µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	50	0.075	2	300	520
050.075CTWX0210_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 2µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	50	0.075	2	1000	520
100.075CTWX0203_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 2µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	100	0.075	2	300	580
100.075CTWX0210_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 2µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	100	0.075	2	1000	580
150.075CTWX0203_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 2µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	150	0.075	2	300	640
150.075CTWX0210_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 2µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	150	0.075	2	1000	640
200.075CTWX0203_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 2µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	200	0.075	2	300	700
200.075CTWX0210_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 2µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	200	0.075	2	1000	700
050.15CTWX0203_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 2µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	50	0.150	2	300	520
050.15CTWX0210_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 2µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	50	0.150	2	1000	520
100.15CTWX0203_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 2µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	100	0.150	2	300	580
100.15CTWX0210_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 2µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	100	0.150	2	1000	580
150.15CTWX0203_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 2µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	150	0.150	2	300	640
150.15CTWX0210_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 2µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	150	0.150	2	1000	640
200.15CTWX0203_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 2µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	200	0.150	2	300	700
200.15CTWX0210_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 2µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	200	0.150	2	1000	700
050.30CTWX0203_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 2µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	50	0.300	2	300	520
050.30CTWX0210_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 2µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	50	0.300	2	1000	520
100.30CTWX0203_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 2µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	100	0.300	2	300	580
100.30CTWX0210_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 2µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	100	0.300	2	1000	580
150.30CTWX0203_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 2µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	150	0.300	2	300	640
150.30CTWX0210_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 2µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	150	0.300	2	1000	640
200.30CTWX0203_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 2µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	200	0.300	2	300	700
200.30CTWX0210_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 2µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	200	0.300	2	1000	700
050.075CTWX0303_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	50	0.075	3	300	520
050.075CTWX0310_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 3µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	50	0.075	3	1000	520
100.075CTWX0303_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	100	0.075	3	300	580
100.075CTWX0310_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 3µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	100	0.075	3	1000	580
150.075CTWX0303_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	150	0.075	3	300	640
150.075CTWX0310_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 3µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	150	0.075	3	1000	640
200.075CTWX0303_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	200	0.075	3	300	700
200.075CTWX0310_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 3µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	200	0.075	3	1000	700
050.15CTWX0303_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	50	0.150	3	300	520
050.15CTWX0310_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 3µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	50	0.150	3	1000	520
100.15CTWX0303_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	100	0.150	3	300	580
100.15CTWX0310_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 3µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	100	0.150	3	1000	580
150.15CTWX0303_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	150	0.150	3	300	640
150.15CTWX0310_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 3µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	150	0.150	3	1000	640
200.15CTWX0303_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	200	0.150	3	300	700
200.15CTWX0310_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 3µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	200	0.150	3	1000	700
050.30CTWX0303_11	PolyCATWAX A® Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	50	0.300	3	300	520

Mixed Bed IEX Capillaries (WCX-WAX) (Available in 2, 3, 5, 12, 20, and 30µm particles with multiple pore sizes.) Custom combinations available upon request info@polylc.com

PN	Product Description	Phase	Type	Hardware	Connection	L (mm)	ID (mm)	Particle (µm)	Pore (Å)	2026 Price (USD)
050.30CTWX0310_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 3µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	50	0.300	3	1000	520
100.30CTWX0303_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	100	0.300	3	300	580
100.30CTWX0310_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 3µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	100	0.300	3	1000	580
150.30CTWX0303_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	150	0.300	3	300	640
150.30CTWX0310_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 3µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	150	0.300	3	1000	640
200.30CTWX0303_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	200	0.300	3	300	700
200.30CTWX0310_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 3µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	200	0.300	3	1000	700
050.075CTWX0503_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 5µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	50	0.075	5	300	520
050.075CTWX0510_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 50mm x 75µm, 5µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	50	0.075	5	1000	520
100.075CTWX0503_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 5µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	100	0.075	5	300	580
100.075CTWX0510_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 100mm x 75µm, 5µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	100	0.075	5	1000	580
150.075CTWX0503_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 5µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	150	0.075	5	300	640
150.075CTWX0510_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 150mm x 75µm, 5µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	150	0.075	5	1000	640
200.075CTWX0503_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 5µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	200	0.075	5	300	700
200.075CTWX0510_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 200mm x 75µm, 5µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	200	0.075	5	1000	700
050.15CTWX0503_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 5µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	50	0.150	5	300	520
050.15CTWX0510_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 50mm x 150µm, 5µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	50	0.150	5	1000	520
100.15CTWX0503_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 5µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	100	0.150	5	300	580
100.15CTWX0510_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 100mm x 150µm, 5µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	100	0.150	5	1000	580
150.15CTWX0503_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 5µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	150	0.150	5	300	640
150.15CTWX0510_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 150mm x 150µm, 5µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	150	0.150	5	1000	640
200.15CTWX0503_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 5µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	200	0.150	5	300	700
200.15CTWX0510_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 200mm x 150µm, 5µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	200	0.150	5	1000	700
050.30CTWX0503_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 5µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	50	0.300	5	300	520
050.30CTWX0510_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 50mm x 300µm, 5µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	50	0.300	5	1000	520
100.30CTWX0503_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 5µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	100	0.300	5	300	580
100.30CTWX0510_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 100mm x 300µm, 5µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	100	0.300	5	1000	580
150.30CTWX0503_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 5µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	150	0.300	5	300	640
150.30CTWX0510_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 150mm x 300µm, 5µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	150	0.300	5	1000	640
200.30CTWX0503_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 5µm,300Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	200	0.300	5	300	700
200.30CTWX0510_11	PolyCATWAX A* Capillary, Glass Lined, PEEK SS, 200mm x 300µm, 5µm,1000Å	MX-IEX	Analytical	PEEK SIL	10-32 coned	200	0.300	5	1000	700

Please contact info@polylc.com with any questions or for information on any product combination not listed in the table above

Mixed-Bed Chromatography

(Mixed-Bed IEX)

2026 PolyLC Inc. Retail Pricelist

PolyLC Inc., 9151 Rumsey Rd., Suite 180, Columbia, MD 21045 410-992-5400 www.polylc.com

